STORZ

KARL STORZ — ENDOSKOPE

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Telescopes and HOPKINS® telescopes

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**Important information for
users of KARL STORZ
devices and instruments**

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**It is recommended to check the suitability
of the products for the intended procedure
prior to use.**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is designed to provide help with the proper application, cleaning, and sterilization of the telescopes manufactured by KARL STORZ. All required details and all actions on your part are clearly presented and explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

The manual also describes the recommended testing and reprocessing methods for telescopes.

WARNING: KARL STORZ telescopes are not sterile when delivered, and must be cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

KARL STORZ is continuously working on the further development of all products. Please understand that changes to the scope of supply, design, equipment, and technology are possible for this reason. Therefore, no claims may be derived from the information, figures, and descriptions in this manual.

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







- ① Screw base for KARL STORZ fiberoptic light cables and Olympus/Winter & Ibe
- ② Plug-in base, 9 mm diameter, for Wolf fiberoptic light cables
- ③ Light inlet piece without adaptor
- ④ Color-code ring to indicate the telescope's direction of view:
 - Green = 0 degrees, straight forward
 - Black = 12 degrees/45 degrees
 - Red = 30 degrees, forward oblique
 - Yellow = 70 degrees, lateral
 - Blue = 90 degrees, lateral
 - White = 120 degrees, rear view
 - Gold = 0° – 120° (ENDOCAMELEON® size 10 mm)
 - Gold = 15° – 90° (ENDOCAMELEON® size 4 mm)
- ⑤ Eyepiece
- ⑥ Sheath
- ⑦ 'Autoclave' marking ring
- ⑧ Type designation and direction of view. The first letter after the number means:
 - A = 0/6 degrees, straight forward
 - B = 30 degrees, forward oblique
 - C = 70 degrees, lateral
 - D = 90 degrees, lateral
 - E = 120 degrees, rear view
 - F = 12 degrees/45 degrees
 If present, the letter W stands for wide angle.
- ⑨ Lens and distal outlet of light cable
- ⑩ Focusing wheel for the eyepiece drive
- ⑪ DCI® connection
- ⑫ Adjustment wheel for direction of view
- ⑬ Rotating wheel for filter replacement (possible with telescopes for fluorescence diagnostics)
- ⑭ Illuminator (possible for VITOM® telescopes)
- ⑮ Irrigation port (LUER-Lock)
- ⑯ Working channel

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1 General information

HOPKINS® rod lens telescopes are reusable and consist of an eyepiece, a connection for light cables with screwed-on adaptors for light cables from other manufacturers, and a sheath made of rust-proof material that encloses the rod lens system and an integrated fiberoptic light cable.

HOPKINS® telescopes are designed to keep the strain on a patient to a minimum. They are put to use hundreds of thousands of times every year.

However, HOPKINS® telescopes are much more sensitive than the devices that a physician would usually use.

They must therefore be handled with care, and checked for correct functioning and for signs of damage before and after every use on a patient. Improper handling or incorrect usage may lead to tissue damage, infection, or instrument breakages.

Observing the following information will ensure that your HOPKINS® telescopes have a long service life.

i **NOTE:** The telescopes referred to in this manual usually have separate instruction manuals. Other general regulations (intended use, user qualification, safety rules etc.) must be taken from these separate instruction manuals and observed!

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2 Safety instructions

2.1 Symbols and hazard information

Please read this manual carefully and follow its instructions exactly. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the following text should be carefully reviewed to ensure the safe and effective use of this instrument. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

2.2 Safety instructions

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces,

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sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light connection, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

Avoid direct tissue contact with the distal end, as well as with the light connection of the telescope and the light cable. Do not rest the endoscope or the light cable on the patient or where it is in direct contact with operating accessories. Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope. Always wear standardized protective equipment when using lasers.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with an applied part of type CF is used. In this case, an endotherapy device with a type CF applied part must also be used in order to minimize the total patient leakage current.

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WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers, it must be ensured that the proximal end (= furthest from the patient) of the telescope is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the setting ring, which can cause cardiac pacemakers or other electronic devices to malfunction.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

NOTE: VITOM® telescopes are not endoscopes and must therefore only be used outside the body. The minimum working distance is 20 cm, and therefore significantly greater than with an endoscope. The minimum distance is set using a spacer bar (209180 20), which is available separately.

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2.3 Intended use

Indication: KARL STORZ HOPKINS® telescopes are used for visualization of the operating field during diagnostic and surgical endoscopic interventions.

VITOM® telescopes are used for the extracorporeal display of the operating field during diagnostic and therapeutic procedures in endoscopic and open surgery.

Contraindication: Use of KARL STORZ telescopes and HOPKINS® telescopes is contraindicated if, in the opinion of the responsible physician, such an application would endanger the patient.

Unauthorized conversions or modifications to the instrument are not permitted for safety reasons.

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2.4 User qualification

KARL STORZ telescopes and HOPKINS® telescopes may only be used by persons with an appropriate medical qualification and who are acquainted with endoscopic techniques. The information given in this manual only serves to instruct in the correct handling, cleaning and sterilization of the telescopes. It is not intended as an introduction to the techniques of endoscopy.

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2.5 Safety precautions when using the telescopes

The telescopes must only be used according to the medical rules and procedures of endoscopy recognized for the method.

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3 Handling**3.1 Unpacking**

Check for missing items and evidence of shipping damage.
File any complaints immediately with KARL STORZ or the supplier.

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3.2 Application

WARNING: The telescope is not sterile when delivered and must therefore be cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

WARNING: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely.

WARNING: Patient leakage currents may accumulate if endoscopes with powered accessories are used. This is particularly relevant when using CF devices. In such a case, the other devices used should also be of type CF in order to minimize the total patient leakage current.

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3.3 Inspection

Inspect the telescope for damage immediately after receipt and before and after each use.

Image waveguide

Examine the optical end faces (distal end and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny.

In order to test the image quality, slowly rotate the telescope while looking through it.

If the image is completely or partially impaired, either the rod lens system in the sheath is broken or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

A cloudy or speckled image can be caused by moisture or disinfectant residue on the optical end faces.

i NOTE: When looking through HOPKINS® PDD telescopes, the image appears yellowish. This effect is normal and is caused by the use of special filters.

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
In these cases, carefully clean the optical end faces with a cotton tip applicator soaked in 70 % alcohol.

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Light cable

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection. Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image.

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However, the light output does significantly deteriorate if a certain number are broken (from approx. 20 % to 25 %), which therefore also results in reduced image quality. The telescope must not be used any longer.

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NOTE: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

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
Sheath

Inspect the entire sheath for dents and scratches. The sheath must not be bent. LUER-Lock stopcocks, telescope locks and other movable parts should be smooth-running and undamaged. Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This applies in particular to seals.

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3.4 Precautions


As with all optical instruments, the HOPKINS® telescopes are a very sensitive part of the endoscopic equipment. In order to ensure the long-term quality of the telescope, the following information must be observed.

- In order to avoid breaks, telescopes must always be held by the eyepiece and never by the distal end alone.
- Do not bend the sheath under any circumstances. This can result in breaks or cracks in the rod lens system.
- Place the telescopes down with care. Hard impacts, particularly to the distal end, can cause damage to or cracks in the luting and allow liquid, steam, etc. to penetrate. This results in cloudy and/or blurred image areas.
- Therefore, during cleaning, disinfection, sterilization and storage, it must be ensured that the telescopes cannot be damaged by other instruments. Ideally, you should store HOPKINS® telescopes individually or use containers in which the telescopes can be secured (see accessories).

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Special protective sheaths are available in different lengths for the HOPKINS® telescopes. These sheaths are suitable for storage and transport, as well as for sterilization (see accessories).

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CAUTION: HOPKINS® telescopes that do not have the *Autoclave* identification ring must not be steam sterilized. This may cause irreparable damage.

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3.5 Instructions for use for accessories

Before use, make sure that the telescope has been correctly cleaned, disinfected, and sterilized where necessary.
Check that the telescope and accessories are all functioning correctly before use.
Operate the video and light source equipment as per the specifications in the relevant manual.

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
To prevent the telescopes from fogging during use, we recommend the use of a warmer for rigid telescopes, Art. no. 10905 B (see page 35).

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
3.6 Connecting to the camera system and light source

The camera heads generally have an integrated instrument coupler.
Squeeze the levers (1) and (2) of the instrument coupler together and insert the endoscope's eyepiece into the coupler. The coupler automatically locks into place when the eyepiece is inserted.
The endoscope can be locked by pressing lever (2) to lever (3).

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Connect the light cable to the outlet of the light source and to the endoscope (twist the knurled screw on the screw base by a quarter turn).

! WARNING: Reduce the intensity of the light source so that optimal illumination of the operating field is still achieved. The higher the intensity of the light source, the greater the heat emission at the distal end of the telescope.

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WARNING: Never place the end of the fiberoptic light cable or a telescope on or underneath a sterile cover cloth. High light intensity can cause burns to the patient and/or the sterile cover cloth.

WARNING: Danger of glare! Never look into the free end of a light cable.

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3.7 Light cable assignment

If there is an identification symbol on the light cable and the endoscope (see table below), assign matching combinations according to the symbols shown*.

Light cable diameter	Endoscope diameter
4.8 – 5.0 mm 	6.6 – 12.0 mm 
3.0 – 3.5 mm 	3.0 – 6.5 mm 
2.0 – 2.5 mm 	0.8 – 2.9 mm 


* Special endoscopes may deviate from this

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
3.8 Focus and image size

Operate the video equipment as per the specifications in the relevant manual. Perform a new white balance as required, or when replacing the telescope.

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3.9 Assembling the sheath and telescope

Sheaths, sheath inserts and telescopes are provided with a 0 mark or a line mark.

Open the telescope lock ① by turning it counterclockwise (as far as it will go, marks in alignment ②).

When assembling, make sure that the marks are aligned with one another (line mark and 0 mark ③ are lined up). Insert the telescope into the sheath lock and close the telescope lock ① by turning it in a clockwise direction.

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4 Reprocessing

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

4.1 Accessories

Accessories required for carrying out reprocessing:

Brushes:

Surface: 27652
Inside lumina 2.5 – 5 mm: 27650 B
Inside lumina 2 – 2.4 mm: 27650 C

4.2 Disassembly

Prior to cleaning and disinfection, the telescope must be separated into its individual components as far as possible. Both adaptors on the light inlet piece (Fig. A), and any adaptors (Fig. B), seals on the working channel and LUER-Lock stopcocks (Fig. C) must be disassembled.

4.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the telescope immediately after use. To this end, preclean the telescope surface by wiping down and rinsing, for example, and rinse the working channel if necessary. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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

C

4.4 Manual precleaning Brushing the surfaces

Depending on the degree of soiling and the degree to which residues have dried, heavy soiling must be cleaned from the surfaces of instruments with the aid of a brush (Art. no. 27652)/sponge under cold running water until any visible contamination has been completely removed (Fig. A).

Brushing the lumina

In order to remove all visible soiling, the working channel (Fig. B), lumina and hollow spaces (e.g. disassembled LUER-Lock stopcocks, Fig. C) of the telescope must be precleaned under cold running water with the aid of a suitable brush. The brush with Art. no. 27650 B is suitable for the LUER-Lock stopcock.

No pretreatment with ultrasound

For technical reasons, the telescope is **not** suitable for ultrasound treatment.

4.5 Manual cleaning

The telescope and its dismantled individual components must be completely immersed in a cleaning solution. By filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with the aid of a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

The dismantled individual components must be completely immersed in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a cloth or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

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4.6 Manual disinfection

The telescope and its dismantled individual components must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations (Fig. A). The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

4.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off afterwards by hand.

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4.8 Connecting

In order to guarantee effective machine cleaning and disinfection, the instrument must be connected up to the washer and disinfector (e.g. via adaptor), or it must be ensured that it is thoroughly rinsed out by means of appropriate rinsing lances (Fig. A).

4.9 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

4.9.1 Removal of deposits

Due to foreign particles (inert gases) in the steam, over time stubborn deposits may form on the cover slips during autoclaving. These deposits can be easily removed using the cleaning paste (Art. no. 27661).

Put a little cleaning paste on a clean, moistened cleaning rod and apply to the cover window (Fig. B). Polishing removes stubborn deposits which can be subsequently rinsed off with water.

In order to achieve a smooth surface, polishing wool (Art. no. 27658) can be used if necessary.

! WARNUNG: Following cleaning, a complete reprocessing cycle must be carried out.

! CAUTION: The optionally available special cleaning paste should only be used for cleaning when the endoscopic image is visibly cloudy; under no circumstances should the paste be used on a regular basis.

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
D

4.9.2 Checking the image waveguides

Examine the optical end faces (distal tip and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny. In order to test the image quality, slowly rotate the telescope while looking through it. If the image is completely or partially impaired, either the rod lens system in the sheath or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

A cloudy or speckled image can be caused by moisture or disinfectant/chemical residue on the optical end faces. In these cases, carefully clean the optical end faces with a cotton tip applicator soaked in 70 % alcohol (Fig. A).

4.9.3 Checking the light cables

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection (Fig. B). Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image. However, the light output does significantly deteriorate if a certain number are broken (from approx. 20 % to 25 %), which therefore also results in reduced image quality. The telescope must not be used any longer.

! WARNING: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

4.9.4 Drying the optical end faces

The fiber areas in the light inlet piece in particular as well as the distal and proximal end faces must be dried carefully with 70 % alcohol (Fig. C and D). Disinfectant and cleaning agent residue, e.g., in the light inlet piece, can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

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- Damaged or corroded medical devices must be withdrawn from use.
- Threads and glide surfaces on the telescope must be lubricated locally with instrument oil (Art. no. 27656 B) and the running surfaces on the cock plug with special grease (Art. no. 27657) (Fig. A).
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

- ❗ **NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).
- ❗ **NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

4.10 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953, Fig. B).

4.11 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective national requirements and in consultation with the device manufacturers.

Telescopes with a ring marked 'AUTOCLAVE' can be steam sterilized without restrictions (Fig. C).

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C – 137 °C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

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CAUTION: Do not open sterilization containers immediately after the sterilization process has finished. Never pour cold sterile water over the telescopes. If the cooling process is accelerated, the telescopes may be damaged.

NOTE: Greased LUER-Lock stopcocks must be dismantled for sterilization so that steam penetration can be ensured.

Steam sterilization is not permitted for the following instrument groups. Only those methods with a maximum temperature of 65 °C may be used, such as gas and plasma sterilization:

- Non-autoclavable telescopes (cleaning up to 93 °C)
- Non-autoclavable miniature endoscopes (cleaning up to 65 °C)
- Magnifiers with drive mechanism
- Non-autoclavable prismatic light deflectors

4.12 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Medical devices must be dismantled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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 **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for telescopes **without** a channel:

- STERRAD® 50, 100S, 200
- STERRAD® NX®
- STERRAD® 100NX®
- STERRAD® 100NX® DUO Cycle

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for telescopes **with** a channel:

- STERRAD® NX®
- STERRAD® 100NX®

4.13 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™1 (only for telescopes without a channel)

Detailed information on the selection of the appropriate cycle for the various device generations is available from the manufacturer STERIS®.

 **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

4.14 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30 – 45 minutes.

 **WARNING:** Medical devices must be dismantled for sterilization.

 **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed for instruments, due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

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4.15 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1® (only for telescopes without a channel)

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

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4.16 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E® (only for telescopes without a channel)

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

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The following procedures are approved in respect of material compatibility, but must be validated by the on-site operator with regard to efficacy:

Low-temperature steam and formaldehyde process (LTSF)

i **NOTE:** Sterilization with FO is tested and approved by KARL STORZ in respect of material compatibility. However, due to the large number of processes used in the global market, validation in respect of efficacy was not carried out. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer. KARL STORZ has had exclusively the aforementioned sterilization methods validated with regard to a reduction in microorganisms of the order of 10^{-6} SAL.

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i **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

4.17 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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- This reprocessing step is part of the validated process
 - This reprocessing step may be carried out additionally in order to aid cleaning
 - x On the grounds of potential damage, this reprocessing step must NEVER be performed
- Fields which are not filled are not part of the performed validation.

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5 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine parts only.

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6 Limitation of liability

KARL STORZ GmbH & Co. KG shall be liable for the failure or a deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to the same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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7 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see Chapter 'Subsidiaries'), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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8 Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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9 Repair program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the replacement plan, you will receive a new replacement telescope and will only be charged the repair costs for the defective telescope.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Service department
Logistics center
Take-Off Gewerbepark 83
78579 Neuhausen ob Eck, Germany

In other countries please contact your local KARL STORZ subsidiary or authorized dealer.

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10 Important information

For hygiene reasons and to prevent infection, telescopes, instruments and devices must be cleaned and disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender. Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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Cleaning, sterilization, storage

Article

Tray holder for cleaning, sterilization, and storage of one rigid endoscope. With holders for light connection adaptors, silicone telescope holders and lid.
External dimensions (W x D x H): 290 x 60 x 52 mm.
For rigid endoscopes with up to 5 mm diameter and 20 cm working length.

Same, external dimensions (W x D x H): 430 x 65 x 52 mm, for use with rigid endoscopes with up to 10 mm diameter and 34 cm working length

Tray holder with lid, for 1 long endoscope, 670 x 80 x 52 mm

Tray holder for cleaning, sterilization, and storage of two rigid endoscopes and one fiberoptic light cable. With holders for adaptors, silicone telescope holder and lid.
External dimensions (W x D x H): 352 x 125 x 54 mm.
For rigid endoscopes with up to 10 mm diameter and 20 cm working length.

Same, external dimensions (W x D x H): 487 x 125 x 54 mm. For rigid endoscopes with up to 10 mm diameter and 32 cm working length.

Tray holder for cleaning, sterilization, and storage of one rigid endoscope with angled eyepiece, including holders for light connection adaptors, silicone telescope holders and lid, external dimensions (W x D x H): 460 x 80 x 52 mm

Tray holder for cleaning, sterilization, and storage of one ENDOCAMELEON® and one fiberoptic light cable, with holder for adaptors, silicone telescope holder and lid.
External dimensions (W x D x H): 480 x 125 x 54 mm

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


KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: lcano@karlstorz.com.mx

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsättrågränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlyoftens 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: info@karlstorz.in

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 Hong Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenheng District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang 1, 102dong, 8F
239 Jeongjail-ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ GmbH & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE

 **KARL STORZ GmbH & Co. KG**

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com