STORZ

KARL STORZ — ENDOSKOPE

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20 1612 20 Kaltlicht-Fontäne LED nova 150

INSTRUCTION MANUAL

20 1612 20 Cold Light Fountain LED nova 150

MANUAL DE INSTRUCCIONES

20 1612 20 Fuente de luz fría LED nova 150

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**1. Important information
for users of KARL STORZ
instruments**

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the cold light fountain LED nova 150. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG

Utilization and reproduction by third parties require express authorization by KARL STORZ SE & Co. KG.

All rights reserved.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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
**3. Controls, displays,
connectors, and
their uses**

- ① Power switch, lights up when in operation
- ② Light outlet point (not exchangeable)
- ③ Brightness control knob
- ④ Air vents
- ⑤ Potential equalization connector
- ⑥ Fuse holder
- ⑦ Power cord socket
- ⑧ to ⑫ - see next page (symbols)

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4. Symbols employed

	Follow instructions for use (see ⑧ on page IV)
	Brightness (see ⑨ on page IV)
	Light outlet (see ⑪ on page IV)
	Type CF applied part according to IEC 60601-1 (Explanation: Limitation of leakage currents; CF = Cardiac Floating. See ⑩ on page IV).
	Potential equalization connector (⑫ on page IV)
	Alternating current
	Electronic information product pollution control (China RoHS). The symbol is found on the identification plate on the underside of the unit.
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE). The symbol is found on the identification plate on the underside of the unit.
	Manufacturer

Symbols on label and packaging of the device

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from www.karlstorz.com.

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1. Important information for users of KARL STORZ instruments.....	III
2. Images of the device	IV
3. Controls, displays, connectors, and their uses	V
4. Symbols employed.....	VI
5. General information	3
5.1 Device description	3
6. Safety instructions.....	4
6.1 Explanation of warnings and cautions.....	4
6.2 Normal use.....	9
6.2.1 Indication	9
6.2.2 Contraindication.....	9
6.3 Intended use.....	9
6.4 User qualification	9
6.5 Patient profile.....	10
6.6 User profile	10
6.7 Intended conditions of use.....	10
6.7.1 Use.....	10
6.7.2 Other intended conditions.....	11
6.8 Safety precautions at the site of installation	11
6.9 Functions and features	12
7. Installation and operating instructions.....	13
7.1 Unpacking the equipment.....	13
7.2 Basic equipment.....	13
7.3 Installation.....	13
7.3.1 Connecting the ground line.....	14
7.3.2 Connecting the power cord	14
7.3.3 Connecting light cable	15
7.3.4 Connecting the light cable to the endoscope or videoscope	15
7.4 Commissioning the LED nova 150 light source.....	16
7.4.1 Setting the brightness.....	16
7.4.2 Decommissioning	16
7.5 Safety functions of the equipment	17

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5. General information

5.1 Device description

The LED nova 150 is a cold light source for diagnostic and therapeutic endoscopic applications.

The brightness of the light is adjusted with the control knob ③ on the front panel.

The main advantages of this LED light source are:

- Efficiency: High light intensity at low power consumption.
- Comfortable working environment: Intelligent cooling concept for minimal noise emission.
- Economic: there are normally no costs for a replacement lamp and no expenses for exchanging the lamp.

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6. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

NOTE: Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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WARNING: Read this instruction manual thoroughly before putting the device into operation. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

WARNING: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or Chap. 16, IEC 60601-1).

WARNING: Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

WARNING: At the end of its useful operating life, dispose of the unit as electronic scrap.

WARNING: Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

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Safety instructions Warnings and cautions

WARNING: Ensure sufficient air circulation for the device.

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WARNING: DO NOT use the device in the presence of flammable anesthetics.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: Test this device prior to each procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

WARNING: Do not use the device if the brightness regulator is not working properly.

WARNING: The device must not be used in oxygenated environments.

WARNING: The device may only be operated within specified environmental conditions.

WARNING: Always disconnect the device from the mains before cleaning.

WARNING: Always pull out the power plug =before any cleaning/disinfection procedures! Make sure the light source system has cooled down.

WARNING: The device must be set up in such a way that the power cord can be simply pulled out at any time.

WARNING: The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

WARNING: The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind.

WARNING: Applied parts of devices used in combination with this device must be of type BF or CF.

WARNING: Only KARL STORZ approved devices/instruments/accessories may be used in combination with the device.

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WARNING: Patient leakage currents may be additive when energized endoscopes are used with energized endotherapy devices. This is particularly important if a CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

WARNING: The compatibility of this device with any accessories and/or energized endotherapy devices should be checked in the individual instructions for use.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Always use the lowest possible settings for your intervention.

WARNING: The light source transmits high-energy light which can lead to high surface temperatures at the light cable and at the tip of the endoscope. Tissue damage is possible.

WARNING: Danger of glare! Never look into the light outlet when the light system is switched on.

WARNING: Insufficient ventilation of the unit can cause an internal build-up of heat which leads to a safety shut-off of the light source.

CAUTION: If a continuous tone is emitted when the unit is switched on, an error has occurred during start-up. Switch the unit off and start it up again. If the error occurs again, there is a general error. Send the unit in for repair.

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WARNING: Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

6.2 Normal use

It is only allowed to use the LED nova 150 light source in combination with approved accessories from KARL STORZ or with equipment that is approved as safety-related uncritical.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

WARNING: No modification of this equipment is allowed.

6.2.1 Indication

The device is intended to create light for endoscopic diagnostic observation and surgical procedures.

6.2.2 Contraindication

The medical device is contraindicated for ophthalmologic interventions. Mandatory protection filters for these kinds of interventions are not implemented for the medical device.

For other interventions the responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

6.3 Intended use

Cold light sources are designed to supply light for endoscopic diagnostic and surgical procedures.

6.4 User qualification

Only qualified doctors and medical assistant staffs are allowed to use the LED nova 150 light source after training. The training has to be conducted by skilled persons only. The area agency or the responsible KARL STORZ sales representative is available for trainings, for advices about further training alternatives.

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6.5 Patient profile

The use of the LED nova 150 is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

6.6 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been thoroughly trained at least once in the use of the device
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

6.7 Intended conditions of use

6.7.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operates properly before using it. During treatment using the LED nova 150 cold light fountain, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

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6.8 Safety precautions at the site of installation

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones.

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This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants. The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

! WARNING: The device should be placed on a flat surface or a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

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
6.9 Functions and features

Following installation and switching on of the unit, the light source emits cold light from its light port.

The brightness of the light is adjusted with the brightness regulator ③ on the front panel.

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
7. Installation and operating instructions

7.1 Unpacking the equipment

Carefully unpack the cold light fountain LED nova 150 and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

7.2 Basic equipment

KARL STORZ Set No.: 20161201
1 LED nova 150 cold light fountain 20161220
1 400 A power cord, length 300 cm
1 Instruction manual 96206245 D
1 Protective wallet for instruction manual

7.3 Installation

i NOTE: The cold light fountain LED nova 150 may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

! WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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Set the device on a horizontal, flat surface (Fig. 1).

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7.3.1 Connecting the ground line

The device is equipped with a connector ⑤ for attaching a ground line.
The device's ground line should be installed by a qualified electrician (Fig. 2).

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7.3.2 Connecting the power cord

Connect power cord and insert the power plug fully into the power cord receptacle ⑦ (Fig. 3).

WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

WARNING: Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

i NOTE: Ensure that access to the power plug is ensured at the site of installation.

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7.3.3 Connecting light cable

Plug the light cable into the socket ② as far as it will go (Fig. 4). Only handle the light cable using the handle piece. Never pull on the cable.

i NOTE: The light outlet ② is provided with an antiglare flap, which does not allow any direct outlet of light. Light is only emitted when the light cable is inserted.

! WARNING: Risk of burning! The light source transmits high energy light which can lead to high surface temperatures. Tissue damage is possible.

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7.3.4 Connecting the light cable to the endoscope or videoscope

Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base (Fig. 5 or 6).

i NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

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General applications

The LED nova 150 is suitable for use in a wide range of applications. The correct light intensity depends on the particular application, the instrument being used, and the connected light cable.

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**7.4 Commissioning the
LED nova 150 light source**

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7.4.1 Adjusting the brightness

The brightness can be adjusted with the control knob ③ (Fig. 8).

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7.4.2 Decommissioning

For decommissioning, switch off the unit at the power switch ① (Fig. 7).

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7.5 Safety functions of the equipment

- The unit includes a range of safety functions. During power-up, the unit performs a self-test. This enables defects or impermissible operating conditions to be detected right from the start and reported once via an acoustic signal.
- If impermissible operating conditions first become apparent during operation, this will also be reported. If the situation deteriorates, the light will switch off automatically. The acoustic signal gives the user the possibility to react to the development, either by resolving the cause of the error or by concluding their treatment in a controlled manner.

Warning signal

The warning signal is a tone lasting 1/2 a second, which is repeated several times.

- During operation, this indicates overheating or a device defect. Light continues to be output until the critical temperature has been reached.

CAUTION: Remedy the cause of overheating (e.g., poor ventilation, effect of heat or light, etc.) or conclude your treatment.

- If the temperature continues to increase, the cold light turns off when the unit reaches the critical temperature. The cold light remains turned off until the device is restarted (restart block), provided that the temperature has dropped below the critical value.
- If the critical temperature is not exceeded or if the unit temperature drops as a result of the measures taken, the unit continues to operate.

i NOTE: If the overheating fuse is activated, turn the device off and allow it to cool down sufficiently.

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Switch off the device and remove the power plug from the electrical outlet (Fig. 9).

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10

Remove the fuse holder ⑥ with a slim screwdriver or other suitable tool (Fig. 10).

WARNING: Only use fuses of the correct rating.

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11

WARNING: The fuse must fulfill the requirements for a component of high integrity characteristics.

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12

Refit the fuse holder ⑥. Reconnect the power cord (Fig. 12). Test the device for proper operation.

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8.2 Reprocessing

8.2.1 General warnings for medical devices

WARNING: Risk of infection: Incorrectly prepared medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: Country-specific laws and regulations must be observed.

i NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

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CAUTION: Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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
8.3 **Maintenance and safety check**

WARNING: Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

8.3.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.3.2 **Safety checks/repeat inspections as per IEC 62353**

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 7.4 and 7.5).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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8.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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8.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

WARNING: At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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8.6 Repair program

For videoscopes and devices, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany you can refer repairs directly to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 NEUHAUSEN

Servicehotline: +49 7461/708 980
E-Mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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8.8 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty. The same is valid for any cases of improper use (i.e. operation outside of the specified operating data).

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- Complete failure of the device.

Possible causes:

- Power supply failure.
- Defective power fuse.
- Power cord connector is not properly connected to device socket.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord connector firmly into device socket.

Symptom:

No light is emitted.

Possible causes:

1. Electronics are defective.
2. Overheating of the LED nova 150 (LED has switched off due to overheating of the unit caused by covered air vent).

Remedy:

1. Notify Service.
2. Remove objects from in front of the air vents. Switch off the LED nova 150, let it cool down (approx. 10 – 15 minutes) and switch it back on. Ensure that it is sufficiently ventilated.

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Symptom:

- No light output, power switch is lit (on).

Possible causes:

- Overheating due to covered air vent.
- Strong impacts or vibrations during operation.
- Power supply or LED is defective.

Remedy:

- Leave the light source to cool down, before switching the device on again.
 - Ensure that there is sufficient ventilation.
- Turn the unit off and then on again.
- Have the unit repaired.

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- Too little light.

Possible causes:

- Ends of light cable and/or endoscope are soiled.
- Light cable or endoscope is defective.
- LED service life has been exceeded.

Remedy:

- Clean ends of light cables, and light entrances and exits on the endoscope.
- Replace light cable or endoscope.
- Have the unit repaired.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

According to IEC 60601-1-2:
Observe the information on electromagnetic compatibility in chapter 11 (p. 30-46).

Directive compliance (for 20161220)
According to the Medical°Device Directive (MDD):
This medical device belongs to Class I
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

Conformidad con la norma (para 20161220)
Según CEI 60601-1, IEC 60601-1-6, ANSI/AAMI ES60601-1, IEC 60601-2-18, UL 60601-1, CAN/CSA-C22.2 No. 60601-1:

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9.3 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Item	Cat. no.
All KARL STORZ light cables are suitable for use	(see catalog; corresponding specialities)

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11. Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The cold light fountain LED nova 150 model 20 161220 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The cold light fountain LED nova 150 model 20 161220 is a Group 1 device (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

i NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic Compatibility (EMC) Information

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the equipment and the system should be observed to verify that they are operating normally.

WARNING: Portable HF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the 20 161220, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING: The use of portable and mobile HF communications equipment may have an impact on this or other pieces of medical equipment.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

WARNING: The use of accessories and cables other than those specified in the KARL °STORZ instruction manual may result in increased emissions or decreased immunity of the device. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

Guidance and manufacturer's declaration – electromagnetic emissions		
<p>The cold light fountain LED nova 150 model 20 161220 is intended for use in the electromagnetic environment specified below. The customer or user of the cold light fountain LED nova 150 model 20 161220 should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	The cold light fountain LED nova 150 model 20 161220 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The cold light fountain LED nova 150 model 20 161220 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. – Professional healthcare facility environment –
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	

Table 4 – ENCLOSURE PORT (according to IEC 60601-1-2)

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 9	
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	

Table 5 – Input a.c. power PORT (according to IEC 60601-1-2)

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
Electrical fast transients/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency	
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV	
Surges Line-to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Voltage dips	IEC 61000-4-11	0 % U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles Single phase: at 0° Grad	
Voltage interruptions	IEC 61000-4-11	0 % U_T ; 250/300 cycles	

Table 8 – Signal input/output parts PORT (according to IEC 60601-1-2)

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Electrical fast transients/bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency	
Surges Line-to-ground	IEC 61000-4-5	±2 kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment

Table 9 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (according to IEC 60601-1-2)						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the LED nova 150 model 20161220 may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skårholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610041, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang 1, 102dong, 8F
239 Jeongjial-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com