STORZ

KARL STORZ — ENDOSKOPE

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20 1614 20-1 Kaltlicht-Fontäne Power LED 175 SCB

INSTRUCTION MANUAL

20 1614 20-1 Cold Light Fountain Power LED 175 SCB

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20 1614 20-1 Koudlichtbron Power LED 175 SCB

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Important information for users of KARL STORZ devices

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It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the 'power LED 175' cold light source. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG. Utilization and reproduction by third parties require express authorization by KARL STORZ GmbH & Co. KG. All rights reserved.

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

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Controls, displays, connectors and their uses

- ① Power switch, lights up when in operation
- ② Light outlet point (not exchangeable)
- ③ Button to increase brightness *
- ④ Button to reduce brightness *
- ⑤ Current light intensity in percent
- ⑥ Standby button **
- ⑦ SCB indicator display ***
- ⑧ to ⑬ – see next page (symbols)
- ⑭ Potential equalization connector
- ⑮ Fuse holder
- ⑯ Power cord socket
- ⑰ Two SCB ports (SCB = KARL STORZ Communication Bus, see next page, top)

*

Keeping the buttons depressed engages continuous adjustment of the light intensity (autorepeat function).

**

- ⑥ Standby button **

** If you press this button once, the brightness is reset to the lowest value. The integrated LED illuminates.

When you press this button again, the brightness is reset to the last valid value. The integrated LED extinguishes.

- If there are several 'power LED 175' light sources available in the same SCB system (see next page), the system assigns each 'power LED 175' light source its own number (1-9). This number is displayed at position ⑦ and makes it clear to the user on which unit he is making the changes when using the remote control via the SCB.

- Example: The user changes the light intensity of LED light source no. 2 from '25%' to '50%' in the camera head's 'SCB menu'. This changes the light intensity of the 'power LED 175' for which a '2' is displayed at position ⑦.

- If only one 'power LED 175' light source is connected to the same SCB system, no number is displayed. Similarly, no number appears if there is no SCB system present.

If the light source is being remote controlled by means of a CCU, an 'A' is shown.

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The KARL STORZ-SCB interface (KARL STORZ Communication Bus) permits remote control of equipment functions, as well as remote display of equipment parameters.

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Read the instructions carefully before operating the equipment! (see ⑧ on page IV)

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Applied part type CF (Explanation: Limitation of leakage currents; CF = Cardiac Floating. See ⑬ on page IV).

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Potential equalization connector (⑭ on page IV)

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Electronic information product pollution control (China RoHS). The symbol is found on the identification plate on the underside of the unit.

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WEEE symbol: At the end of its service life, dispose of the device as electronic waste (see also p. 27). The symbol is found on the identification plate on the underside of the unit.

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Important information for users of KARL STORZ devices	III
Property rights	III
Photographs of the equipment.....	IV
Controls, displays, connectors and their uses	V
Symbols employed	VI
General information/Safety information	2
Device description	2
Warnings and cautions	2
Definitions.....	3
Intended use	8
User qualification	10
User profile	10
Patient profile	11
Affected body part or tissue type.....	11
Use	11
Function and features.....	12
Safety precautions at the site of installation	13
Safety precautions when operating the unit	13
Installation and operating instructions	14
Unpacking the equipment	14
Basic equipment	14
Installing and connecting the system	14
Creating SCB connections	16
Commissioning the SCB system	16
Operating the cold light source.....	17
Switching the unit on.....	19
Setting the brightness manually.....	20
Safety functions of the equipment	21
Maintenance	23
Fuse replacement.....	23
Reprocessing	24
Maintenance and safety check	26
Servicing and repair	27
Disposal	27
Repair program	28
Important information	28
Limitation of liability.....	29
Manufacturer's warranty	29
Technical description.....	30
Troubleshooting	30
Technical data	32
Standard compliance	33
Technical documentation.....	34
Spare parts, recommended accessories.....	35
List of spare parts	35
Accessories.....	35
Appendix	36
Electromagnetic Compatibility (EMC) Information	36

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Device description

The 'power LED 175' is designed to supply light for endoscopic diagnostic and surgical procedures.

The device can be controlled via buttons on the front panel or via SCB (SCB = KARL STORZ Communication Bus, see 17 on page IV, or the application example in the last three paragraphs on page V).

The main advantages of this LED light source are:

- Efficiency: High light intensity at low power consumption.
- Comfortable working environment: Intelligent cooling concept for minimal noise emission.

WARNING: The 'power LED 175' may be subject to sudden failures. The operator is bound to keep a spare cold light system available during endoscopic interventions which require the permanent availability of cold light.

- Economic: There are normally no costs for a replacement lamp and no expenses for exchanging the lamp.

Safety information

Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device.

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Warnings and cautions

Definitions

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A Note indicates special information about operating the unit.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

WARNING: Additional equipment which is connected to medical electrical equipment must be proven to comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (MDD: Appendix 1: 13.6.c, IEC 60601-1: 6.8.2.c, 19.2.b, 19.2.c).

WARNING: Do not use the device in case of visible damages.

WARNING: Test this equipment prior to and after each surgical procedure to ensure that it functions correctly.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Do not use in areas with flammable/explosive and/or oxidising gases.

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Safety instructions

Warnings and cautions

WARNING: Danger of glare! Never look directly into the light port or into the end of a light cable connected to the light port.

WARNING: Never look into the light beam using optical instruments (e.g. magnifying glasses).

WARNING: Do not store liquids on the device.

WARNING: If liquid enters the device, provide adequate time for evaporation.

WARNING: Danger of an electrical shock. The electrical installations in the room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Danger of an electrical shock:

- Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' mains supply.
- Routinely inspect power cord and outlet. Do not use the device if the inspection reveals damage.

WARNING: Danger of electric shock! Do not open the LED light source! Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee.

WARNING: Always unplug the unit before doing any maintenance work or performing any cleaning/disinfection procedures.

WARNING: If the light switches itself off during operation, an error might have occurred. Turn the device off and switch to your spare device. See details on page 21.

WARNING: The light port and the housing may be hot.

WARNING: Please observe the installation guidance in this instruction manual.

WARNING: Please observe cleaning/disinfection guidelines in this instruction manual.

WARNING: Ensure that air is able to circulate freely around the unit. Cooling inlets and outlets must not be covered.

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Safety instructions**Warnings and cautions**

WARNING: When changing the fuse, only use the fuses specified on the fuse holder.

WARNING: When using accessories, please refer to the respective instruction manuals.

WARNING: The device must be placed on a position which allows the operator to identify visual information signals.

WARNING: Do not use the unit if the displays on the front of the unit are not in perfect working order.

WARNING: Always supervise the start up sequence of the device in order to detect possible key or display errors.

WARNING: Do not use the unit if the entry buttons (operating buttons on the front or rear) are not in perfect working order.

WARNING: The KARL STORZ 'power LED 175' may be subject to sudden failures. The operator is bound to keep a spare cold light system available during endoscopic interventions which require permanent availability of cold light.

WARNING: The device may only be operated at the supply voltage stated on the identification plate or in the instruction manual.

WARNING: Do not use the device with any accessories or cables other than those specified in the instruction manual.

WARNING: If the light source stops and the display indicates overheating, allow sufficient time for the device to cool down before reactivating it.

WARNING: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix during installation and commissioning.

WARNING: The device does not provide a life supporting function.

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Safety instructions

Warnings and cautions

WARNING: The light intensity must be set to a value which is reasonable for the specific intervention.

WARNING: Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's wellbeing is maintained in the event of such interference.

WARNING: The housing may be subject to a maximum load of 20 kg.

WARNING: When endoscopes are used with energized endoscopically-used accessories, the patient leakage currents may be additive. This is particularly important if a Type CF endoscope is used, in which case a Type CF endoscopically-used accessory should be used in order to minimize the total patient leakage current.

WARNING: Risk of burns to the patient and user: The temperature of the endoscope's light exit window can exceed 41°C. Moreover, high amounts of radiant energy are emitted from the light exit window. This can, for example, lead to irreversible tissue damage and coagulations, regardless of whether there is direct tissue contact or not. The user must therefore ensure

- that sufficient distance is always maintained between the tissue and the endoscope's light exit window;
- that the light energy emitted from the device is reduced as much as is possible;
- that the tissue is exposed to the light energy for as short an amount of time as is possible.

CAUTION: Acclimatization time!

If the unit is stored/transported at temperatures less than +10 °C before being put into operation, condensation may build up on or in the unit. In this case, it is essential to allow for an acclimatization time of approx. 3 hours before use.

CAUTION: Turn the cold light source off each time it has been used.

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CAUTION: Only operate the unit with the voltage stated on the unit's identification plate.

CAUTION: When changing the unit's fuses, only use fuses with the specified parameters (see Technical data on page 32).

CAUTION: Do not allow the unit to overheat. It is irrelevant whether the cause of overheating is attributed to just the external circumstances (e.g., storage) or operation in improper circumstances. In cases of overheating the light source does not work and switches off.

CAUTION: If a continuous tone is emitted when the unit is switched on, an error has occurred during start-up. Switch the unit off and start it up again. If the error occurs again, there is a general error. Send the unit in for repair.

NOTE: If the unit is already too hot when it is switched on, it will not start and no light is emitted. If the overheating occurs during operation, optical and acoustic warnings are given and the unit will continue to work for a while. However, when the critical temperature is exceeded, the light output switches off. (Further information on possible errors can be found from page 30 onwards.)

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Light cable

WARNING: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely.

WARNING: Observe the instruction manuals, warnings, and cautions of all the products used in combination.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended Use

Cold light sources are designed to supply light for endoscopic diagnostic and surgical procedures.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The 'power LED 175' light source must only be operated under the following conditions:

- Use for medical and veterinary purposes.
- Used exclusively with accessories which are offered, recommended and designated as suitable by KARL STORZ (see page 35).
- Used exclusively in locations appropriate for medical applications (operating theaters) including inside the patient area.

WARNING: Operate device only under standard environmental conditions (see Technical data page 32).

- The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.
- Operate only in properly ventilated areas. Do not place the light source on its side, front, rear, or upside down (correct orientation is shown on page 15).
- No additional energy supply (apart from the electrical energy used to operate the device) is possible/required.

WARNING: It must not be used in the proximity of flammable liquids or gases or in areas of increased oxygen levels under any circumstances.

- It must not be used in areas with high dust levels under any circumstances.
- It must not be used in areas of increased air pressures (diving stations) or low air pressures (at elevations above 3000 m above sea level).
- It must not come into contact with aggressive liquids (acids, alkalis, solvents) or aggressive gases.
- It must not be operated outdoors or in areas where condensation forms (in humid environments).
- It must not be installed or operated in vehicles (on land, water or in the air) or in mobile units (containers). An exception is made here for mobile equipment trolleys.

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- It must not be installed in the proximity of radiators or other locations where the unit would be exposed to high levels of heat or light.
- It must not be installed or operated in areas with strong electromagnetic fields (electric motors, transformers, RF generators, X-ray units, MRI machines, etc.).
- It must not be installed or operated in areas of high radiation levels (nuclear diagnosis and treatment units).

WARNING: The housing may be subjected to a maximum load of 20 kg. Never climb onto the housing.

User qualification

WARNING: The 'power LED 175' may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the instrument. Such instruction should only be provided by persons equipped to do so on the basis of their skills.

Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and information on further training alternatives.

User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

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Patient profile

The use of 'power LED 175' is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

Affected body part or tissue type

The field of application is not confined to certain parts of the body or types of tissue.

Use

The device may only be used in hospitals and doctors' offices in the ambient conditions given in the technical data.

Other conditions of use:

- Frequency of use up to several times a day
- Length of use from a few minutes to several hours a day
- Position it on a vibration-free, fixed base
- Can be moved only by installing it on an equipment trolley
- Can be swiveled only by installing it on ceiling supply units

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Functions and features

Following installation and the switching on of the unit, the light source emits cold light from its light port. The intensity of the light can be regulated by pressing the + and – buttons ③ and ④. The current intensity is displayed by the green LED. When the unit is switched on, the lowest setting is always active.

i **NOTE:** The light intensity can also be controlled by remote control or automatically, provided the technical requirements for this are satisfied (see following paragraph).

i **NOTE:** Various light sources (up to 9) can be integrated in one SCB installation. The units can then only be differentiated on the remote station (monitor) by their SCB address. For further information, please see the last three bullet points (•) on page V.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The unit is equipped with a connector for attaching a ground line. It should be connected before power is applied to the unit in compliance with current national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

Also observe the safety instructions on pages 3-12.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit. During treatment using the cold light source 'power LED 175' the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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Unpacking the equipment

Carefully unpack the cold light source and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Basic equipment

- 1 Cold light source 'power LED 175'
- 1 Power cord 400 A
- 1 Instruction manual 96206197NL
- 1 Wallet for instruction manual
- 1 SCB connection cable 100 cm 20090170

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Installing and connecting the system

- i** **NOTE:** The cold light source 'power LED 175' may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.
- !** **WARNING:** This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Installation and operating instructions

CAUTION: Only operate the unit on a flat, horizontal surface.

CAUTION: Danger of overheating. It is important to ensure free circulation of air, above, behind, and around the sides (particularly the **right** side) (rear distance to wall min. 10 cm).

Place the unit on its **horizontal underside** onto a base with the following characteristics:

- sufficiently large and stable;
- flat and horizontal;
- fixed, non-slip surface.

This can also be the top of another KARL STORZ unit.

Position the unit so that you can operate it and see the displays from the place where you are working.

NOTE: For the safety of the application it is important that the unit can be easily operated from the normal working position and the displays are clearly visible. The set-up position must be selected accordingly.

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The unit is equipped with a connector for attaching a ground line (14). The unit's ground line should be installed by a qualified electrician.

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Creating SCB connections

- i** **NOTE:** The term ‘SCB’ is explained on page VI. The bottom of page V also provides an example of use of the SCB function (see last three bullet points (•)).

If you wish to integrate the ‘power LED 175’ into an existing SCB system, connect the SCB cable 200901 70 up to the SCB port (17) on the rear of the unit.

- i** **NOTE:** The SCB ‘power LED 175’ system can be operated with any camera control unit (CCU) which has an SCB.

On the rear of every SCB unit there are 2 SCB ports (17). Use the supplied SCB cable 200901 70 to connect the ‘power LED 175’ with the camera control unit (e.g., IMAGE 1 HUB™ HD) – there is no difference between the two SCB ports.

- i** **NOTE:** Ensure that the no SCB cable is only inserted at one end as this can cause transmission problems within the system. In cases of malfunction, please consult the Service department.

- i** **NOTE:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed. Please read the instruction manual for the KARL STORZ-SCB system carefully and follow the instructions exactly.

Commissioning the SCB system

The units of the SCB system can be turned on in any order. When the camera control unit (e.g., IMAGE 1 HUB™ HD) is turned on, the following SCB message appears: New Device

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**Operating the
cold light source**

CAUTION: Use only with electrical supply lines that agree with the line voltage stated on the unit's identification plate.!

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Insert the power cord into the unit's receptacle ⑯, pushing it firmly into place.

WARNING: Danger of explosion and fire hazard: Only insert the power cord into and remove it from electrical outlets located outside areas subject to explosion hazards.

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Plug the light cable into the socket as far as it will go. Our photo shows a light cable for an endoscope. Only handle the light cable using the handle piece. Never pull on the cable.

ⓘ NOTE: The light outlet ② is provided with an antiglare flap, which does not allow any direct outlet of light. Light is only emitted when the light cable is inserted.

⚠ WARNING: Risk of burns: The light outlet may heat up to over 41°C.

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Connecting the light cable to the endoscope or videoscope


- Photo 1: Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base.
- Photo 2: Connect the light cable to the videoscope (quarter turn of the knurled screw on the screw base).

ⓘ NOTE: The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

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
Switching the unit on

Switch the unit on with switch ①.

i **NOTE:** If no error is noted, the unit's self test is ended with a short tone. The light source is then at the lowest intensity setting (approx. 5%). However, the intensity display ⑤ shows the last setting used before the unit was switched off. That means: the segment which was the highest illuminated segment prior to the unit being switched remains permanently lit.

Choosing the last setting

The unit starts in standby mode (at the lowest setting). If the standby button is then pressed, the last valid setting is activated.

Choosing the lowest setting

The unit starts in standby mode and is already automatically at the lowest setting. Therefore, no further action is necessary.

Choosing the individual setting

The unit starts in standby mode and is automatically at the lowest setting. To select the desired setting, press the + button (repeatedly, if necessary).

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Setting the brightness manually

The brightness can be set with the buttons ③ and ④.
To increase the brightness press the button ③.
To decrease the brightness press the button ④.
The current brightness setting is shown by the LED display ⑥.

i NOTE: If the buttons ③ and ④ are held down, the autorepeat function is activated. The value then continues to change in the selected direction until the minimum or maximum setting is reached.

Standby button ⑥

Pressing this button once resets the brightness to the lowest value. The integrated LED illuminates. Pressing this button again resets the brightness to the last valid value. The integrated LED extinguishes.
If the plus button is pressed when the device is in standby mode, the light intensity returns to the lowest value (5 %).

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Safety functions of the equipment

- The unit includes a range of safety functions. During power-up, the unit performs a self-test. This enables errors or impermissible operating conditions to be detected, right from the start. In the case of errors, the unit enters safety mode immediately. No light is emitted.
- If impermissible operating conditions first become apparent during operation, this will initially be signaled. If the situation deteriorates, the light will switch off automatically. This gives the user the possibility to react to the development either by resolving the cause of the error or by concluding his treatment in a controlled manner.

The following error displays may be shown:

a) Constant tone after turning on and dark intensity display:

There are two possible causes of this:

- A hardware problem was discovered. The light source is not working properly. This can also be assumed if the error is signaled constantly even when the unit is cold. The unit must be sent in for repair.
- The unit is too hot:
It is irrelevant whether the unit has become too hot during operation or solely as a result of the external circumstances (storage). As long as it is too hot, it will not start.
The unit must be switched off and allowed to cool.

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**b) Warning sound (two short tones) every
15 seconds; simultaneously, the fields of the
display which show the current light intensity
flash**

• This warning can only occur during operation. The same is true of the cause. The unit has become too hot during operation. If the signal sounds, the unit can still be used for a while.

CAUTION: Remove the cause of the overheating (e.g., poor convection, effect of heat or light, etc.) or reduce the light intensity (output). If this is not possible, conclude your treatment.

- If the temperature continues to increase, the cold light turns off when the unit reaches the critical temperature. The optical/acoustic signal remains on. It turns off once the temperature drops back down to an admissible value. However, the cold light remains turned off (restart block).
- If the critical temperature is not exceeded or if the unit temperature drops as a result of the measures taken, the unit continues to operate. As soon as the temperature reaches admissible values, the signal (optical/acoustic) is deleted.

Notes on errors a) and b)

NOTE: If the overheating fuse has responded and the unit is subsequently turned off and immediately turned back on, error a) is signaled (constant alarm). The unit must firstly be allowed to cool down to permissible temperatures.

- As soon as the unit has cooled down, turn the unit back on. As long as there are no other problems (e.g., blocked fan), the unit starts again normally.

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
Fuse replacement

Switch off the unit and remove the power cord from the electrical outlet.

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
Release the line fuse holder (15) using a slim screwdriver (tip approx. 6 mm wide) or another suitable tool.

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Insert new fuses.

! **WARNING:** Use only fuses with the correct ratings.

	20161420-1
	100 – 240 VAC
Power fuse	2 x T1.6AH250V

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Refit the power fuse holder (15).
Reconnect the power cord.
Test the instrument for proper operation.

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Reprocessing

General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **CAUTION:** National laws and regulations must be observed.

 **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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WARNING: Always pull out the power plug before any cleaning/disinfection procedures! Make sure the light source system has cooled down.

CAUTION: Avoid allowing fluids to enter the unit.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

The chemical manufacturer's specifications concerning material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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
Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety checks/repeat inspections as per IEC 62353

WARNING: *Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.*

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section "Operation").

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty. The same is valid for any cases of improper use (i.e. operation outside of the specified operating data).

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Troubleshooting

WARNING: Always unplug the unit before
all maintenance work!

Symptom:

Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction
manual. Make sure you use the correct type of
fuse.

Symptom:

No light output, power switch is lit (on).

Possible causes:

- Strong impacts or vibrations during operation.
- Power supply or LED is defective.

Remedy:

- Turn the unit off and then on again.
- Have the unit repaired.

Symptom:

Too little light.

Possible causes:

- Ends of light cable and/or endoscope are
soiled.
- Light cable or endoscope is defective.
- LED service life has been exceeded.

Remedy:

- Clean ends of light cables, and light entrances
and exits on the endoscope.
- Replace light cable or endoscope.
- Have the unit repaired.

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Symptom:

No light output, power switch is lit (on).
LED-brightness indicator is flashing.

Possible causes:

- Defective LED.
- Defective fan.
- Defective temperature sensor.
- Overheating of the unit.

Remedy:

- Switch off the unit and allow to cool down.
- Have the unit repaired.

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- Led defect.
- Ventilator defect.
- Temperatuursensor defekt.
- Apparaat oververhit.

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‘power LED 175’	20 1614 20-1
Power supply voltage	100-240 VAC (± 10 %)
Power frequency	50/60 Hz
Power consumption	1.1 A
Power fuses	2 x T1.6AH250V
Operating temperature	10 °C...40 °C
Dimensions	305 mm x 110 mm x 233 mm (w x h x d)
Weight	4.6 kg
Storage/transport conditions:	
Air humidity (RH, non-condensing)	5 %...95 %
Temperature	-10 °C...60 °C
Abs. air pressure	+700 hPa...+1060 hPa (= -600 ... +3000 m above sea level)

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- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type CF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility
Information in the appendix (p. 36-47).

**Directive compliance
(for 20 1614 20-1)**

This medical device belongs to Class I

This medical device bears the CE mark in
accordance with the Medical Device Directive
(MDD) 93/42/EEC. A code number after the CE
mark indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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List of spare parts

Item	Cat. no.
Power fuses T1.6 AH250V, pack of 10	1973290
Power cord (grounded)	400 A
'Hospital Grade' Power cord (USA)	400 B
Instruction manual	96206197NL

Accessories

All KARL STORZ light cables currently available on
the market are suitable for use.

For detailed information please contact:
KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

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WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix during installation and commissioning.

The 'power LED 175' model 20 1614 20-1 corresponds to IEC 60601-1-2:2007 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The 'power LED 175' model 20 1614 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The 'power LED 175' model 20161420-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

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Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The 'power LED 175' model 20161420-1 is intended for use in the electromagnetic environment specified below. The customer or user of the 'power LED 175' model 20161420-1 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	The 'power LED 175' model 20161420-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The 'power LED 175' model 20161420-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	

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Table 2			
Guidance and manufacturer's declaration – electromagnetic immunity			
<p>The 'power LED 175' model 201614 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the 'power LED 175' model 201614 20-1 should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Complies <5 % U_T * (>95 % dip in U_T) for 0.5 cycle Complies 40 % U_T (60 % dip in U_T) for 5 cycles Complies 70 % U_T (30 % dip in U_T) for 25 cycles Complies <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The 'power LED 175' model 201614 20-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the 'power LED 175' model 201614 20-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF acc. to IEC 61000-4-6</p> <p>Radiated RF acc. to IEC 61000-4-3</p>	<p>3 V_{rms} 150 kHz to < 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V_{rms}</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the 'power LED 175' model 201614 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = [3.5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 6			
Recommended separation distances between portable and mobile RF communications equipment and the 'power LED 175' model 201614 20-1			
The 'power LED 175' model 201614 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseam.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6. Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0) 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, Kumho Asiana Plaza
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang Intenational Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 028, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.jp

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen
Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com