STORZ

KARL STORZ — ENDOSKOPE

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IMAGE1 S H3 Kameraköpfe

TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P, TH 104 H3-ZA, TH 106 H3-M COVIEW®

INSTRUCTION MANUAL

IMAGE1 S H3 camera heads

TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P, TH 104 H3-ZA, TH 106 H3-M COVIEW®

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IMAGE1 S H3-camerakoppen

TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P, TH 104 H3-ZA, TH 106 H3-M COVIEW®

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**1 Important note for users of
KARL STORZ devices and
instruments**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ name. Like all of our other
products, this product is the result of years of
experience and great care in manufacture. You
and your organization have decided in favor of
a modern, high quality piece of equipment from
KARL STORZ.

1.1 Intended use

The camera head is used together with the
camera control unit and is connected to a
telescope during endoscopic or microscopic
interventions.

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2 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

2.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

This instruction manual is intended to serve only as a supplement to the KARL STORZ video systems instruction manual, which was enclosed with your camera control unit. To ensure the safe and efficient use of this KARL STORZ product, it is essential that you read the entire video system instruction manual before connecting the camera system.

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2.2 General safety information

WARNING: Read this instruction manual thoroughly before using the equipment.

WARNING: Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.

WARNING: Perceived electric shocks may be experienced when using the camera system in conjunction with electro-surgical units.

WARNING: Electrostatic discharge (ESD) may be perceived as an electric shock when touching the camera system for the first time. In the event of image loss due to ESD, switch the camera control unit off and back on again to restore the image.

WARNING: When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.

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WARNING: All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC/EN 60601-2-18, the standard for medical electrical equipment – particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery. HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

WARNING: Risk of blindness. Never look into the open end of a light cable.

WARNING: Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.

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



- ① Instrument coupling, C-MOUNT connection with H3-M
- ② Focusing ring/knob
- ③ Zoom ring
- ④ User control button: for scrolling (downwards) through the menu or activating a camera function
- ⑤ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑥ User control button: for accessing the menu/ making a selection

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 Number of products in the product packaging

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5 Operating instructions

5.1 Frequency of the High Definition camera head (50/60 Hz)

The IMAGE1 S H3 camera heads are suitable for use with frequencies of both 50 Hz and 60 Hz. The frequency can be adapted to suit local conditions.

i NOTE: The method for switching frequency is described in the instruction manual IMAGE1 S – TC 200, TC 300, TC 301, TC 302 (96206286NL).

5.2 Using endoscopy camera heads (TH 100, TH 101, TH 103, TH 104)

5.2.1 Connecting an endoscope and light cable

IMAGE1 S H3 camera heads have an integrated instrument coupler and telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the instrument coupler clockwise and insert the eyepiece of the endoscope into the coupler. When connecting the endoscope to a Pendulum camera head, the Pendulum lock must click into place before the outer ring of the grasping mechanism is turned. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the eyepiece of the endoscope by tightening the knurled screw on the light cable through a quarter turn.

! **CAUTION:** Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

i NOTE: For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (96206342ZS), which can be downloaded or requested at www.karlstorz.com. KSEA uses a different fiberoptic light cable manual. The part number for KSEA is PI-000011.

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5.2.2 Adjusting the focus

The focus can be adjusted by turning the focusing ring/knob on the camera lens.

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5.2.3 Adjusting the image size with zoom camera heads

The image size can be adjusted by turning the zoom ring.

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5.2.4 Operating the user control buttons

- ④ Scrolling downwards: Scroll downwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑤ Scrolling upwards: Scroll upwards through menus, options, etc. Can also be assigned a programmable function/feature.
- ⑥ Menu/select: Access the menu. Once the menu is displayed, this button can be used to select the highlighted option (identical to Enter on a computer keyboard). Hold this button down for seven seconds to call up the 'Options' menu.

For detailed instructions on the programming and use of the camera head buttons, please refer to the KARL STORZ IMAGE1 S (Art. No.: 96206286NL) and KARL STORZ IMAGE1 (Art. No.: 96206114D) instruction manuals.

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5.3 Using the microscope camera head TH 106 IMAGE1 S H3-M COVIEW®

The IMAGE1 S H3-M COVIEW® microscope camera head is connected to a suitable interface on the surgical microscope using a TV adaptor. Suitable interfaces are optical beam splitters and co-observer connectors. These options are available for all surgical microscopes.

Some surgical microscopes are equipped with a C-MOUNT connection that is directly coupled with the IMAGE1 S H3-M COVIEW® camera head.

5.3.1 Detaching the camera head cable

The IMAGE1 S H3-M COVIEW® camera head is delivered with a 9 m cable as standard. Detaching the camera cable from the camera head is recommended in order to make installation easier. To detach the camera cable, unscrew the cable connector from the camera head by turning in the counterclockwise direction and carefully pull out the cable.

5.3.2 Connecting the camera head to the TV adaptor and surgical microscope

To connect the the IMAGE1 S H3-M COVIEW® camera head to the surgical microscope, a TV adaptor is required (e.g. KARL STORZ QUINTUS® TV adaptor). Screw the camera head carefully onto the C-MOUNT connection of the TV adaptor and then secure the camera head and TV adaptor to the co-observer interface of the surgical microscope. For more detailed installation instructions, please refer to the instruction leaflet provided with the QUINTUS® TV adaptor.

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5.3.3 Connecting the camera cable

To connect the camera head cable to the camera head, align the tip of the cable plug carefully with the recess provided and press the cable carefully onto the camera head. Screw the cable connector to the camera head by turning in the clockwise direction.

CAUTION: To avoid damage, the camera cable should not be connected to the Camera Control Unit (CCU) during installation!

CAUTION: The camera cable should NOT be connected to/disconnected from the camera head when the Camera Control Unit (CCU) is switched on.

5.3.4 Operating the user control buttons

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For detailed instructions on the programming and operation of the camera head buttons, please refer to the KARL STORZ IMAGE1 instruction manual.

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6.2 Storage and Operating conditions

Refer to the compatible camera control unit (CCU) manual for relevant storage and operating conditions. The CCU and camera head were evaluated as a system for IEC 60601-1.

6.3 Standard compliance

According to IEC 60601-1; IEC 60601-2-18; UL 60601-1:

- Protection of applied part against electric shock: Type CF 

6.4 Directive compliance

According to the Medical Device Directive (MDD):

This medical device belongs to Class I

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

6.5 Information on electromagnetic compatibility (EMC)

Refer to the Electromagnetic Compatibility (EMC) information Appendix of the compatible camera control unit (CCU) manual for relevant electromagnetic compatibility (EMC) information.

The EMC warning statements, caution, notes and Emission/Immunity limits specified in the EMC Tables of the compatible CCU apply to the camera head. The CCU and camera head were tested as a system.

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7 Reprocessing

The following chapter provides information on the care and cleaning of the IMAGE1 S H3 camera heads.

7.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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7.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile upon delivery. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, national regulations regarding personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ.

CAUTION: National laws and regulations must be observed.

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i **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

! **CAUTION:** Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be processed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

! **CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.

! **CAUTION:** The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.

! **CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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7.3 Preparation for cleaning and disinfection

Immerse the camera head in a suitable container containing a pH-neutral (pH value of 6.0 – 8.0) enzymatic cleaning solution as per the manufacturer's instructions immediately after use to prevent blood, protein and other contaminants from drying onto the equipment.

CAUTION: Inspect the camera head cable for breaks and cuts. Camera heads with damaged cables must not be immersed. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

Heavy soiling must be removed from the camera head immediately after use. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.3.1 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixing effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material incompatibility must be observed.

CAUTION: The TH 106 IMAGE1 S H3-M COVIEW® microscope camera head may only be cleaned by wipe-down disinfection. Other reprocessing or sterilization methods are not permitted.

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7.3.2 Manual cleaning

The camera head and camera head cable must be completely immersed, preferably in a pH-neutral enzymatic cleaning solution as per the manufacturer's instructions and microbiologically pure/completely demineralized water. It must be ensured that the surfaces are covered and that no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with a soft brush (Art. No. 27652) or sponge. Finally, it must be rinsed with microbiologically pure, sterile or completely demineralized water to ensure neutralization. Thorough rinsing of the camera head is necessary for removing any soiling and cleaning agent residues which could interfere with disinfection.

7.3.3 Manual disinfection

The camera head and the camera head cable must be completely immersed in a disinfectant solution. It must be ensured that the surfaces are covered and that no air bubbles are present. Following exposure, rinse the camera head three times with microbiologically pure water, sterile water, or completely demineralized water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Thorough rinsing of the camera head is necessary for removing any residual disinfectant that could interfere with sterilization.

All the surfaces are then dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose. Wipe the camera head connector with 70% isopropyl alcohol to remove any residual disinfectant.

Do not allow exposed glass windows to air dry. 70% isopropyl alcohol may be applied to glass surfaces with a cotton tip applicator to prevent streaks and spots. Dry the surfaces thoroughly with a cotton tip applicator after applying the alcohol.

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7.4 Inspection and care

The cleaned camera head and cleaned camera head cable must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the camera head and the camera head cable must be manually cleaned and subjected to a full cleaning procedure once more.
- Camera heads with damaged cables must be withdrawn from use.

CAUTION: Inspect the camera head cable for breaks and cuts. Camera heads with damaged cables must not be sterilized or disinfected. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

- Afterwards, a functional check must be carried out.

Before sterilization and/or disinfection, coil the camera head cable into loops with a diameter of at least 15 cm (6"). In doing so, do not kink or twist the cable.

7.5 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953). National regulations must be observed.

7.6 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S H3 camera heads:

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7.6.1 Steam sterilization in the fractionated prevacuum procedure only suitable for autoclavable camera heads

CAUTION: Of the camera heads compatible with IMAGE 1 HUB™ HD and IMAGE1 S only, IMAGE1 S H3-ZA (TH 104) is autoclavable.

NOTE: Place the camera head in a suitable sterilization container (e.g. 39301 Z3TS).

The autoclavable camera head and the camera head cable must be sterilized using a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 °C – 137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

It is recommended to leave the camera heads to dry for at least 4 minutes, so that the card-edge connector is completely dry before its insertion.* After steam sterilization, the camera heads should not be touched until adequately cooled.

CAUTION: Do not immerse the camera head in cold liquid after removing it from the steam sterilizer. Due to the ‘temperature shock’, this may result in damage to the camera lens and invalidation of the warranty.

7.6.2 Steam sterilization using the gravitation method only suitable for autoclavable camera heads

CAUTION: Of the camera heads compatible with IMAGE 1 HUB™ HD and IMAGE1 S, only IMAGE1 S H3-ZA (TH 104) is autoclavable.

NOTE: Place the camera head in a suitable sterilization container (e.g. 39301 Z3TS).

The autoclavable camera head and the camera head cable must be sterilized using a gravitation procedure at 121 °C – 122 °C with an exposure time of 30 minutes.

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CAUTION: Do not immerse the camera head in cold liquid after removing it from the steam sterilizer. Due to the 'temperature shock', this may result in damage to the camera lens and invalidation of the warranty.

* When using steam sterilizers without an integrated drying cycle, the card-edge connector can be dried off with a lint-free, sterile cloth once it has cooled down. Before insertion, all excess moisture must be removed from the card-edge connector.

7.6.3 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

NOTE: To ensure that the medical device concerned can be sterilized using the different STERRAD® units, refer to the 'STERRAD® Sterility Guide' (www.sterradsterilityguide.com).

CAUTION: STERRAD® sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

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The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S H3 camera heads:

- STERRAD® 100S
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
‘short’ cycle
- TH 103 H3-P
‘long’ cycle
- STERRAD® NX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
‘standard’ cycle
- STERRAD® 100NX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 S H3-ZA
‘standard’ cycle
‘DUO’ cycle

7.6.4 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO®

More detailed information regarding the cycle used for various device generations can be obtained from the manufacturer STERIS®.

CAUTION: V-PRO™ sterilization may cause cosmetic changes to the camera head that do not usually impact the camera’s functionality.

The following STERIS® AMSCO® sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S H3 camera heads:

- V-PRO™ 1
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
Equivalent to ‘lumen’ Cycle
- V-PRO™ 1 Plus
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
‘lumen’ cycle
- V-PRO™ maX
TH 100 H3-Z, TH 101 H3-ZI, TH 103 H3-P,
TH 104 H3-ZA
‘lumen’ cycle

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7.6.5 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30-45 minutes.

 **NOTE:** Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

 **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated for the materials by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing procedure of the EO sterilizer (concentration, process control).

7.6.6 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1 (except U. S.)

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS®.

7.6.7 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E (only valid in U. S.)

More detailed information regarding sterilization parameters can be obtained from the manufacturer STERIS®.

7.6.8 High-Level Disinfection

 **WARNING:** High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

 **CAUTION:** Any deviations from the recommended parameters for disinfection must be validated by the user.

 **CAUTION:** Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

 **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life.

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The following high-level disinfection methods have been validated and approved by KARL STORZ for the IMAGE1 S H3 camera heads:

7.6.9 CIDEX (2.4% glutaraldehyde solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4% glutaraldehyde solution for a minimum of 45 minutes at 25 °C (77°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

CAUTION: Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head and the card-edge connector with sterile water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

7.6.10 Resert XL HLD (2.0% hydrogen peroxide solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

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To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head with sterile water is very important for removing any residual disinfectant.

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7.7 Reprocessing limit

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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8 Maintenance**8.1 Maintenance and safety check****8.1.1 Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

8.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.

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Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see Chapter 7.3.12 of the IMAGE1 S CCU instruction manual no. 96206286 NL).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

8.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

8.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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9 Repairs and warranty

9.1 Repair program

CAUTION: Opening the camera head will result in damage and void the warranty.

To bridge the repair period, you will, as a rule, receive a unit on loan which must then be returned to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel).

We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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9.2 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

9.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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i NOTE: The S technologies are available with IMAGE1 S.

* IMAGE1 camera control units are available with various configurations and image outputs. These camera control units are compatible with IMAGE1 S camera heads.

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**11 Spare parts,
recommended accessories**

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Extension cable 5.2 m	22 2200 72
Plastic Container for Sterilization and Storage of camera heads IMAGE1 S H3-Z (TH 100) and IMAGE1 S H3-ZA (TH 104), autoclavable	39301 Z3TS
H3-M COVIEW® Camera Head Cable 9 m	22 2200 71
Instruction Manual	96206511NL
Plastic Container for Sterilization and Storage of camera heads IMAGE1 S H3-ZI (TH 101) and IMAGE1 S H3-P (TH 103)	39301 PHTS

* For safe and comfortable working with KARL STORZ cameras we recommend the use of sterile drapes.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: ausra.kukeniene@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlyoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: + 386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6. Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0) 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, Kumho Asiana Plaza
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang Intenational Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 028, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.jp

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

 **KARL STORZ GmbH & Co. KG**

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com