STORZ

KARL STORZ—ENDOSKOPE

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INSTRUCTION MANUAL
IMAGE1 S™ 4U camera head TH 120 4U

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Important note for users of KARL STORZ devices and instruments

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1 Important note for users of KARL STORZ devices and instruments

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual, in combination with the IMAGE1 S™ instruction manual (96206286D), is intended to aid you in the proper installation, connection, and operation of this device. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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1 Important note for users of KARL STORZ devices and instruments	IV
2 Safety instructions	4
2.1 Explanation of warnings and cautions	4
2.2 General safety information	5
3 Intended Use Camera Heads	7
3.1 Contraindications	7
3.2 User qualification	7
4 Images of the equipment	8
5 Symbols employed	9
5.1 Symbols for operation	9
5.2 Label and packaging symbols	9
6 Operating instructions for the IMAGE1 S™ 4U camera head	10
6.1 Frequency of the IMAGE1 S™ 4U camera head (50/60 Hz)	10
6.2 Operating the IMAGE1 S™ 4U camera head	10
6.2.1 Connecting an endoscope and light cable	10
6.2.2 Adjusting the focus	11
7 Technical description	12
7.1 Specifications	12
7.2 Storage and Operating conditions	12
7.3 Standard compliance	12
7.4 Directive compliance	12
8 Reprocessing (not validated for US customers)	13
8.1 Care and handling	13
8.2 Reprocessing	14
8.2.1 Reference chart	14
8.2.2 General warnings	15
8.3 Accessories	17
8.4 Preparation for cleaning and disinfection	17
8.4.1 Manual wipe-down disinfection	17
8.4.2 Manual precleaning	18
8.4.2.1 Brushing the surfaces	18
8.4.2.2 Auxiliary pretreatment with ultrasound	18
8.4.3 Manual cleaning	18
8.4.4 Manual disinfection	18

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8.5	Assembly, inspection and care	19
8.6	Packaging systems	20
8.7	Sterilization	21
8.7.1	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	21
8.7.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™1	22
8.7.3	Ethylene oxide sterilization (EO)	22
8.7.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®	23
8.7.5	High-Level Disinfection	23
8.7.6	CIDEX (2,4 % glutaraldehyde solution)	23
8.7.7	Resert XL HLD (2,0 % hydrogen peroxide solution)	24
8.8	Reprocessing limits	25
9	Reprocessing (for US customers only) ...	26
9.1	Care and handling	26
9.2	Reprocessing	27
9.2.1	Reference chart	27
9.2.2	General warnings	28
9.3	Accessories	30
9.4	Water Quality Requirements	30
9.5	Pre-cleaning	30
9.5.1	Manual cleaning	31
9.6	Assembly, inspection and care	32
9.7	Packaging systems	33
9.8	Sterilization	33
9.8.1	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	33
9.8.2	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™1	34
9.8.3	Ethylene oxide sterilization (EO)	35
9.8.4	Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E® (only valid in USA)	35
9.8.5	High-Level Disinfection	36
9.8.6	CIDEX (2,4 % glutaraldehyde solution)	36
9.8.7	Resert XL HLD (2,0 % hydrogen peroxide solution)	37
9.9	Limits of reprocessing	37

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10 Maintenance.....	38
10.1 Maintenance and safety check	38
10.1.1 Maintenance	38
10.1.2 Safety check.....	38
10.2 Servicing and repair	39
10.3 Disposal.....	39
10.4 Repair program.....	40
10.5 Limitation of liability.....	41
10.6 Warranty	41
11 Compatible camera control unit (CCU)	42
12 Spare parts, recommended accessories.....	43
13 Information on electromagnetic compatibility (EMC).....	44
14 Subsidiaries	56

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2 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

2.1 Explanation of warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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2.2 General safety information

-  **WARNING:** Read this instruction manual thoroughly before using the equipment.
-  **WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.
-  **WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.
-  **WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.
-  **WARNING:** Patient leakage When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.
-  **WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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WARNING: All endoscopic devices (e.g., HF surgical probes) used in conjunction with KARL STORZ endoscopic cameras must comply with IEC 60601-2-18, the standard for medical electrical equipment – particular requirements for the safety of endoscopic equipment. The insulation of the endoscopic devices must be designed to withstand the highest HF voltage encountered during surgery.

HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.

WARNING: Risk of blindness. Never look into the open end of a light cable.

WARNING: Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.

WARNING: Test all equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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3 Intended Use Camera Heads

The camera head is used together with the camera control unit (CCU) and is attached to a telescope during endoscopic or microscopic procedures.

3.1 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

3.2 User qualification

The IMAGE1 S™ 4U camera head may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the IMAGE1 S™ 4U camera head.

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
- ① Instrument coupler
- ② Focusing ring/knob
- ③ User control button: for scrolling (downwards) through the menu or activating a camera function
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑤ User control button: for accessing the menu/ making a selection

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Number of products in the product packaging

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Federal (USA) law restricts this device to sale by or on the order of a physician

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6 Operating instructions for the IMAGE1 S™ 4U camera head

6.1 Frequency of the IMAGE1 S™ 4U camera head (50/60 Hz)

The IMAGE1 S™ 4U camera head is suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

i **NOTE:** The method for switching frequency is described in the instruction manual IMAGE1 S™ – TC 200, TC 300, TC 301 (96206286D).

6.2 Operating the IMAGE1 S™ 4U camera head

6.2.1 Connecting an endoscope and light cable

The IMAGE1 S™ 4U camera head has an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

! **CAUTION:** Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.

i **NOTE:** For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (97000213), which can be downloaded or requested at www.karlstorz.com.

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6.2.2 Adjusting the focus

The focus on the IMAGE1 S™ 4U camera head can be adjusted by turning the focusing ring/knob ② on the camera lens.

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7.2 Storage and Operating conditions

Refer to the compatible camera control unit (CCU) manual for relevant storage and operating conditions. The CCU and camera head were evaluated as a system for IEC 60601-1.

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7.3 Standard compliance According to IEC 60601-1; EC 60601-2-18.

Provides Type CF protection when used with KARL STORZ Camera Control Units that have CF symbol on connection port.

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7.4 Directive compliance According to the Medical Device Directive (MDD):

This medical device belongs to Class I
This medical device bears a CE mark in accordance with MDD 93/42/EEC.

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8 Reprocessing (not validated for US customers)

In the following chapter, you will find information on care and cleaning of the IMAGE1 S™ 4U camera head.

8.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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- Sterilization method has been validated.
- X On the grounds of potential damage, this reprocessing step must not be performed.

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**Reprocessing
(not validated for US
customers)**

8.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile upon delivery. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, national regulations regarding personal safety must be observed.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: Medical devices may be damaged by using chemicals that have not been approved by KARL STORZ.

CAUTION: National laws and regulations must be observed.

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i **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

! **CAUTION:** Inspect the camera head cable for breaks and cracks before and after reprocessing. Camera heads with damaged cables must NOT be reprocessed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

! **CAUTION:** Chemical residue on the camera head plug can cause failure of the camera electronics.

! **CAUTION:** The camera head must be thoroughly cleaned and disinfected/sterilized using validated infection control procedures before use and before every subsequent application. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.

! **CAUTION:** The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

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8.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

8.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

CAUTION: Inspect the camera head cable for breaks or cracks. Camera heads with damaged cables must not be sterilized or disinfected. Following wipe-down disinfection, return camera heads with damaged cables to KARL STORZ for repair.

8.4.1 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixing effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surfaces with a dry low lint cloth.

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8.4.2 Manual precleaning

8.4.2.1 Brushing the surfaces

Depending on the degree of soiling and the degree to which residues have dried, heavy soiling must be cleaned from the surfaces of instruments with the aid of a brush/sponge under cold running water until any visible contamination has been completely removed.

8.4.2.2 Auxiliary pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

8.4.3 Manual cleaning

The medical device must be completely immersed in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a brush or sponge. The device must then be rinsed with cold water to ensure neutralization.

8.4.4 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces and joints are dried completely with (preferably sterile) compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Wipe the camera head connector with a cotton tip applicator soaked in 70 % isopropyl alcohol to remove any residual disinfectant. Do not allow exposed glass windows to air dry. Glass surfaces can be cleaned using a cotton tip applicator soaked in 70 % isopropyl alcohol to prevent streaks and spots.
- Afterwards, a functional check must be carried out.

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i NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

8.6 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2 – 10, EN ISO 11607 Parts 1 + 2, DIN 58953). National regulations must be observed.

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8.7 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for the IMAGE1 S™ 4U camera head:

8.7.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

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NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S
- STERRAD® NX®
'standard' cycle
- STERRAD® 100NX®
'standard' cycle
'DUO' cycle

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8.7.2 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™1

Detailed information on the selection of the appropriate cycle for the various device generations is available from the STERIS® manufacturer.

CAUTION: V-PRO™ sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

8.7.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30 – 45 minutes.

WARNING: Medical devices must be dismantled for sterilization.

NOTE: Place the camera head in a sterilization tray. Do not place any instruments on top of the camera.

CAUTION: For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

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8.7.4 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

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CAUTION: Any deviations from the recommended parameters for disinfection must be validated by the user.

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CAUTION: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life.

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8.7.6 CIDEX (2.4 % glutaraldehyde solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4 % glutaraldehyde solution for a minimum of 45 minutes at 25 °C (77 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

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To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water

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Reprocessing (not validated for US customers)

for each rinse. Thorough rinsing of the camera head and the card-edge connector with sterile water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

8.7.7 Resert XL HLD (2.0 % hydrogen peroxide solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0 % hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

CAUTION: Do not leave the camera head in the solution for any longer than 60 minutes.

To finish, rinse the camera head three times with sterile water. Each rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head with sterile water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

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- Sterilization method has been validated.
- X** On the grounds of potential damage, this reprocessing step must not be performed.

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9.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: When preparing the cleaning and high-level disinfectant solutions, follow the manufacturer's instructions for dilution to proper concentration and temperature.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices.

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CAUTION: Inspect the camera head cable for breaks and cracks before and after preparation. Camera heads with damaged cables must NOT be processed. A break or cut in the cable will allow fluids to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

CAUTION: Chemical residue on the camera head plug can cause failure of the camera electronics.

CAUTION: The camera head must be thoroughly cleaned and high-level disinfected or sterilized using validated infection control procedures before use and before every subsequent application.

CAUTION: The use of sterilization/disinfection methods other than those recommended by KARL STORZ may have adverse effects on the materials of the camera and endoscopic adaptor. Any resulting damage will not be covered by the warranty.

CAUTION: Use only one method of high-level disinfection or sterilization per reprocessing cycle.

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9.4 Water Quality Requirements

Utility water:

Water as it comes from the tap that meets specifications as defined per AAMI TIR 34: 2014 Water for Reprocessing of Medical Devices. Utility water is acceptable for pre-cleaning and cleaning of the instruments.

Critical water:

Water that is extensively treated (usually by a multistep treatment process that could include a carbon bed, softening, DI and RO, or distillation) to ensure that the microorganisms and the inorganic and organic material are removed from the water; a final submicron filtration could also be part of the treatment process. Critical water is recommended for rinsing after high-level disinfection.

9.5 Pre-cleaning

Immediately after a procedure, at the point of use (e. g. operating room), wipe the camera head with a soft, lint-free disposable cloth moistened with water or diluted mild/neutral pH enzymatic cleaning solution to remove all heavy soiling.

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9.5.1 Manual cleaning

1. Thoroughly rinse the camera head using a commercial overhead sprayer nozzle with cold water to remove all gross debris. Rinse the camera head unit for a minimum of 1.5 minutes, paying special attention to the knobs, crevices, and seams. Then rinse the connecting cable for a minimum of 30 seconds.
2. Completely immerse the camera head in a diluted mild/neutral pH enzymatic cleaning solution (e.g. Enzol). Keep immersed for at least 1 minute.
3. While immersed, use a soft bristle brush (P/N: 27652/3) to brush all the surfaces of the camera head unit. Pay close attention to the seams between the instrument coupler, focusing ring and zoom ring. Do not brush the glass window.
4. Brush for a minimum of 1 minute.
5. Then, use a low-lint wipe to wipe the entire length of the connecting cable and the glass window of the camera head.
6. Finish by brushing the entire camera head, including the entire silicon cable, card-edge connector, paying special attention to the operating buttons, screws, knobs, and grasping mechanism. Rinse the bristle of the brush if heavy debris collects on them. Keep brushing until there is no visible debris on the bristles of the brush (minimum of 1 minute).
7. After the immersion and brushing period, remove the camera head from the cleaning solution and rinse using a commercial overhead sprayer nozzle with cold water. Rinse the camera head unit for a minimum of 1.5 minutes, paying special attention to the knobs, crevices and seams. Then rinse the connecting cable for a minimum of 30 seconds.

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Reprocessing (for US customers only)

9.6 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning procedure once more.
- The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.
- Wipe the camera head connector with a cotton tip applicator soaked in 70% isopropyl alcohol to remove any residual disinfectant. Do not allow exposed glass windows to air dry. Glass surfaces can be cleaned using a cotton tip applicator soaked in 70% isopropyl alcohol to prevent streaks and spots.

Before sterilization and/or disinfection, coil the camera head cable into loops with a diameter of at least 6 inches (15 cm). In doing so, do not kink or twist the cable.

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**Reprocessing
(for US customers only)**

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9.8 Sterilization

The following sterilization methods have been validated and approved by KARL STORZ for IMAGE1 S™ 4U:

9.8.1 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

CAUTION: STERRAD® sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

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9.8.3 Ethylene oxide sterilization (EO)

1. Place the camera head in an FDA-cleared sterilization tray.
2. Double wrap the sterilization tray with polypropylene wrap.
3. Start sterilization under the following parameters:

Conditioning parameters

Temperature: 55 °C (131 °F)
Humidity: 70 % RH
Conditioning dwell time: 30 minutes

Sterilization parameters

Sterilant: 100% ethylene oxide
Temperature: 55°C (131°F)
Gas concentration: 735±30 mg/l
Exposure time: 180 minutes

Aeration parameters

Time: 12 hours, minimum
Temperature: 55°C (131°F)

9.8.4 Chemical low-temperature sterilization with peracetic acid – STERIS® system 1E® (only valid for USA)

1. Place the camera head in one of the following trays:
 - C1200E SS1E General Processing Tray/ Container
 - C1220E SS1E Directed Flow Processing Tray/ Container
 - C1140E SS1E Flexible Endoscope Processing Tray/Container
 - C1160E SS1E Universal Flexible Processing (Open) Tray
2. Load the tray into the SS1E and start the sterilizer.

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9.8.5 High-level disinfection

WARNING: High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

CAUTION: Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life.

The following high-level disinfection methods have been validated and approved by KARL STORZ:

9.8.6 CIDEX (2.4 % glutaraldehyde solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.4% glutaraldehyde solution for a minimum of 45 minutes at 25°C (77°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

CAUTION: Do not leave the camera head in the solution for longer than 60 minutes.

To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head and the card-edge connector with water is very important for removing any residual disinfectant.

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9.8.7 Resert XL HLD (2.0 % hydrogen peroxide solution)

In order to prevent potential damage, do not place camera head in a container with other instruments. Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20°C (68°F). It must be ensured that the surfaces are covered and that no air bubbles are present.

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To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head with water is very important for removing any residual disinfectant.

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9.9 Limits of reprocessing

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

Always perform a functional test with the full IMAGE1 S™ system prior to using the camera head. The system should not be used if any of the accessories are damaged or non-functional.

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10.1 Maintenance and safety check

10.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

10.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the inscriptions relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents measured according to IEC 62353 (current edition), and checked for compliance with the limits specified in the standard.

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Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual (see Chapter 7.3.12 of the IMAGE1 S™ CCU instruction manual no. 96206286 D).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

10.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

10.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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10.4 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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10.5 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

10.6 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**Kompatible Kamera-
Kontrolleinheit**

**Compatible camera control
unit (CCU)**

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**11 Compatible camera control
unit (CCU)**

Item	Cat. No.
IMAGE1 S™ 4U-LINK in combination with IMAGE1 S Connect™	TC 304
or IMAGE1 S Connect™ II	TC 200
	TC 201

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**12 Spare parts,
recommended accessories**

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533 TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530 BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96286003D
Instruction manual IMAGE1 S™	96206286D

* For safe and comfortable working with
KARL STORZ cameras, we recommend the use
of sterile drapes.

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Information on electromagnetic compatibility (EMC)

13 Information on electromagnetic compatibility (EMC)

The CCU and camera head were tested as a system.

WARNING: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Please observe the EMC instructions in this instruction manual during installation and operation.

WARNING: The IMAGE1 S™ 4U camera head is only intended for use by physicians and medical assistants. This is a CISPR Class A medical equipment/system. In a domestic environment, IMAGE1 S™ 4U camera head may cause radio interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the unit or filtering the connection to the public mains network.

WARNING: The use of portable and mobile RF equipment may have impact on this and other pieces of medical equipment.

WARNING: The use of accessories or cables with the IMAGE1 S™ 4U camera head other than those specified in this manual may result in increased emissions or decreased immunity of the IMAGE1 S™ 4U camera head. When using accessories or cables with the 4U camera head other than those specified in this manual, it becomes the responsibility of the user of the 4U camera head to determine compliance with IEC 60601-1-2 when using this item.

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Table 1		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The IMAGE1 S™ 4U camera head is intended for operation in an electromagnetic environment as defined below. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic enforcement – guidance
RF emissions acc. to CISPR 11	Group 1	The IMAGE1 S™ 4U camera head uses RF energy only for its internal function. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The IMAGE1 S™ 4U camera head is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ 4U camera head is intended for use in the electromagnetic environment specified below. The user of the IMAGE1 S™ 4U camera head should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Air discharge	± 6 kV Contact discharge ± 8 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV/1 kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV Differential mode ± 2 kV Common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	<5 % U_T * (>95 % dip of the U_T) for ½ cycle 40 % U_T (60 % dip of the U_T) for 5 cycles 70 % U_T (30 % dip of the U_T) for 25 cycles <5 % U_T (>95 % dip of the U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ 4U camera head require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ 4U camera head be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the IMAGE1 S™ 4U camera head further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
 *NOTE: U_T is the mains alternating voltage before applying the test levels.			

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Table 4			
Guidance and manufacturer's declaration – electromagnetic immunity			
The IMAGE1 S™ 4U camera head is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ 4U camera head should ensure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™ 4U camera head, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
	NOTE: At 80 MHz and 800 MHz, the higher frequency range applies.		
	NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.		
^a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.			
^b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.			
^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 4U camera head is used exceeds the applicable RF compliance level above, the 4U camera head should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 4U camera head.			
^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6 Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™ 4U camera head			
The IMAGE1 S™ 4U camera head is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the IMAGE1 S™ 4U camera head can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ 4U camera head as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to < 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to < 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: lcano@karlstorz.com.mx

KARL STORZ Endoscopia México S.A de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ GmbH & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jiangmenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ GmbH & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com