# STORZ

**KARL STORZ—ENDOSKOPE**

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**1 Important information for  
users of KARL STORZ  
devices**

**It is recommended to check the suitability of  
the product for the intended procedure prior  
to use.**

Thank you for your expression of confidence in  
the KARL STORZ brand name. Like all of our  
products, this product is also the result of years  
of experience and careful manufacturing. You  
and your organization have decided in favor of a  
modern high quality product from KARL STORZ.

This instruction manual is intended to serve as an  
aid in the proper setup, installation, and operation  
of the ENDOFLATOR® 40. All of the necessary  
details and actions are clearly explained.

Please read these instructions carefully. Keep this  
manual in a convenient and conspicuous location  
and in its protective casing close to the device.

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**3 Controls, displays,  
connectors,  
and their uses**

- ① Power switch ("O" = off)
- ② Touch screen
- ③ Insufflation connection to patient
- ④ Gas connection (American connection)
- ⑤ SCB connectors
- ⑥ Service interface (concealed)
- ⑦ Ethernet interface (concealed)
- ⑧ Potential equalization connector
- ⑨ Power cord socket
- ⑩ Line fuse holder
- ⑪ Holder for CO<sub>2</sub> bottle  
(optional – cat. no. UI005)

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	Follow instructions for use
	Power on
○	OFF
	Potential equalization connector
	Applied part type CF
~	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).
	Manufacturer

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**Symbols on label and packaging of the device**

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. This can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).

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
	Settings
	Working area
	Start/Stop
	Sensitive Mode
	High Flow Mode
	Gas bottle
	House connection
	Flow
	Pressure
	Audio signals on
	Audio signals off
	House connection, input pressure acceptable
	House connection, input pressure too high or too low
	Gas bottle; input pressure above 40 bar
	Gas bottle; input pressure between 40 and 30 bar
	Gas bottle; input pressure between 30 and 10 bar
	Gas bottle; input pressure below 10 bar
	Gas bottle; input pressure too high
	Display of gas consumption / Reset the display by clicking on the -0- button

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
4.3 Gas filter label symbols

	Manufacturer
	Article catalog no.
	Product must not be reused
	Product must not be resterilized
	CE mark
	Do not use the product if the packaging has been damaged
	Consult instructions for use
	In the USA, federal law restricts this device to sale or use by or on the order of a physician
	Store in a dry place
	Store out of direct sunlight
	Sterilized with ethylene oxide
	Use by date of the product
	Production batch label

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7.3.6	CO <sub>2</sub> bottles equipped with German-standard or ISO connection.....	21
7.3.7	CO <sub>2</sub> bottles equipped with PIN-Index connection .....	21
7.3.8	Opening the valve of the CO <sub>2</sub> bottle, if necessary .....	22
7.3.9	Connecting to the central gas supply.....	22
7.4	Commissioning .....	22
7.4.1	Commissioning for the first time.....	22
7.4.2	Normal commissioning.....	25
7.4.3	Test for proper functioning .....	27
7.4.4	Preparing insufflation .....	28
7.4.5	Performing CO <sub>2</sub> insufflation .....	29
7.4.6	Decommissioning.....	30
7.5	Settings .....	31
7.5.1	Device settings .....	31
7.5.2	Service.....	31
7.5.3	System Log .....	31
7.5.4	Administration .....	32
7.5.5	Device information .....	33
<b>8</b>	<b>Maintenance .....</b>	<b>34</b>
8.1	Fuse replacement .....	34
8.2	Reprocessing .....	35
8.2.1	Reprocessing of the ENDOFLATOR® 40.....	36
8.2.2	Reprocessing of insufflation tube set with gas filter for single use, sterile .....	37
8.2.3	Reprocessing of reusable insufflation tube.....	37
8.2.3.1	Preparation for cleaning and disinfection .....	37
8.2.3.2	Manual precleaning .....	37
8.2.3.3	Manual cleaning .....	38
8.2.3.4	Manual disinfection.....	38
8.2.3.5	Machine cleaning and disinfection.....	39
8.2.3.6	Assembly, inspection and care.....	39
8.2.3.7	Packaging systems .....	41
8.2.3.8	Sterilization.....	41
8.2.3.9	Reprocessing limits .....	41
8.3	Maintenance and safety check.....	42
8.3.1	Maintenance .....	42

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8.3.2	Safety checks/ repeat inspections as per IEC 62353.....	42
8.4	Servicing and repair .....	43
8.5	Disposal.....	43
8.6	Repair program.....	44
8.7	Important information.....	44
8.8	Limitation of liability .....	45
8.9	Warranty .....	45
<b>9</b>	<b>Technical description.....</b>	<b>46</b>
9.1	Signals.....	46
9.1.1	Visual signals .....	46
9.1.2	Acoustic signals.....	47
9.2	Troubleshooting .....	49
9.3	Technical data.....	51
9.3.1	Standard compliance (for UI400).....	52
9.3.2	Directive compliance (for UI400).....	52
9.4	Technical documentation .....	53
<b>10</b>	<b>Spare parts, recommended accessories.....</b>	<b>54</b>
10.1	Spare parts/accessories .....	54
10.2	Accessories .....	55
<b>11</b>	<b>Electromagnetic Compatibility (EMC) information .....</b>	<b>56</b>
<b>12</b>	<b>Subsidiaries .....</b>	<b>70</b>

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## 5 General information

### 5.1 Description of the device

The ENDOFLATOR® 40 is an insufflation device for universal application in laparoscopic and thoracoscopic examinations and operations as well as endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy). In addition, the ENDOFLATOR® 40 is designed to expulse ambient air in open and endoscopically assisted cardiac surgery as well as endoscopic vessel harvesting. The very latest technologies for pressure and flow measurement and control allow for various operating modes, including modes tailored to specific situations.

The 'SENSITIVE' and 'HIGH FLOW' modes enable various applications to be accommodated. Easy handling and clearly arranged controls together with several safety circuits ensure maximum patient safety.

For particularly sensitive applications, the device features a 'SENSITIVE' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min).

To quickly compensate for the considerable loss of gas occurring in complex laparoscopic operations, the ENDOFLATOR® 40 has been designed for a high flow rate of up to 40 l/min.

The ENDOFLATOR® 40 is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface. The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

Thus the current patient pressure, current flow and the amount of gas used during the current intervention are displayed simultaneously during operation.

#### 5.1.1 Property rights

This product is protected in the USA by (at least one of the following) US Patent(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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## 6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 6.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this instruction manual, the following text must be carefully reviewed to ensure the safe and effective operation of this device. In addition, to make the signal words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** A Note provides special information regarding the operation of the product, or clarifies important information.

**WARNING:** Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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**Safety instructions**  
Warnings and cautions

-  **WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.
-  **WARNING:** Set up the device out of reach of patients.
-  **WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if
  - they are identified as such in the respective instruction manuals or
  - the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1).
-  **WARNING:** Grounding reliability can only be achieved when the device is connected to a 'Hospital only' or 'Hospital Grade' outlet. Check the plugs and cables regularly and do not use if damaged.
-  **WARNING:** To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.
-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.
-  **WARNING:** For safety reasons, do not work with a leaky device. There is a risk that the intraabdominal pressure will increase in an uncontrollable manner.
-  **WARNING:** Always rinse out the system with 1l CO<sub>2</sub> before use.

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**Safety instructions**  
Warnings and cautions

-  **WARNING:** Make sure there are no external forces acting on the cavity. This can increase the intracavitary pressure and/or cause pressure fluctuations.
-  **WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.
-  **WARNING:** Always use the lowest possible pressure and flow for insufflation.
-  **WARNING:** Always check which operating mode is set. For safety reasons, use the 'Sensitive' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.
-  **WARNING:** Any treatment may only be performed if visual observation of the action of the device is ensured.
-  **WARNING:** After 24 h of continuous operation, a self-test must be performed. Take note of the warnings to ensure that the device is fully functional.
-  **WARNING:** The device must not be used if either the touch screen or display is defective.
-  **WARNING:** If the touch screen does not react when pressed, you must turn the device off.
-  **WARNING:** Do not press several points on the touch screen simultaneously.
-  **WARNING:** Only sterile CO<sub>2</sub> gas may be used.
-  **WARNING:** You must secure CO<sub>2</sub> bottles connected up to the device to prevent them from falling over.
-  **WARNING:** Follow the safety instructions when changing the CO<sub>2</sub> bottle.
-  **WARNING:** The CO<sub>2</sub> bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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**Safety instructions**  
Warnings and cautions

-  **WARNING:** Always unplug the device from the power supply before carrying out any cleaning and maintenance work.
-  **WARNING:** Danger of electric shock! Do not open the device. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unauthorized persons will void warranty claims.
-  **WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.
-  **WARNING:** Only operate the device under the specified environmental conditions.
-  **WARNING:** Always keep a spare device to hand in case the first device should fail.
-  **WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.
-  **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.
-  **WARNING:** Always use a hydrophobic bacterial filter to prevent cross-contamination of patients.
-  **WARNING:** You must insert a sterile CO<sub>2</sub> filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. Replace the CO<sub>2</sub> filter after every use.

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 **WARNING:** Ensure that liquid cannot flow back into the device. Position the device higher than the patient and use a hydrophobic filter between the insufflation connection and insufflation tube.

 **WARNING:** Check the packaging of the tube set for damage! If the packaging is damaged, the sterility of the tube set may be impaired and the tube set must not be used.

 **WARNING:** Follow the country-specific regulations/laws for the disposal of accessories.

 **WARNING:** Avoid allowing fluids to enter the housing. Do not store liquids on or directly above the device. If, despite the precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies for the formation of condensation).

 **WARNING:** The device may only be operated at the voltage stated on the rating plate.

 **WARNING:** When replacing fuses, use only fuses of the correct rating.

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Before using the device on the patient for the first time it is imperative that you be acquainted with how the device operates and is controlled.

## 6.2 Intended use

CO<sub>2</sub> insufflators and their accessories are used to create and maintain a cavity during diagnostic or therapeutic endoscopic interventions as well as to expulse ambient air during open, endoscopically assisted or endoscopic interventions.

### 6.2.1 Indication

The insufflator is used to create a cavity in the following diagnostic and therapeutic interventions:

- Laparoscopy
- Thoracoscopy
- Endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy)

In addition, the ENDOFLATOR® 40 is also designed to expulse ambient air in open and endoscopically assisted cardiac surgery and endoscopic vessel harvesting.

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Use of the device in fields other than those indicated above is not allowed for safety reasons.

The ENDOFLATOR® 40 may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

**WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.

**WARNING:** Unauthorized conversions or modifications to the device are not allowed for safety reasons.

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**WARNING:** Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

### 6.2.2 Contraindications

No contraindications relating directly to the products used for CO<sub>2</sub> insufflation are currently known.

The use of an insufflator is contraindicated if, in the opinion of a qualified physician, such an application would pose risks to the patient, e.g. due to the patient's general condition.

The device is contraindicated for hysteroscopic applications.

### CO<sub>2</sub> absorption

CO<sub>2</sub> absorption occurs during insufflation. The body absorbs a portion of the CO<sub>2</sub> used for the insufflation. A CO<sub>2</sub> concentration in the blood or in the respiratory system which is too high can lead to the death of the patient in isolated cases. Observe vital functions throughout the entire procedure and ensure sufficient patient respiration. High pressure or high gas flow encourage CO<sub>2</sub> absorption.

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**Changes to circulatory, metabolic and respiratory functions**

The following symptoms can occur during CO<sub>2</sub> insufflation:

- Hypercapnia
- Acidosis
- Decrease in cardiac output
- Decrease in venous return

**Embolism**

Incorrect positioning of the insufflation instrument or a high intraabdominal pressure can lead to a CO<sub>2</sub> embolism. Avoid high pressures and verify the correct positioning of the insufflation instrument.

**Hypothermia**

The CO<sub>2</sub> insufflation can lead to a decrease in body temperature. Monitor the body temperature throughout the intervention.

**Idiosyncratic reactions**

The risk of a metabolic imbalance due to increased CO<sub>2</sub> absorption is increased for patients with sickle cell anemia or pulmonary insufficiency.

**Dehydration**

For longer operations with high gas levels, insufflation can lead to a drying out of the tissues or cause tissue damage. Unnecessary leaks must be avoided for this reason.

The decisive factor for the incidence of these complications is the length of the pneumoperitoneum, the level of the intraabdominal pressure and the CO<sub>2</sub> absorption associated with it. Increased pressure and flow values must thus be avoided. Pressures in excess of 15 mmHg are only required in isolated cases.

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### 6.3 User qualification

The ENDOFLATOR® 40 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

### 6.4 User profile of physician and assistants

Recognized medical skills in the field of application (specialist physician, qualified medical staff).

Adequate powers of comprehension to rationally assess the situation in hand.

Familiarity with the use of technical apparatus.

Suitable command of one of the languages used by the device and in the instruction manual.

Must have been trained thoroughly at least once in the use of the device.

No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

### 6.5 Intended conditions of use

#### 6.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

#### 6.5.2 Other intended conditions

- Frequency of use: one or more times a day.
- Length of use: from a few minutes to several hours a day.
- Place of installation: positioning on a level, vibration-free surface.
- Mobility: can be moved if positioned on a cart.
- Combination: can be used on the patient at the same time as other equipment required for the operations.

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- Control: can be controlled via the KARL STORZ SCB.
- The device is placed outside the sterile area.

**6.5.3 User position**

- To operate the device, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$ , at a distance of approx. 30-70 cm to the front panel.
- For observation of the actual values during the application, a visual distance from the device of 2 m is assumed.

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
### 6.6 Safety precautions at the installation site

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

**! WARNING:** DO NOT use the device in the presence of flammable anesthetics.

**! WARNING:** The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

**! WARNING:** The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

### 6.7 Safety precautions when operating the device

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

This device has been tested and found to comply with the EMC limits of the Medical Device Directive 93/42/EEC (EN 55011 Class B and IEC 60601-1-2).

Please observe the Electromagnetic Compatibility information in section 11.

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## 6.8 Safety features

The ENDOFLATOR® 40 has the following safety features:

### 6.8.1 Self-test

A self-check procedure is performed each time the device is switched on (for this the device must be supplied with gas). If the test is negative, a corresponding error message is shown and the device can no longer be used as normal.

The following components are tested:

- Leakproofness of the system (valves)
- Checking pressure measurement
- Checking flow measurement

### 6.8.2 Monitoring during operation

- **Overpressure:** If, during insufflation, the intracorporeal pressure rises above the set setpoint value by a limit dependent on the setpoint value for more than 3 seconds, an acoustic information signal sounds (see chapter 9.1.2) and the gas supply is stopped. After a time period which can be set (5...99 seconds), the overpressure is actively reduced via a blow-off valve.
- **Mode-specific pressure monitoring**
- **Monitoring of gas supply**  
If the system recognizes a lack of CO<sub>2</sub> gas, an optical signal is displayed.
- **Monitoring of gas preheating**  
The internal gas preheating is monitored by a second circuit, independently of the control. If a fault arises the heating is switched off.

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- Passive safety measures
  - The pressure and flow measurement systems have a redundant stand-by.
  - Both pressure control stages and the flow measurement are equipped with safety pressure release valves.
  - The function of the mechanical high pressure regulator is monitored.

**i** **NOTE:** If the monitoring system recognizes an error which renders safe functioning seemingly impossible, the device is transferred into a "safe state". From this state, the device can only be reset by turning off the power supply.

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**7 Installation and operating  
instructions**

**7.1 Unpacking the equipment**

Carefully remove the ENDOFLATOR® 40 and the accessories from their packaging. Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

**7.2 Basic equipment  
KARL STORZ Set no.: UI400S1**

- 1 ENDOFLATOR® 40 UI400
- 1 Power cord 400A
- 5 Insufflation tube sets with gas filter for single use, sterile
- 1 Universal wrench 20 4000 30
- 1 SCB connecting cable 20 0901 70, length 100 cm
- 1 Instruction manual
  
- 1 CO<sub>2</sub> supply tube (customized)

**7.3 Installing and connecting the  
device**

**i NOTE:** The ENDOFLATOR® 40 including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

**! WARNING:** Position the device higher than the patient.

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1. Place the device on a flat surface.  
**! WARNING:** Ensure that there is sufficient air circulation for the device.  
**! WARNING:** The device must be set up in such a way that the power cord can be simply pulled out at any time.

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- 7.3.1 Connecting the ground line**
- ! WARNING:** The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.
  2. Connect the potential equalization cable to the potential equalization connector.
  3. Where necessary, have grounding performed by qualified personnel.

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- 7.3.2 Connecting the power cord**
- ! WARNING:** The device may only be operated at the voltage stated on the rating plate.
  4. Push the power cord into the power socket as far as it will go.  
**! WARNING:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.  
**! WARNING:** The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.  
**! WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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### 7.3.3 KARL STORZ SCB

**i** **NOTE:** To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets. Connect the other end of the cord to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB devices (see KARL STORZ SCB control NEO System Instruction Manual).

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**7.3.4 Affixing the CO<sub>2</sub> bottle to the device (optional)**

Attach the optionally available holder to the rear of the device. Use the enclosed screws of length 6 mm.

Insert the CO<sub>2</sub> bottle vertically into the holder and close the locking clip.

**WARNING:** Only sterile CO<sub>2</sub> gas may be used.

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**7.3.5 Connecting the CO<sub>2</sub> bottle**

Connect the gas outlet of the CO<sub>2</sub> bottle to the device's gas connection using the high-pressure tube supplied.

**WARNING:** CO<sub>2</sub> bottles connected up to the device must be secured to prevent them from falling over.

**WARNING:** Follow the safety instructions when changing the CO<sub>2</sub> bottle.

**WARNING:** The CO<sub>2</sub> bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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**7.3.6 CO<sub>2</sub> bottles equipped with German-standard or ISO connection**

Place the high-pressure tube (e.g. 20 4002 22; see the section 'Accessories') onto the outlet of the CO<sub>2</sub> bottle and screw tight using the universal wrench provided.

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**7.3.7 CO<sub>2</sub> bottles equipped with PIN-index connection**


Place the connecting piece of the high-pressure tube (e.g., 20 4000 22; see the section 'Accessories') onto the CO<sub>2</sub> bottle and screw this tight.

Use the wrench provided (see the illustration) to open the CO<sub>2</sub> bottle.

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**7.3.8 Opening the valve of the CO<sub>2</sub> bottle, if necessary**

Slowly open the valve stopcock of the CO<sub>2</sub> bottle by turning it approx. one-half turn counterclockwise.

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**7.3.9 Connecting to the central gas supply**

Connect the low-pressure tube (e.g., UI001; see the section 'Accessories') to the device's gas connection and connect the other end directly to the wall outlet of the central gas supply (3.3...7 bar).

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**7.4 Commissioning**

**i NOTE:** The gas supply must be connected when the device is switched on, otherwise an error message is emitted.

**7.4.1 Commissioning for the first time**

The ENDOFLATOR® 40 is delivered as follows:

- (Screen) Language: English
- High Flow mode (pressure in mmHg)

**i NOTE:** Changing/adjustment of the settings for language, audio and the gas source are described over the next pages.

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1. Press the power switch ①.  
➤ Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo. After this, a white screen appears.

The device is now running a self-test. Once the test has been successfully completed, an acoustic availability signal sounds and a touch screen appears.

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- If no gas is connected or if the input pressure is too low, the warning 'Input pressure too low' appears.  
Now check whether the bottle is open.

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- If the gas supply was secured with a delay, the switch-on test is delayed and the user is notified via message 110 that the tests are still active.

**ⓘ NOTE:** Normal use is only possible after successful conclusion of the tests.

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If the self-test fails, the screen **Safe State** appears with an information message and an acoustic information signal sounds. Verify that the device is properly connected and that the preparations have been made according to the instructions.

Further information on troubleshooting can be found in the chapter 'Technical description - Troubleshooting'.

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The **Working screen** is displayed after the switch-on message.


**NOTE:** In routine operation, insufflation can be started and the device settings changed from here.

3. Press the button ① (wrench) to change to the settings mode.

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### Changing the language and other presets

**NOTE:** The screen language is preset as English.

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
The **Device Settings** screen appears with the list of available user settings.

- Language
- Audio
- Source of gas supply

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5. Select your preferred language under *Language*.

Confirm the selection with ① or reject the changes with ②.

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6. Select the desired volume for warning signals (information signals) and key tones (from 0 ... 4) under *Audio*.

Confirm the selection with ① or reject the changes with ②.

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7. Select the desired gas supply source under *Source of gas supply*.

Confirm the selection with ① or reject the changes with ②.

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#### 7.4.2 Normal commissioning

After successfully starting (see section 'Commissioning for the first time'), the *Working screen* appears.

**i NOTE:** In routine operation, insufflation can be started and the device settings changed from the *Working screen*.

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1. Attach the sterile CO<sub>2</sub> gas filter to the insufflation connection ③.

**WARNING:** You must insert a sterile CO<sub>2</sub> filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. The CO<sub>2</sub> filter must be replaced after each use.

2. Attach the insufflation tube to the CO<sub>2</sub> gas filter and connect the other end of the insufflation tube to the VERESS needle.

**NOTE:** In the case of the disposable insufflation tube set, the filter is firmly connected to the tube.

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#### 7.4.3 Test for proper functioning

**WARNING:** Test the device prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

#### Pressure stability and leakproofness

1. Set the mode to 'Sensitive'.
2. Use 15 mmHg and 2.0 l/min as the setpoint values.
3. Check the connection of the patient tube to the device.
4. Close the stopcock on the instrument.
5. Start insufflation.

➡ Pressure then builds up in the tube which increases above the setpoint value and which, after a few seconds, is released via the overpressure valve (blow-off valve). This is repeated cyclically.

6. Observe the instrument's conduct over a 30 second period.

**NOTE:** The flow display may only show values greater than zero for short periods. If the set pressure is only reached slowly even when the instrument's stopcock is closed or if you note stationary flow rates of 0.2 l/min or more, there is a leak in the system filter – patient tube – instrument.

7. In this case, check the connections and, if necessary, change the tube set.

**WARNING:** Do not work with a leaky device. There is a risk that the intra-abdominal pressure will increase in an uncontrollable manner.

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### 7.4.4 Preparing insufflation Enter the operating parameters for the insufflation mode\*, pressure and flow.

1. In the *Working screen* click on the 'High Flow' or 'Sensitive' button to select the mode.
2. Set the setpoint values for the pressure and flow using the ± buttons.

\*The ENDOFLATOR® 40 is equipped with two different insufflation modes:

- High Flow ('man' symbol)
- Sensitive ('baby' symbol)

**i NOTE:** To compensate for the considerable loss of gas occurring in complex operations, the High-Flow mode has been fitted with a flow rate of up to 40 l/min. For particularly sensitive applications, the device features a 'SENSITIVE' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min).

**! WARNING:** Always check the mode set. For safety reasons, use the 'Sensitive' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.

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**7.4.5 Performing CO<sub>2</sub> insufflation**

**WARNING:** Before starting insufflation rinse the entire system with 1 l CO<sub>2</sub>.

**WARNING:** Always use the lowest possible pressure and flow for insufflation.

3. Reset the volume display using the '0' button.
4. Insert the VERESS needle and open the locking lever. Press the START/STOP button to initiate the insufflation process.

**NOTE:** In order to use the maximum device flow in the High Flow mode, HICAP® instruments with a suitable connection must be used. The LUER-lock connection fitted to the end of the insufflation tube must be removed in this case.

**NOTE:** In the case of the disposable tubing set, cut the tube off at the mark behind the LUER-lock.

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### Displays

The insufflation pressure, flow and volume of the gas used are continuously displayed in the *Working screen*.

If the patient pressure rises above the preselected value, an acoustic information signal will be sounded. After the set delay (see *Settings* – Administration screen), the pressure is released via the overpressure valve.

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
### 7.4.6 Decommissioning

Close the valve on the CO<sub>2</sub> bottle on ending the intervention (when using a high-pressure gas supply) and remove the insufflation tube (incl. filter) used from the device.

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### 7.5 Settings

- Device Settings
- Service
- System Log
- Administration
- Device information

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#### 7.5.1 Device settings

The Device Settings are described in the section on commissioning for the first time.

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#### 7.5.2 Service

The **Service** screen is protected by a password.

**NOTE:** The *service programs* are only available to *service personnel*. The *settings* are described in the *Service Manual*.

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#### 7.5.3 System Log

The System Log contains the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/info information ID.

The most recent entry is allocated the number 000.

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**7.5.4 Administration**

- Export system log
- Unit pressure
- Date and Time
- Overpressure valve time delay
- Change administration password

**i NOTE:** The Administration screen is protected by a password. On delivery the password is '2132'.

It is the responsibility of the user to change this if necessary.

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- Export system log
- For service purposes, data can be exported via the service interface ⑥ using the USB stick.

1. Open the cover on the rear of the device with a Torx screwdriver (size T10).
  2. Insert a USB stick into the USB interface.
  3. Select Export system log on the Administration screen.
- Exporting of the system log starts.

**i NOTE:** The current status is shown via a progress bar. This progress bar shows that the export is underway, however it does not show how far it has advanced.

4. Screw the cover back on once export is complete.

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

- Unit pressure
- Select the desired unit (mmHg or hPa) under Unit pressure.

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- Date and Time  
The setting of date and time is self-explanatory.  
Confirm the selection with ① or reject the changes  
with ②.

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
- Overpressure valve time delay  
Select the desired time (5...99s) under  
**Overpressure valve time delay**.  
Confirm the selection with ① or reject the changes  
with ②.

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**i NOTE:** If there is overpressure for longer  
than the set time, this is reduced. If the  
overpressure cannot be reduced, the device  
emits a warning.

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### 7.5.5 Device information

The following elements are displayed:

- Device serial number
- Operating hours
- Software versions
- Software license information

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**8 Maintenance**

**8.1 Fuse replacement**

1. Turn the device off and disconnect from the power supply.

2. Remove the screw bits on the line fuse holder ⑩ using a screwdriver or another suitable tool.

**! WARNING:** Only use fuses of the correct rating.

3. Insert new fuses.

	100...240 V
Line fuse	2 x T2.0AH250V

4. Now put the line fuse holder ⑩ on again.
5. Restore the power supply.
6. Switch the device on and ensure that the start screen appears.

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## 8.2 Reprocessing

- UI 400 – ENDOFLATOR® 40
- UI004 Insufflation tube, can be sterilized

### General warnings

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

**WARNING:** Risk of infection: The reusable insufflation tube is not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect the reusable insufflation tube for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the reusable insufflation tube before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**WARNING:** Risk of infection/product damage: KARL STORZ has validated the procedures given in the instruction manual. Only use these procedures. The reprocessing method must be validated by the user on site.

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**WARNING:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**CAUTION:** Risk of damaging the medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at [www.karlstorz.com](http://www.karlstorz.com).

**WARNING:** National laws and regulations must be observed.

**NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

### 8.2.1 Reprocessing of the ENDOFLATOR® 40

**WARNING:** Always disconnect the device from the mains before cleaning!

**WARNING:** Liquid must not be allowed to enter the housing under any circumstances.

### Manual wipe-down disinfection of device and supply tube

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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**8.2.2 Reprocessing of insufflation tube set with gas filter for single use, sterile**

**WARNING:** Insufflation tube set with gas filter for single use, locks and hydrophobic bacteria filters must not be reused.

Sterile disposable products are not approved for reprocessing.

**8.2.3 Reprocessing of reusable insufflation tube**

**WARNING:** The reusable insufflation tube must be cleaned, disinfected and sterilized prior to use and subsequent reuse using validated reprocessing procedures.

**Disassembly of the connectors**

As far as possible, separate the reusable insufflation tube into its individual components prior to cleaning and disinfection.

**8.2.3.1 Preparation for cleaning and disinfection**

The reusable insufflation tube must be kept moist for transportation to the reprocessing station in order to prevent contaminations from drying on or reduce this as far as possible.

**8.2.3.2 Manual precleaning Immersion in cold water for 5 minutes**

To dissolve dried organic residues, immerse the reusable insufflation tube in cold water (10-20 °C) for 5 minutes.

**Brushing the surfaces**

Any visible contamination or heavy soiling must be completely removed from all external surfaces by cleaning under cold running water with the aid of a brush (Art. no. 27652) or a sponge.

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**Brushing the lumina/hollow spaces**

In order to remove all visible contamination, the lumina of the connectors and the tube ends of the reusable insufflation tube must be precleaned under cold running water with the aid of a suitable brush (Art. no. 27650A/27650B).

**Rinsing out with a pressure water gun**

The lumina of the reusable insufflation tube must be rinsed out using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 10 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) with the special attachment (Art. no. 27660F) is ideal for this purpose.

**8.2.3.3 Manual cleaning**

The reusable insufflation tube must be completely immersed in a cleaning solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically using a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

**8.2.3.4 Manual disinfection**

The reusable insufflation tube must be completely immersed in a disinfectant solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. Following the exposure time, the reusable insufflation tube and connectors must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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### 8.2.3.5 Machine cleaning and disinfection

Prior to machine cleaning and disinfection, carry out the steps of manual precleaning.

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

#### Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A<sub>0</sub> value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

#### Connecting

In order to guarantee effective machine cleaning and disinfection, connect the LUER-Lock to the washer and disinfectant using the LUER connector provided. Place the reusable insufflation tube on the rinsing lance in order to ensure that it is rinsed out completely. In addition, put the plastic connector with the tube adapter side to the top on the rinsing lance.

**i** **NOTE:** If necessary, the reusable insufflation tube must be dried off afterwards by hand.

### 8.2.3.6 Assembly, inspection and care

**! WARNING:** A visual and functional check of the reusable insufflation tube must be carried out!

**! WARNING:** Check the reusable insufflation tube for leaks after cleaning and disinfection. By doing this it is possible to detect a leak early, through which gas could escape.

**i** **NOTE:** Discolorations impede visual inspections for possible contamination.

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The cleaned and disinfected reusable insufflation tube must be visually inspected to ensure that it is clean, complete, dry and has not been damaged:

- If residues or contamination are still present, the reusable insufflation tube must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged medical devices must be withdrawn from use.
- There should be no holes, cuts, dents or scratches on the tube surface. If the reusable insufflation tube is damaged as described above, it must be replaced. If the reusable insufflation tube is discolored as a result of the reprocessing process, it must be removed immediately from the sterile supply cycle.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

The integrity of the reusable insufflation tube can be tested in the following way:

Attach a syringe to one end of the tube and kink the other end. Place the insufflation tube in completely demineralized water or microbiologically pure, sterile water. Force air into the insufflation tube using the syringe. If any bubbles appear in the water while forcing air through the reusable insufflation tube, the tube must not be reused.

**i NOTE:** The leakage tester 13242 XL may also be used to check for leaks.

To check for leakage use leak tester 13242 XL, for example. To do so, connect one tube end to the leakage tester and block off the other end, for example, with your finger. If there is a leak in the tube set, a pressure drop will be observed on the pressure gauge of the leakage tester.

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### 8.2.3.8 Sterilization

**WARNING:** Prior to the sterilization of the reusable insufflation tube, particularly in cases of manual reprocessing, it must be ensured that all chemical residue has been rinsed off or removed.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method is to be selected taking into account the specific national requirements and in consultation with the device manufacturers.

The following sterilization method has been validated and approved by KARL STORZ for this medical device:

#### Steam sterilization in the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 – 137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

### 8.2.3.9 Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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
### 8.3 Maintenance and safety check

#### 8.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Information regarding maintenance services can be obtained from your local representative or from the manufacturer.

#### 8.3.2 Safety checks/repeat inspections as per IEC 62353

**WARNING:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

#### Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

#### Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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**Test for proper functioning**

Carry out the functional check as per the instruction manual (see section 7.4.3).

**Documentation**

The safety check/repeat inspection and results must be documented.

**WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**NOTE:** For detailed information on the range and implementation of the safety check/ repeat inspection, please refer to the latest issue of the relevant Service Manual.

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**8.4 Servicing and repair**

**WARNING:** Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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**8.5 Disposal**

This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).

**WARNING:** At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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Servicehotline: +49 7461/708 980

E-Mail: [technicalsupport@karlstorz.com](mailto:technicalsupport@karlstorz.com)

**8.6 Repair program**

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen

Service hotline: +49 7461/708 980

E-mail: [technicalsupport@karlstorz.com](mailto:technicalsupport@karlstorz.com)

In other countries please contact your local KARL STORZ branch or authorized dealer.

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**8.7 Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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**8.8 Limitation of liability**

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the device has been used in accordance with its operating instructions at all times.

**8.9 Warranty**

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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## 9 Technical description

### 9.1 Signals

All of the signals created, with the exception of those which display the cause of an inoperable device, are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by an alarm signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized last will be signaled, i.e., the previous signal is overwritten.

#### 9.1.1 Visual signals

Signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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**9.1.2 Acoustic signals**


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

**i NOTE:** An acoustic signal is always issued as a complete sequence (regardless of how long the signal conditions apply).

The acoustic signal can be switched on or off in the 'Device settings' menu.

The acoustic signal can, however, only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

• Information signal (general)

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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

• **Availability signal**

On completion of the successful self-test after switching on the device, an availability signal with the following characteristics sounds:

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

• **Low gas supply**

When using a high-pressure supply, the drop in supply pressure shows that there is no more liquid CO<sub>2</sub> in the reservoir.

In addition to the visual messages described in the 'Visual signals' section, when the pressure levels 40 bar – 30 bar – 20 bar are fallen short of, a double gong is produced as an acoustic signal.

This is a sampled sequence. It lasts for 1.5 seconds.

If a minimum supply pressure is fallen short of, the information signal is output as specified both with a high-pressure supply and low pressure supply.

• **Key click**

Clicking on a key on the touch screen is indicated by a short beep. The volume of this beep can be set in the 'Device Settings' menu independently of all other information signals. The key click can also be permanently deactivated.

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**9.2 Troubleshooting**

**WARNING:** Always unplug the device before all maintenance work.

**Description of error:**

Complete failure of the device.

**Possible causes:**

1. Power supply failure.
2. Defective line fuse.

**Measure:**

1. Check that there is electricity to the wall outlet.
2. Change fuses as described in the instruction manual. Use only fuses of the correct rating.

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246: Contact KS service

The System detected a serious failure and established a safe state. You have limited access to settings menu. You need to restart the device.

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Device defective.

**Measure:**

Switch unit off and back on. If the warning message still appears after restarting, contact Service. Prior to this, open the System Log under Settings and note the ID, Word1, Word2 and information text of the last 5 entries.

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No gas flow.

**Measure:**

Observe the warnings (device temperature exceeded, temperature of high pressure regulator too low, gas supply not sufficient).

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No pressure build-up.

**Possible causes:**

1. Excessively large leak.
2. Excessively high instrument resistance.
3. Leaky insufflation system.

**Measure:**

1. Increase the flow setpoint value (for information: if the gas outflow from the cavity with the desired intraabdominal pressure is larger than the set flow setpoint value, the insufflation regulation goes to flow limitation and the cavity pressure setpoint value is not achieved).
2. Check whether it is possible to change the instrument. (for information: there are situations in which the insufflation pressure provided is not sufficient to build up the desired cavity pressure due to a high level of instrument resistance and simultaneous gas outflow. Depending on the instrument and leak conditions, the current values for pressure and flow are below the setpoint values.)
3.
  - a) Check the connections of the accessories and the leakproofness of the trocar. If the test for proper functioning is passed, the likelihood of a loose connection is very high.
  - b) Change the tube/tube set and perform a test for proper functioning.
  - c) If necessary, change the filter and repeat the test for proper functioning.
  - d) Contact service if the leak remains.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

- Type of IP protection: IP21

**According to IEC 60601-1-2, FCC-B:**

Observe the information on electromagnetic compatibility in chapter 11 (pages 56-69).

**9.3.2 Directive compliance  
(for UI400)**

**According to the Medical Device Directive (MDD):**

This medical product belongs to Class IIb.

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

**NOTE:** The code number after the CE mark indicates the responsible notified body.

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#### 9.4 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

**i NOTE:** We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Article	Cat. no.
<b>High-pressure tube</b> , American connector/ German connector, length: 55 cm	20 4000 21
<b>Same</b> , length 102 cm	20 4000 27
<b>High-pressure tube</b> , American connector/ PIN-Index connector, length: 55 cm	20 4000 22
<b>Same</b> , length 102 cm	20 4000 28
<b>High-pressure tube</b> , American connector/ American connector, length 102 cm	20 4000 25
<b>Same</b> , length 55 cm	20 4000 29
<b>High-pressure tube</b> , American connector/ ISO connector, length 102 cm	20 4002 22
<b>CO<sub>2</sub> bottle</b> , 1 l, empty, with German connector	26 4000 90
<b>CO<sub>2</sub> bottle</b> , 1 l, empty, with PIN-Index connector	26 4000 91
<b>Low-pressure tube</b> , for the central CO <sub>2</sub> gas supply, length 150 cm	UI001
<b>Same</b> , length 300 cm	UI002
<b>Same</b> , length 600 cm	UI003
<b>HD-Inline gas filter</b>	20 4000 32
<b>Universal wrench</b>	20 4000 30
<b>Line fuses</b> T2.0 AH 250 V, pack of 10	1858390
<b>Power cord</b> (grounded)	400 A
<b>Power cord</b> 'Hospital Grade' (USA)	400 B
<b>Instruction manual</b>	96116044 D

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**10.2 Accessories**

Article	Cat. no.
<b>Insufflation tube</b> , sterilizable, length 2.5 m, ID 9 mm	UI004
<b>Gas filter</b> , disposable, sterile, pack of 25	031122-25
<b>Insufflation tube set with gas filter</b> for single use, sterile, pack of 10	031200-10
<b>VERESS pneumoperitoneum needle</b> with spring-loaded blunt stylet, LUER lock, length 10 cm	26120 J
<b>Same</b> , length 12 cm	26120 JL
<b>Same</b> , length 7 cm	26120 JK
<b>Same</b> , length 15 cm for creating pneumoperitonea by means of DOUGLAS puncture with adipose patients	26120 JF
<b>Bottle holder</b> for the KARL STORZ CO <sub>2</sub> bottle	UI005
<b>Leakage tester</b>	13242 XL

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**11 Electromagnetic Compatibility (EMC) information**

**WARNING:** Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The KARL STORZ ENDOFLATOR® 40 Model UI400 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

The limits used are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The KARL STORZ ENDOFLATOR® 40 Model UI400 is a Group 1 device (as per CISPR 11).

Group 1 contains all of the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the individual devices
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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**Electromagnetic  
Compatibility (EMC)  
information**

**WARNING:** The device should not be used adjacent to or stacked with other equipment. If use adjacent to or stacked with other devices is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

**WARNING:** The use of accessories / transducers and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the device. The accessories / transducers and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories / transducers and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Cables which have been shown to comply with IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
SCB Cable	Yes	≥ 1	No	connecting SCB network, other SCB units
Power Cord	No	3	No	Power supply

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The ENDOFLATOR® 40 model UI400 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 40 model UI400 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ENDOFLATOR® 40 model UI400 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The ENDOFLATOR® 40 model UI400 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

<b>Table 2</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The ENDOFLATOR® 40 model UI400 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 40 model UI400 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T^*$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 seconds	<5 % $U_T^*$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

<p align="center"><b>Table 4</b> <b>Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting</b></p>			
<p align="center">The ENDOFLATOR® 40 model UI400 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 40 model UI400 should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF disturbances IEC 61000-4-6</p> <p>Radiated RF disturbances IEC 61000-4-3</p>	<p>3 V<sub>rms</sub> 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V<sub>rms</sub></p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the ENDOFLATOR® 40 model UI400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> $d = [3.5/3]\sqrt{P}$ <p><math>d = [3.5/3]\sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = [7/3]\sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

**Table 6**

**Recommended separation distances between portable and mobile RF communications equipment and the ENDOFLATOR® 40 model UI400**

The ENDOFLATOR® 40 model UI400 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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