# STORZ

## KARL STORZ — ENDOSKOPE

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**1. Important information  
for users of KARL STORZ  
instruments**

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the ENDOMAT® SELECT. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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



**3. Controls, displays,  
connectors,  
and their uses**

- ① Power switch
- ② TFT touch screen
- ③ Pump rollers
- ④ Light barrier
- ⑤ Pump lever (cartridge locking lever)
- ⑥ Connection socket (example) for footswitch
- ⑦ SCB connectors\*
- ⑧ Service interface (Ethernet)\*\*
- ⑨ Service interface (USB)\*\*
- ⑩ Power cord receptacle
- ⑪ Potential equalization connector
- ⑫ Disposable tubing cartridge

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\* The KARL STORZ-SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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**\*\*WARNING:** The ENDOMAT® SELECT is not intended to be connected to other equipment via the USB or Ethernet interface. These interfaces are exclusively reserved for KARL STORZ service. Connecting the ENDOMAT® SELECT to an IT network that includes other equipment could result in previously unidentified risks to patients, operators or third parties. Since the framework condition, installation, and operation of an IT network are the responsibility of the operator, KARL STORZ cannot guarantee failure-free operation. The operator must ensure the protection, safety, and reliability of the IT network through risk management in accordance with IEC 80001.

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

	Follow instructions for use
	ON/OFF (Standby)
	Applied part of type CF
	Connection socket (example) for footswitch
	Potential equalization connector
	Ethernet
	USB
	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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


**4.2 Symbols on label and packaging of the device**

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For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).

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4.3 User interface symbols

1	Name of discipline or field of application
2	Parameters which can be set (depending on the tubing cartridge inserted – prerequisite: installed software packages) <ul style="list-style-type: none"> <li>• Irrigation (ml/min) – LAP</li> <li>• Irrigation (mmHg) – ART</li> <li>• Suction (ml/min) – IBS® Shaver</li> </ul>
3	Increase/decrease the setpoint
4	Irrigation/suction display activated Actual value – white bar display Setpoint – orange marks and 'number'
5	Start/stop irrigation (flow or pressure) Standby / pump activated
6	Start/stop suction using the footswitch
7	only available in the disciplines ART/URO Boost pressure increase (via button or footswitch – in the example 'footswitch'). (in 10 % increments – from 10 % to 50 %)
8	Flow or pressure limitation (depending on the specialization and installed additional package (Advanced)) <ul style="list-style-type: none"> <li>• Flow limitation in ml/min</li> <li>• Pressure limitation in mmHg</li> </ul>
9	Alarm audio paused (30 s)
'no image'	
	Open the menu
	Cancel (reject changes)
	Confirm (accept changes)
	Menu back
	= scrolling in the menu
	Scroll through pages (symbol rotated through 90°)

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**4.4 Abbreviation used**

ART	Arthroscopy
BS	Bottle suction
CV	CLEARVISION®
CYST	Cystoscopy
DS	Direct suction
ENT	Ear Nose Throat
FC	Flow-controlled
FESS	Functional Endoscopic Sinus Surgery
GI	Gastrointestinal
Gyn	Gynecology
HYS	Hysteroscopy
IBS®	Intrauterine BIGATTI Shaver
LAP	Laparoscopy
NEURO	Neurosurgery
PC	Pressure-controlled
PCN	Percutaneous Nephroscopy
PRO	Proctology
RES	Resection
SPINE	Vertebral column
SURG	Surgery
THOR	Thoracoscopy
URO	Urology
URS	Ureterorenoscopy
VET	Veterinary Medicine

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ART	Arthroscopie
BS	Afzuiging via fles
CV	CLEARVISION®
CYST	Cystoscopie
DS	Afzuiging direct
ENT	HNO
FC	Door stroom geregeld
FESS	Functional Endoscopic Sinus Surgery (functionele endoscopische chirurgie aan de neusbijholten)
GI	Gastro-intestinaal
GYN	Gynaecologie
HYS	Hysteroscopie
IBS®	Intra-uteriene BIGATTI Shaver
LAP	Laparoscopie
NEURO	Neurochirurgie
PC	Druk geregeld
PCN	Percutane nefroscopie
PRO	Proctologie
RES	Resectie
SPINE	Wervelkolom
SURG	Chirurgie
THOR	Thoracoscopie
URO	Urologie
URS	Ureterorenoscopie
VET	Diergeneeskunde

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
<b>1. Important information for users of KARL STORZ instruments</b> .....	<b>II</b>
<b>2. Images of the device</b> .....	<b>III</b>
2.1 Operation as irrigation pump .....	III
2.2 Operation as suction pump .....	III
<b>3. Controls, displays, connectors, and their uses</b> .....	<b>V</b>
<b>4. Symbols employed</b> .....	<b>VI</b>
4.1 Symbols on the device .....	VI
4.2 Symbols on label and packaging of the device .....	VI
4.3 User interface symbols .....	VII
4.4 Abbreviation used .....	VIII
<b>5. General information</b> .....	<b>4</b>
5.1 Device description .....	4
5.2 Property rights .....	5
<b>6. Safety instructions</b> .....	<b>6</b>
6.1 Explanation of warnings and cautions .....	6
6.2 Intended use .....	10
6.3 User qualification .....	13
6.4 Training in the operation and function of the device .....	13
6.5 Patient profile .....	13
6.5.1 SURG, URO, ART, SPINE, ENT/NEURO, GI, VET .....	13
6.5.2 GYN .....	13
6.6 Intended region of treatment on the patient .....	14
6.7 User profile of physician and assistants .....	15
6.8 Intended conditions of use .....	15
6.8.1 Use .....	15
6.8.2 Other intended conditions .....	15
6.9 User position .....	16
6.10 Safety precautions at the site of installation .....	16
6.11 Safety precautions when operating the device .....	17
6.12 Safety features .....	17
6.12.1 Self-test .....	17
6.12.2 Testing the touch display .....	17
6.12.3 Monitoring the irrigation pressure .....	17

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
## 5. General information

### 5.1 Device description

The ENDOMAT® SELECT is a new roller pump from KARL STORZ. Depending on the application, it can be used for irrigation or suctioning away liquids during interventions in a range of medical fields.

The device automatically adapts to suit the type of operation being performed by providing optimal operating parameters when a field-specific tubing cartridge is inserted, assuming the device is activated for the field of application in question.

This concept has the advantage that the device can be configured specifically for the individual requirements of the user. For example, only applications which he wishes to use are displayed. If necessary, disciplines from the field of human medicine can be added at a later date (see Fig. 1). An 'Advanced software package' with extended setting options (see pages 26, 28, 29 and 34) can also be purchased.

Additional functions:

- BOOST: With certain URO/ARTHRO applications, the irrigation pressure can be increased briefly using the BOOST function.
- Pressure/flow limitation: With certain applications, in addition to the freely selectable main setting, e.g., the irrigation flow rate, the irrigation pressure can also be limited or vice versa.

The ENDOMAT® SELECT can also be configured as the model for veterinary medicine (Fig. 2).

However, it is not possible to combine disciplines from the field of veterinary medicine with disciplines from the field of human medicine.

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The device is operated and monitored via a touch display. The display for the setpoint and actual values of the irrigation pressure/flow enables monitoring of the current operating status. An electronic safety circuit interrupts irrigation or suction if device performance persistently deviates from the setpoint, and emits an audible signal to inform the user. The high safety standard of ENDOMAT® SELECT is further enhanced by an electronic auto-check system, which tests the system components each time the device is switched on, and notifies the operator of any failures detected.

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## 5.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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## 6. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**NOTE:** A Note indicates special information about operating the product, or clarifies important information.

**WARNING:** Read this instruction for use thoroughly before putting the device into operation. Read the section on safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**NOTE:** Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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## Safety instructions

### Warnings and cautions

**WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

**WARNING:** Set up the device out of reach of patients.

**WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.

**WARNING:** The instructions for use and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**WARNING:** Combinations of medical devices are only assured to be safe if:

- they are identified as such in the particular instructions for uses or
- the intended use and interface specifications of the devices used in combination permit this.

**WARNING:** The device may only be connected to the devices listed below, which have been shown to comply with the requirements of IEC 60601-1:

- CALCUSON
- UNIDRIVE® S III
- UNIDRIVE® S III ARTHRO

**WARNING:** Reel pump. Do not touch the pump reels when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

**WARNING:** The device controls the pressure measured at the outlet of the irrigation pump. The intracavitary pressure depends, among other things, on the height difference between the device and the cavity. Try to keep the difference as low as possible.

**WARNING:** Avoid allowing fluids to enter the housing. Do not store liquids on or directly above the device. If, despite the precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies for the formation of condensation).

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**Safety instructions****Warnings and cautions**

**WARNING:** Before you start the intervention, make sure that the tubing set used is not leaking or kinked and that the ends of the tubes are properly secured.

**WARNING:** Be carefully when inserting the tubing cartridge. A squeezed pump hose can cause the safety system to stop the pump.

**WARNING:** Servicing may only be performed by specialists authorized by KARL STORZ.

**WARNING:** Test this device prior to each surgical procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

**WARNING:** Follow the country-specific regulations/laws for the disposal of accessories.

**WARNING:** The device may only be operated within specified environmental conditions.

**WARNING:** The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind if the intervention requires controlled fluid transportation.

**WARNING:** To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Any opening of the device by unauthorized persons will void warranty claims.

**WARNING:** To avoid the risk of an electric shock, this equipment must only be connected to a supply mains with protective earth.

**WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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**WARNING:** Do not use in a magnetic resonance (MR) environment. This device is MR Unsafe.

**WARNING:** The ENDOMAT® SELECT is not to be used where MRI, CT, Diathermy, Electrocautery or RFID equipment may cause Electromagnetic interference (EMI). If the user is unsure or suspects that such equipment may be causing technical problems, the clinical/biomedical engineer of the facility should be consulted.

**WARNING:** To prevent electrical medical device equipment from being adversely effected by unknown sources of electromagnetic interference, it is highly recommended that Ad hoc radiated RF testing be conducted by a Clinical/Biomedical engineer or an EMC professional.

**WARNING:** Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

**WARNING:** The correct configuration of the device must be checked after each software update.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

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## 6.2 Intended use

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs, joints and operating fields as well as the suctioning off of irrigation fluids and bodily fluids, secretions, tissue and gas during diagnostic or therapeutic interventions.

### Indication

Suction and irrigation pumps with the corresponding accessories are suitable for all diagnostic or therapeutic interventions which require the energy and substances provided by the device as per section 9.6 'Technical data' for successful implementation.

The device provides irrigation and suction functions for a range of indications:

- Urology
- Gynecology
- Surgery (Thoracoscopy, Laparoscopy and Proctology)
- Arthroscopy
- Gastroenterology
- Lense Cleaning e.g. during FESS
- Spine surgery

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**Contraindication**

Use is contraindicated if, in the opinion of the attending physician, the device is not compatible with successful completion of the planned intervention due to its technical design.

The accessories must only be used with the intended devices. In general, medical devices must not be used on patients who are not part of the defined patient group, or if the operation itself is contraindicated.

The use of pump systems is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated, or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Pump systems must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

Suction and/or irrigation pumps without pressure monitoring must not be used for dilations of hollow organs or joints.

The device must not be used for ERCP and cholangioscopy (gastro), the lumbar spine (interlaminar approach), cervical spine as well as all applications which require a non-pulsatile flow.

Likewise, the device must not be used to administer medication. Nor may it be used for thoracic drainage.

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## Safety instructions

### Intended use

The ENDOMAT® SELECT is a roller pump which is intended for use either as an irrigation or suction pump in diagnostic or therapeutic interventions. The application must be carried out using appropriate suction/irrigation instruments. Use of the device in fields other than those indicated above is not allowed for safety reasons. The ENDOMAT® SELECT may only be operated with tubing cartridges or tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

**WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing system other than those specified are used.

The following serious defects may arise in particular:

- Inadequate grip on the tube connectors means a loss of sterility and may lead to the device being damaged by ingress of liquid.
- Premature breakage of the pump tube segment and consequent loss of sterility and damage to the device caused by ingress of liquid.
- Display of an incorrect actual irrigation rate due to deviations in tube diameter.

**WARNING:** No modification of this equipment is allowed.

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**6.3 User qualification**

The ENDOMAT® SELECT may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

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**6.4 Training in the operation and function of the device**

Training should only be provided by persons equipped to do so on the basis of their skills.  
Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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**6.5 Patient profile**

**6.5.1 SURG, URO, ART, SPINE, ENT/NEURO, GI, VET**

Gender: irrelevant  
Age: irrelevant  
Weight: irrelevant  
Nationality: irrelevant  
State of health: suitable for endoscopic interventions, in the physician's judgment

**6.5.2 GYN**

Gender: female  
Age: irrelevant  
Weight: irrelevant  
Nationality: irrelevant  
State of health: suitable for endoscopic interventions, in the physician's judgment

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**6.6 Intended region of treatment on the patient**

SURG: Abdomen, thorax, pelvis  
GYN: Uterus  
URO: Upper and lower urinary tract  
ART: Joints in the foot, knee, hips, shoulder, hand and finger  
SPINE: Thoracic and lumbar vertebral column  
ENT/NEURO: Lens cleaning with endoscopically assisted FESS surgery as well as with endoscopically assisted microneurosurgery and transnasal pituitary gland surgery  
GI: Upper and lower gastrointestinal tract

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**6.7 User profile of physician and assistants**

Training: recognized medical skills in the field of application (specialist physician, qualified medical staff)  
Intelligence: adequate powers of comprehension to rationally assess the situation in hand  
Ability to act: familiarity with the use of technical apparatus  
Language: Suitable command of one of the languages used by the device and in the instruction manual.  
Experience: at least one comprehensive session of training in the use of the device.  
Impairments: no physical impairment which makes it difficult or impossible to perceive acoustic or visual signals.

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**6.8 Intended conditions of use****6.8.1 Use**

The device is intended for use in clinics and practices. The technical data and ambient conditions are described in the instruction manual.

**6.8.2 Other intended conditions**

Frequency of use: one or more times a day  
Length of use: from a few minutes to several hours a day.  
Place of installation: Positioning on a level, vibration-free surface  
Mobility: can be moved if positioned on a cart  
Combination: can be used on the patient at the same time as other equipment required for the operations  
Control: can be controlled via the KARL STORZ-SCB.

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**6.9 User position**


For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$  at a distance of approx. 30 – 70 cm to the front panel.

For observation of the actual values during the application, a visual distance from the device of 2 m is assumed (the tube length is 2 m).

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**6.10 Safety precautions at the site of installation**

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

**! WARNING:** DO NOT use the device in the presence of flammable anesthetics.

**! WARNING:** The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use (Fig. 3).

**! WARNING:** The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

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**6.11 Safety precautions when operating the device**

It is the user's responsibility to make sure the device is safe and operates properly before using it.

During treatment with the ENDOMAT® SELECT the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic effects.

Any treatment may only be performed if visual observation of the device is ensured.

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**6.12 Safety features**

The ENDOMAT® SELECT is equipped with the following safety features:

**6.12.1 Self-test**

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

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**6.12.2 Testing the touch display**

After switching on the supply voltage, a KARL STORZ logo with the device name is displayed for a few seconds, depending on the loading time of the operating system (Fig. 4).

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**6.12.3 Monitoring the irrigation pressure**

An electronic safety circuit blocks the supply if the maximum adjustable target value is persistently exceeded by 20 mmHg.

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## 7. Installation and operating instructions

### 7.1 Unpacking the equipment

Carefully unpack the ENDOMAT® SELECT and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

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**i NOTE:** Depending on the installed software package, the range of functions can vary. Necessary accessories can be looked up in chapter 10.2.

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### 7.3 Installation

**i NOTE:** The ENDOMAT® SELECT including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

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1. Set the device on a horizontal, flat surface (Fig. 5).

**! WARNUNG:** The device must be set up in such a way that the power cord can be simply pulled out at any time.

**! WARNUNG:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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## Installation and operating instructions

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**! WARNING:** The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

- Connect the potential equalization cable to the potential equalization connector ⑪.
- Where necessary, have grounding performed by qualified personnel (Fig. 6).

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**! WARNING:** The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

Connect power cord and insert the power plug fully into the power cord receptacle ⑩ (Fig. 7).

**! WARNING:** Only operate the device with the voltage stated on the manufacturer's identification plate.

**! WARNING:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

**! WARNING:** Grounding reliability can only be achieved when the device is connected to a properly installed 'Hospital only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room).

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**! WARNING:** The ENDOMAT® SELECT is not intended to be connected to other equipment via the USB or Ethernet interface. These interfaces are exclusively reserved for KARL STORZ service.

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**7.4 KARL STORZ-SCB**

**!** **NOTE:** To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device. Therefore, to unplug the SCB cable, pull out the connector – do not pull directly on the cable!

Connect the SCB cable 20 0901 70 to an (equivalent) SCB socket ⑦ on the back of the device (Fig. 9). Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ-SCB control NEO System Instruction Manual).

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**7.5 Operating the device**

Switch on the device at its power switch ① (green indicator light illuminates) (Fig. 10).

➔ The device will then conduct a self-test, checking all system components for proper operation.

Case A: Once charging is complete, a start screen will appear, with the button: 'Please press here to continue' (Fig. 11).

Case B: Once charging is complete, a start screen will appear, with the button: 'Check height setting' (Fig. 12).

In Case B, the user presses the button to confirm that the position of the device matches the height difference setting (see Section 7.15.1 'Settings – height difference').

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
**!** **NOTE:** When the button on the start screen is pressed (Fig. 11/12), an acoustic ready signal will sound in order to check the alarm system.


**!** **WARNING:** For safety reasons, the device may only ever be used if, immediately after 'first click', the acoustic availability signal sounded!

**!** **WARNING:** Using advanced software package provides a setting for height difference between the device and the cavity. This setting corrects the pressure measurement and control. Check the setting after power on to avoid unintended over pressure.

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


## 7.6 Discipline/procedure

### 7.6.1 Human medicine

- After starting up the device, a screen is displayed for the user to select the discipline (Fig. 1).
- After selecting the discipline URO, for example, the screen 'Select procedure for URO' appears (Fig. 13).

### 7.6.2 Veterinary medicine

- After starting up the device, a screen is displayed for the user to select the discipline (Fig. 2).
- After selecting the discipline VET ART, for example, the screen 'Select procedure for VET ART' appears (Fig. 14).

**i NOTE:** If a tubing set was selected which then only supports the light procedures, these are grayed out.

Following selection of the desired procedure, an animation explains how to insert the tubing cartridge (see Fig. 15, top of page 22).

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
**7.7 Inserting the tubing cartridge**

**!** **CAUTION:** Only use the tubing cartridges or tubing sets designed for the ENDOMAT® SELECT.

☛ If the inserted cartridge and selected procedure do not match, a corresponding message is displayed (Fig. 16).

**7.7.1 Inserting the tubing cartridge for  
the first time/startup**

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**7.7.2 Inserting the tubing cartridge for  
the second and subsequent times**

☛ The procedure last selected appears. If the matching cartridge is now inserted, the user is taken to the operating menu. Should you wish to change, press the button to return to the medical specialization selection.

☛ If the inserted cartridge and selected procedure match, the user is taken to the appropriate working level (if the pump lever is closed).

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**Disposable tubing cartridge**

**!** **WARNING:** Check the expiration date and the package for possible damage prior to use. Products past expiration date or with damaged packaging must not be used.

**!** **WARNING:** Disposable cartridges and tubes are intended for single use and must not be resterilized. Discard after use.

Insert the tubing cartridge and turn the lever counterclockwise (Fig. 17, 18).

**i** **NOTE:** Reusable tubing cartridges are also available. The instruction manual supplied with the reusable tubing cartridges must be observed.

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
Connect the correct end of the length of tubing to the irrigation-liquid bag (puncture needle) or to the irrigation port on the instrument (LUER-Lock) (Fig. 19).

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The preselected procedure (in this case "Ureterorenoscopy" ❶) is shown on the touch display (Fig. 20).

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**7.8 Removing the tubing cartridge**

Turn the lever clockwise to the 9 o'clock position and remove the tubing cartridge (Fig. 21).

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### 7.9 Connecting the tubing cartridge for suction

**Using the ENDOMAT® SELECT as a suction  
pump (IBS®/RES/CALCUSON)**

**WARNING:** Use of the **KARL STORZ pump together with a KARL STORZ morcellator or shaver is not approved for use in the United States.**

**Alternative A – indirect suction  
(suction for CALCUSON)**

When there are solid components in the suctioned liquid such as stone concretions (urology), there is a risk of perforating the pump head tubing. In this case, it is recommended to connect a suction bottle (see accessories) between the suction instrument and the ENDOMAT® SELECT (Fig. 22). Use disposable tubing set 031647-10 for indirect suction.

**CAUTION:** When there are solid components in the suctioned liquid, the pump head tubing can be perforated and the ENDOMAT® SELECT can be contaminated.

**NOTE:** When the suction bottle is filled with water, the air volume of the suction bottle is reduced. One advantage of this is that the liquid is suctioned faster on starting suction on the ENDOMAT® SELECT (smaller dead volume). A second advantage is that almost no fluid will be suctioned after the suction on the ENDOMAT® SELECT has been stopped.

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
1. Attach the inflow tube on the tubing cartridge (tube on top; see arrow, Fig. 23) with tube connector onto the connector on the cap where the wire tray is located (Fig. 24).
2. Place the outflow tube on the tubing cartridge (tube at bottom) into a collection container (Fig. 27).
3. Attach the blue end of the 2nd tube in the tube set (without cartridge) onto the free bottle connection tube (Fig. 25).
4. Connect the orange end of the tube to the instrument.

**i NOTE:** The suction bottle must be placed at the height of the patient to prevent it from filling up due to the height difference between the patient and the suction bottle.

**Alternative B – direct suction  
(IBS® and RES)**

Use the disposable tubing set 030647-10 for direct suction (see Fig. 26).

1. Connect the inflow tube (marked orange, see Fig. 26) with the connector on the instrument.
2. Connect the outflow tube (marked blue, see Fig. 27) to the collection container.

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## 7.10 Using the ENDOMAT® SELECT as an irrigation pump

Depending on the field of application (determined by the tubing cartridge), either the irrigation flow or irrigation pressure can be freely selected. For the other parameters (pressure or flow), a 3-level limitation can be specified.\*

**i** **\*NOTE:** If some setting parameters are not visible on your system, it is possible that the 'Advanced' package has not been enabled.

### 7.10.1 Specifying the irrigation flow (SURG, GI, VET SURG)

In the display field (4) the setpoint (digital numerical value with orange limitation triangles) and the actual value (white bar display) are shown for the irrigation flow (in ml/min) (Fig. 28).

If the Advanced package is installed, the limitation of pressure (in mmHg) is shown in display field (8).

The roller pump is activated and deactivated with the associated Start/Stop button (5).

**Setting the flow:** Set the setpoint in 100 ml/min increments using the ± buttons (3). This is not affected by the pump status.

SURG: 100...2,500 ml/min  
(ADVANCED: 100...3,500 ml/min)

VET SURG: 100...3,500 ml/min

GI: 100...1 000 ml/min

**Set the pressure limit (prerequisite: installed Advanced software package or VET software package):** Set the pressure limitation using the button ► (range 100 mmHg --> 300 mmHg --> 500 mmHg --> 100 mmHg...).

**i** **NOTE:** If the Advanced package is not installed, the default value for pressure limitation is 500 mmHg.

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ENT/NEURO: 50 – 65 – 80 – 95 – 110 – 130 ml/min

**7. 10.2 Specifying the irrigation flow  
(ENT/NEURO)**

In combination with irrigation sheaths, this roller pump is used to clean the distal lens (endoscope window) during the CLEARVISION® procedure. The pump delivers liquid via the distal lens for as long as it is activated by the footswitch (see Accessories). The continual oscillation of the pump head removes drops of liquid from the lens. When the footswitch is released, an ongoing cycle will be fully completed, meaning that liquid will still be suctioned into the sheath.

The display field ④ shows the target value (digital figure with orange limit markings) and the actual value (white bar display) for the irrigation flow (in ml/min) (Fig. 29).

The roller pump is activated and deactivated by the associated footswitch ⑥.

**Setting the flow:** Adjust the target value in the following steps using the ± buttons ③.

This is not dependent on the pump status.

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**7.10.3 Specifying the irrigation pressure  
(HYS, URO, ART, SPINE, VET ART)**

**! WARNING:** Always start with the lowest possible pressure needed to achieve the desired distension. Increase the distension pressure until there is a clear, liquid medium.

In the display field (4) the setpoint (digital numerical value with orange limitation marks) and the actual value (white bar display) are shown for the irrigation pressure (in mmHg).

**URO/ART/VET ART:** In the display field (7), the pressure increase (Boost mode) (in %) is shown (Fig. 30).

**URO/GYN:** In the display field (4), the safety threshold at 100 mmHg can be seen by means of a vertical white line and a 'locked' padlock symbol. In order to set values between 100 mmHg and 150 mmHg, the + button (3) will need to be pressed until the padlock symbol has completely filled up and is shown as 'opened' (Fig. 31).

**Advanced software package or VET software package installed (Fig. 32, 33):**

In the display field (8), flow limitation is shown (in ml/min).

The roller pump is activated and deactivated with the associated Start/Stop button (5) or footswitch button (6).

**Setting the pressure:** Set the setpoint in increments of 10 mmHg using the ± buttons (3) (range 20...150 mmHg). The range of 20...150 mmHg applies for all pressure-regulated procedures apart from VET ART 'Large Animal', for which a maximum pressure of 400 mmHg can be set (Fig. 33).

This is not affected by the pump status.

**Boost:** Using the button ►, set the (brief\*) increase in pressure as a percentage (range 10% - 20% - 30% - 40% - 50% -> 10%...).

**i \*NOTE:** The Boost run-on is 2 seconds. If the Advanced software package is installed, the run-on can be extended by up to 60 seconds.

**i NOTE:** When the ENDOMAT® SELECT is being used in combination with the UNIDRIVE® S III via control cable UP 006, the Boost function is activated as soon as the Shaver function is activated.

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**Set the flow limit ⑧ (prerequisite: installed Advanced software package or VET ART software package):**

Set the flow limitation using the button ► (Fig. 34).

HYS, URO, SPINE: 200 ml/min – 400 ml/min – 600 ml/min --> 200 ml/min...

ART, VET ART: 1,500 ml/min – 2,000 ml/min – 2,500 ml/min --> 1,500 ml/min...

**① NOTE:** If the Advanced package is not installed, the default values for flow limitation are: 400 ml/min (HYS & URO), 200 ml/min (SPINE) and 1,500 ml/min (ART).

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
**7.10.4 Starting the irrigation pump**

If the irrigation pump is activated, the white bar in the display field ④ shows the actual value of the irrigation flow (see Fig. 35) or the actual value of the irrigation pressure (see Fig. 36). The current setpoint value is indicated by two orange triangles in the bar display. The bar display reflects the entire setting range.

**① NOTE:** The settings are retained in the device's memory after the device is switched off.

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## 7.11 Using the ENDOMAT® SELECT as a suction pump

### 7.11.1 Specifying the suction flow (IBS® shaver, RES, CALCUSON VET SURG)

**! WARNING:** Use of the **KARL STORZ pump together with a KARL STORZ morcellator or shaver is not approved for use in the United States.**

In the display field **4**, the setpoint (digital numerical value with orange limitation marks) and the actual value (white bar display) are shown for the suction flow (in ml/min) (Fig. 37, 38).

The roller pump is activated and deactivated with the associated Start/Stop button **5** or footswitch button **6**.

**Setting the flow:** Set the setpoint in 20 ml/min increments (IBS®) or 100 ml/min increments (RES/CALCUSON) using the ± buttons **3**.

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300 ml/min – 1,000 ml/min

### 7.11.2 Starting the suction pump

If the pump is activated, the white bar in the display field **4** shows the actual value of the suction flow (see Fig. 37, 38). The current setpoint value is indicated by two orange marks in the bar display. The bar display reflects the entire setting range.

**i NOTE:** The settings are retained in the device's memory after the device is switched off.

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**7.12 Test for proper operation and air-relief of the tubing system**

**WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.


**NOTE:** Have a collecting container ready to catch irrigation liquid as it runs out.

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Press the Start/Stop button **5** and let the irrigation pump run until all of the air has been released from the tubing system (i.e. no bubbles).


**SURG and GI**

The irrigation flow rate indicated on the bar display **4** must agree with the flow-rate setpoint (see orange marks). As soon as the inflow valve on the instrument is closed, the roller pump must stop irrigation.

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**HYS/URO/ART/SPINE**

Slowly close the inflow valve on the instrument. As soon as the pressure indicated on the bar display **4** rises above the setpoint (see orange marks), the roller pump must stop irrigation.

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
**With suction**

Actuate the footswitch. Close off the suction tube (instrument tube) with your finger and check whether a vacuum forms (Fig. 41).

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### 7.13 Initiating irrigation

Press the Start/Stop button 5 to initiate irrigation. Depending on the application, the actual irrigation flow or irrigation pressure is shown continuously on the bar display 4 (Fig. 42, 43).

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### 7.14 Initiating suction

Press the footswitch 6 to initiate suction. The actual suction flow rate is shown continuously on the bar display 4 (Fig. 44).

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**7.15 Device menu**

☰ Press the 'Menu' symbol to call up the device menu (see Fig. 45).

ⓘ **NOTE:** This screen appears when the lever is in the 9 o'clock position.

In the 'Menu' window, you can select various setting and management fields.

☛ The *Menu* screen appears (see Fig. 46).

There are submenus in the *Menu* screen:

- Settings
- Event protocol
- Information
- Service (only authorized service technicians)

**7.15.1 Settings**

Press the word 'Settings' (see Fig. 46).

☛ The *Settings* screen with the item Language, Misc. Settings\* and Height Difference\* appears (see Fig. 47).

\*these items require the Advanced software package (see page 34).

• **Language**

Press the term 'Language' in the *Settings* window (see Fig. 47).

☛ The *Language* screen appears (see Fig. 48).

In this window, you can select the screen language. The following languages are available (Czech, German, English, Spanish, French, Italian, Portuguese, Russian, Polish, Swedish, Chinese, Japanese and Thai). If necessary scroll to page 2 or 3. If you do not choose a new language, only the Back button (arrow pointing left) is available. Pressing the Back button returns you to the previous screen.

Select the required language (see Fig. 48). Clicking on X cancels the action or rejects changes, clicking on ✓ confirms changes – the previous *Settings* screen appears.

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• **Misc. settings**  
(prerequisite: installed Advanced or VET software package)

Press the term 'Misc. Settings' in the *Settings* window (see Fig. 47).

- The *Misc. settings* screen appears with the following setting options (see Fig. 49):

- Brightness of the display (adjustment range: 1 ... 8)
- Volume of the button tones (adjustment range: 0 ... 5; 0 = 'mute')
- 'Boost time' – delay time/duration of increased pressure (adjustment range: 2 seconds ... 60 seconds; up to 10 seconds in 2-second increments, thereafter in 5-second increments)
- Pressure unit 'mmHg' or 'cmH2O'

If you do not make a new selection, only the 'Back' button is available.

• **Height difference**  
(prerequisite: installed Advanced or VET software package)

**! WARNING:** Using advanced software package provides a setting for height difference between the device and the cavity. This setting corrects the pressure measurement and control. Check the setting after power on to avoid unintended over pressure.

Press the term 'Height difference' in the *Settings* window (see Fig. 47).

- The *Height difference* screen appears (see Fig. 50).

**i NOTE:** To ensure precise pressure measurement, the device must be positioned level with the patient. If the device is above or below the patient, incorrect measurements must be anticipated due to the pressure of the water column.

Enter the height difference between the device and the patient (adjustment range: +6 ... -6; the difference corresponds to a device height of 110 mm).

If you do not make a new selection, only the 'Back' button is available.

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If the height difference is set to a positive value, the device will factor in the additional hydrostatic pressure at this height. This can be identified by means of an enduring actual value (white bar display), which can also be seen even when the pump is not activated.

If a negative height difference is set, the device will factor in the missing hydrostatic pressure and will add this to the measured pressure value. The target values are still limited to the maximum values with this correction.

If the set height difference does not match the actual height difference between the device and the patient, then the displayed pressure will not match the actual pressure.

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### 7. 15.2 Event protocol

Press the term 'Event Protocol' in the *Menu* window (see Fig. 46).

☛ The *Event protocol* screen appears (see Fig. 51).

The saved alarms and information messages can be found here. The maximum number of displayed entries is 200. In the scroll bar, the current page number out of the maximum pages present is shown (e.g., 1/34).

Each entry occupies one line in the list and consists of the date, time, and info ID. The most recent entry is in the top line on page 1. Tapping on a line opens the corresponding event log entry (Fig. 52).

The event log will be backed up in the event of voltage drops and when the device is switched off, and contains entries relating to switch-on and switch-off times.

The entire event log has a capacity of 50,000 entries. If the maximum number of entries is exceeded, the oldest entries will be overwritten by new ones.

All alarms and user information are saved in the event log together with the time of entry, and can be exported to the Service area.

### 7. 15.3 Device information

Press the term 'Device Info' in the *Menu* window (see Fig. 46).

☛ The *Device info* screen appears (see Fig. 53).

In this window, you will find information about your device (e.g., the serial number, software version, and operating hours).

### 7. 15.4 Service

Press the term 'Service' in the *Menu* window (see Fig. 46).

☛ The *Service Password* screen appears (see Fig. 54).

The Service menu is reserved for the authorized Service and, for this reason, access is password-protected.

'Cancel' calls up the 'Settings' window again.

The settings that are possible in the Service Menu are described in the Service Manual.

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### 8.1 Reprocessing

#### 8.1.1 Wipe-down disinfection of device

**! WARNING:** Always disconnect the device from the mains before cleaning!!

**! CAUTION:** Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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#### 8.1.2 General warnings for reusable medical devices

**! WARNING:** Reusable accessory must be cleaned, disinfected and sterilized before use and subsequent reuse using validated reprocessing procedures.

**! WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" instructions and the accompanying documentation.

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**Maintenance****Reprocessing**

**WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**WARNING:** Risk of infection/Danger of damage to medical devices: The methods listed in instructions for use have been validated by KARL STORZ. You must use only these methods. The reprocessing process has to be validated by the responsible organization on site.

**WARNING:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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**CAUTION:** Risk of damaging the medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A helpful guide is available on the Internet at [www.karlstorz.com](http://www.karlstorz.com).

**WARNING:** National laws and regulations must be observed.

**NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D) can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

### 8.1.3 Suction bottle and cap

Carefully pull the tubes from the suction lances and take the bottle out of the holder.

**NOTE:** Prompt transportation to the reprocessing station and a subsequent decontamination process must be ensured.

#### 8.1.3.1 Accessories

Accessories required for carrying out reprocessing: Brushes: 27652

#### 8.1.3.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

#### 8.1.3.3 Manual precleaning

##### Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

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**Brushing the surfaces**

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

**8.1.3.4 Manual cleaning**

The medical device must be completely immersed in a cleaning solution. To ensure bubble-free wetting, the lumens must be filled. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or a sponge. Finally, it must be rinsed with cold water to ensure neutralization.

**8.1.3.5 Manual disinfection**

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumens must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

**8.1.3.6 Machine cleaning and disinfection**

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D):

**Machine cleaning/thermal disinfection**

Thermal disinfection is preferred. The relevant national requirements and the  $A_0$  value must be taken into account when using this method.

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The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

**i NOTE:** If necessary, the instrument must be dried off afterwards by hand.

**Connecting**

In order to ensure effective machine cleaning and disinfection, as well as adequate rinsing out of the instrument, it must be connected up to the washer and disinfectant.

**8.1.3.7 Assembly, inspection and care**

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

**i NOTE:** During care procedures, use items from the catalog "Care, Sterilization and Storage Techniques".

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**8.1.3.9 Sterilization**

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

**The following sterilization methods have been validated and approved by KARL STORZ for this medical device:**

**Steam sterilization using the fractionated prevacuum procedure**

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C – 137 °C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

**i NOTE:** *Greased components must be disassembled for sterilization, so that steam penetration is ensured.*

**8.1.3.10 Limits of reprocessing**

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

**8.1.3.11 Storage**

Store disinfected suction bottles and caps in a closed cupboard to protect them from dust, moisture and temperature fluctuations.

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
## 8.2 **Maintenance and safety check**

**WARNING:** Always unplug the device before carrying out any cleaning and maintenance work.

### 8.2.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

### 8.2.2 **Safety check/repeat inspection according to IEC 62353**

**WARNING:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

### **Visual inspection**

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

### **Electric measurements**

- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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**Test for proper functioning**

A functional inspection must be carried out in accordance with the instruction manual (see section 7.12).

**Documentation**

The safety check/repeat inspection and results must be documented.

**WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**NOTE:** For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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**8.3 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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**8.4 Disposal**

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

**WARNING:** At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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**8.5 Repair program**

For devices, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In the case of repairs, please contact your local KARL STORZ subsidiary or authorized dealer.

In Germany you can refer repairs directly to

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen  
Service hotline: +49 7461/708 980  
E-mail: technicalsupport@karlstorz.com

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**Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### 8.6 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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### 8.7 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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## 9. Technical description

### 9.1 Alarm specification

The device issues a low-priority alarm '301: Maximum pressure' (Fig. 55).

The alarm signals are not self-retaining. This means that the alarm signal is only issued for the duration that the signal condition is present. The alarms are displayed for a minimum of 5 seconds, and at least one tone sequence is output.

The alarm thresholds and delays are pre-programmed as fixed settings.

The excess pressure alarm '301: Maximum pressure' (Fig. 55) indicates that excessive pressure was measured at the device output. The excess pressure alarm is triggered when the pressure exceeds 150 mmHg for a maximum of 2 seconds during pressure-regulated irrigation applications in urology and gynecology.

**WARNING:** Continued use of the pump without eliminating the condition that has lead to the excess pressure alarm can lead to serious undesirable effects for the patient.

#### 9.1.1 Visual signaling

The alarm is shown in the title line.

The alarm is shown as blue text on a cyan-colored background.

The alarm will overwrite any other text messages.

#### 9.1.2 Acoustic signaling

The alarm signal is a burst of 2 tones (tone 1 – frequency 320 Hz, tone 2 – frequency 254 Hz) lasting for 170 ms each, with a pause of 240 ms. The burst is repeated every 16 seconds. Four harmonic waves are generated (Fig. 57).

If an alarm condition is present, at least one complete tone sequence will be output. The audio signal of the alarm can be deactivated for 30 seconds by pressing the 'bell' (pause alarm audio).

The bell will be shown in yellow during this time (Fig. 56).

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### 9.1.3 Verification of the functioning of the alarm system

The functioning of the alarm system is verified by the sounding of the ready signal when the button is pressed on the start screen.

The individual alarm conditions can be checked as described in Section 9.3.

## 9.2 Information signals

The information signals are self-explanatory messages that explain the device behavior, provide the user with support with the individual operating functions and thus contribute to the usability of the device, or help the service technicians with troubleshooting.

Apart from those indicating cause for an inoperable device ('Safe State'), the information messages are not self-retaining.

This means that the signal is only issued for the duration that the signal condition is present. For very short-term signal conditions, the minimum display duration is 5 seconds.

The signals are ordered according to priority. Alarms have higher priority than information signals. Information signals are split into five different priorities.

A high-priority signal will overwrite an existing lower-priority signal, and a lower-priority signal will be suppressed for as long as there is a higher-priority signal.

In the event of multiple signal conditions with the same priority, the most recently detected condition will overwrite earlier conditions in the title line on the display.

### 9.2.1 Visual signaling

Information signals are also displayed in the title line.

An information signal is displayed as blue text on a white background (the opposite way round to a normal display).

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
**9.2.2 Acoustic information signal**

The acoustic information signal is a double tone c-e (263 Hz – 330 Hz) lasting 300 ms. Depending on the priority of the message, the tone sequence will either be output once, or will be output repeatedly. The volume is 63 dBA (Fig. 58).

In the event of messages indicating an inoperable device, or if the message '300' (high pressure) appears, then the tone sequence will be repeated every 20 seconds. In all other cases, it will only be output once.

• **Ready signal**

Once a self-test has been successfully carried out and the button on the start screen has been pressed, a ready signal with the following characteristics will sound:

The tone pitch is modulated by  $\pm 2$  Hz with a 1.5 Hz modulation frequency. Five harmonic waves are generated (Fig. 59).

• **Button press**

The pressing of a button on the touchscreen is indicated by a short beep. The volume of this beep can be adjusted in the 'Settings' menu independently of the volume of all other information signals. The button press tone can also be permanently deactivated.

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### 9.3 Test conditions for alarm and information signals

#### 9.3.1 Excess pressure alarm

In order to check the excess pressure alarm, which only exists in the disciplines of urology and gynecology, proceed as follows:

1. Switch on the device
2. Correctly insert the pressure-regulated tubing cartridge 031523-01
3. Choose one of the following procedures if available:  
Uro – CYST, PCN or URS  
Gyn – HYS
4. Set the target value to 50 mmHg
5. Securely attach the leakage tester (Fig. 60) or pressure cuff to the bottom tube (Fig. 61).
6. Pump to generate a pressure of up to 170 mmHg at the leakage tester

☛ The visual and acoustic pressure alarm is output.

#### 9.3.2 Information signals

To check the information signal, proceed as follows:

1. Switch on the device
2. Correctly insert the pressure-regulated tubing cartridge 031523-01
3. Choose one of the following procedures if available:  
Uro – CYST, Gyn – HYS, SPINE – LUMBAR, ART – Knee or VET ART – small animal
4. Set the target value to 50 mmHg
5. Securely attach the leakage tester (Fig. 60) or pressure cuff to the bottom tube (Fig. 61).
6. Generate pressure at the leakage tester
  - a. If using the procedures Uro – CYST or Gyn – HYS, pump to generate pressure of up to 60 mmHg
  - b. If using the procedures SPINE – LUMBAR, ART – Knee or VET ART – small animal, pump to generate pressure of up to 330 mmHg

☛ The visual and acoustic excess pressure warning is output.

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**9.4 Information messages**

<b>Brief description (1 line with approx. 40 characters possible)</b> <b>Possible cause</b>	<b>Response</b>
102: Position the tubing set again Last sensor test more than 24 h ago, deviation in pressure values	Remove the tubing cartridge briefly and then start the pump again.
150: Pump stopped Communication or pressure measurement disrupted	Device function interrupted, pump will need to be re-activated if error message is automatically reset.
180: Paused Transient discrepancy between redundant pressure measurement values	Transient sensor discrepancy, pump continues automatically
190: Cartridge detection error Cartridge detection electronic error	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
191: Cartridge locking error Cartridge locking electronic error	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
200: Hardware initialization error Device defect	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
203: C2B error Self-test: Communication error	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
20C: Restart necessary in < 4 h The device has been in continuous operation for more than 20 h.	After 24 h of continuous operation, the pump can no longer be restarted.
20D: Sensor test overdue The device has been in continuous operation for more than 24 h.	Turn device off and on again.
251: C2B error Internal communication error	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
252: Electronics error Internal communication error	Device function interrupted, pump will need to be re-activated if error message is automatically reset.
254: Default restored File error	Information message, work can continue, if necessary check the time setting
255: Electronics error IO error	Information message, it may be possible to continue work
256: Date reset RTC Battery empty	Information message, work can continue
258: USB overcurrent FE: USB power too high	Information message, work can continue, disconnect USB connection
259: Electronics error BE: Internal error	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.

**Information messages**

<b>Brief description (1 line with approx. 40 characters possible)</b>	<b>Response</b>
<b>Possible cause</b>	
300: High pressure Stopcock on instrument closed	Observe surgical field, ensure pressure reduction
301: Maximum pressure Excess pressure alarm in URO and HYS application if 150 mmHg is exceeded	Observe surgical field, ensure pressure reduction
500: Main functions not active FE: A touch screen activation lasting more than 25 s has been detected	Acoustic and visual signal – device can still be operated: A touchscreen calibration can be started by pressing the touchscreen for 25 seconds. This text message will appear after 25 seconds. Remedy: Check touchscreen and clean if necessary. This text message may also appear in the event of a touchscreen short circuit. In this case: Turn the device off and on again, and contact KARL STORZ Service if the error occurs again.
501: Touch calibration FE: A touch screen activation lasting more than 30 s has been detected	This text message is part of the touchscreen calibration procedure (see text message 500) and is displayed after the touchscreen is touched continuously for 30 seconds. The screen must be released and touched again within 5 seconds in order to start calibration.
0x750: Defective electronics Electronics error in the cartridge module	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x760: Defective drive Electronics error in the drive control unit	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x761: Defective control unit Electronics error in the general control unit	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x762: Defective H bridge Electronics error in the H bridge	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x763: Motor current offset Electronics error in the motor current measurement	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x764: Peak current Electronics error, overcurrent	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x765: Index missing Electronics error in the speedometer	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x766: Index moved Electronics error in the speedometer	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x767: Motor PWM Control error in the motor	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.

**Information messages**

<b>Brief description (1 line with approx. 40 characters possible)</b>	<b>Response</b>
<b>Possible cause</b>	
0x768: Direction change Electronics error in the speedometer	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x770: Defective electronics (PS) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x771: SW error (PS) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x772: PS error (PUST) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x773: PS error (ST) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x774: PS error (H off) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x775: PS error (H on) Error in the monitoring system	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x780: Voltage error Operating voltage outside of permissible range	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.
0x790: WD task error BE: Program crash	Turn the device off and on again, contact KARL STORZ Service if the error reoccurs.

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- Complete failure of the device.

*Possible causes:*

- Power cord connector is not properly connected to device socket.
- Power supply failure.
- Defective internal device fuse.

**WARNING:** The user can not replace a defective internal device fuse.

*Remedy:*

- Push power cord connector firmly into device socket.
- Check that there is electricity to the wall outlet.
- Contact Service.

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- Poor suction.

*Possible causes:*

- Leakage occurring in suction system.

*Remedy:*

- Check all suction lines, and replace any exhibiting evidence of leakage.
- Check the seating of the sealing cap.

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- No suction.

*Possible causes:*

- Float blocking suction inlet.
- Bacterial filter on suction bottle is moist and therefore impermeable.

*Remedy:*

- Check the fluid level in the vessel and empty, if required.
- Clean float and check for freedom of movement.
- Change the bacterial filter.

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- No irrigation pressure build-up.

*Possible causes:*

- Leakage in tubing or tubing incorrectly connected.
- Failure in device's control circuitry.

*Remedy:*

- Check all tubing, in particular the connections and replace if necessary.
- Send device in for repair.

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**Standard compliance** (for UP 210)

**According to IEC 60601-1, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:**

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF 
- Type of protection against moisture: Drip-water protection as per IP 21

**According to IEC 60601-1-2:**

Observe the information on electromagnetic compatibility in chapter 11 (p. 67–83).

**Directive compliance** (for UP 210)

**According to Medical Device Directive (MDD):**

This medical device belongs to Class IIb  
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

 **NOTE:** The code number after the CE mark indicates the responsible notified body.

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### 9.7 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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**10. Spare parts,  
recommended accessories**

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Item	Cat. no.
Tubing Set, Irrigation, PC, for single use, sterile, package of 10	031523-10*
Tubing Set, Suction, BS, for single use, sterile, package of 10	031647-10*
Tubing Set, Suction, DS, for single use, sterile, package of 10	030647-10*
Tubing Set, Irrigation, FC, for single use, sterile, package of 10	031524-10*
Tubing Set, Irrigation, CV, for single use, sterile, package of 10	031529-10*
Tubing Set, Suction, BS, reusable, sterilizable	UP 010
Tubing Set, Suction, DS, reusable, sterilizable	UP 009
Tubing Set, Suction, FC, reusable, sterilizable	UP 007
Tubing Set, Suction, PC, reusable, sterilizable	UP 008

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Item	Cat. no.
Leakage tester	13242 XL
<b>For IBS® shaver and RES recommended accessories</b>	
Suction bottle, 5 l, sterilizable	20 3000 50
Cap for suction bottle 20 3000 50, 5 l, sterilizable	20 3000 34
Metal filter	20 3000 38
Bottle stand for suction bottle, 5 l	20 3000 32
Bottle stand holder	20 3000 33
<b>Accessories required for operation with CALCUSON</b>	
Suction bottle, 0,5 l, sterilizable	20 3000 51
Cap, for suction bottle	20 3000 39
Bottle stand, for suction bottle	20 3002 31
Control cable	20 7010 70
Control cable for combination with UNIDRIVE® S III ARTHRO	UP 006
Control cable for combination with UNIDRIVE® S III or CALCUSON	20 7010 70
<b>Accessories required for ENT/NEURO</b>	
Single-pedal footswitch	20 0142 30
<b>Accessories required for GI / accessories recommended for boost confirmation</b>	
Single-pedal footswitch	20 0141 30
<b>For »Advanced« Package recommended accessories</b>	
SCB connecting cord, length 100 cm	20 0901 70

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## Electromagnetic Compatibility (EMC) Information

### 11. Electromagnetic Compatibility (EMC) Information

**WARNING:** Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and commissioning.

The ENDOMAT® SELECT model UP 210 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The ENDOMAT® SELECT model UP 210 is a Group 1 device (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**i NOTE:** The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**Electromagnetic  
Compatibility  
(EMC) Information**

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ENDOMAT® SELECT UP 210, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

Guidance and manufacturer's declaration – electromagnetic emissions		
<p>ENDOMAT® SELECT model UP 210 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOMAT® SELECT model UP 210 should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	ENDOMAT® SELECT model UP 210 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	<p>The ENDOMAT® SELECT model UP 210 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p> <p><b>– Professional healthcare facility environment –</b></p>
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions acc. to IEC 61000-3-3	Complies	

**Table 4 – ENCLOSURE PORT** (according to IEC 60601-1-2)

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 9	
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	

**Table 5 – Input a.c. power PORT (according to IEC 60601-1-2)**

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
Electrical fast transients/bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency	
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV	
Surges Line-to-ground	IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment
Voltage dips	IEC 61000-4-11	0 % $U_T$ ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
		0 % $U_T$ ; 1 cycle and 70 % $U_T$ ; 25/30 cycles Single phase: at 0° Grad	
Voltage interruptions	IEC 61000-4-11	0 % $U_T$ ; 250/300 cycles	

**Table 8 – Signal input/output parts PORT (according to IEC 60601-1-2)**

Phenomenon	Basic EMC standard	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	
Electrical fast transients/bursts	IEC 61000-4-4	±1 kV 100 kHz repetition frequency	
Surges Line-to-ground	IEC 61000-4-5	±2 kV	
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable since released only for professional healthcare facility environment

**Table 9 – Test specifications for ENCLOSURE PORT IMMUNITY to  
RF wireless communications equipment (according to IEC 60601-1-2)**

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ENDOMAT® SELECT model UP 210 may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga  
Ontario, L5N 3R3, Canada  
Phone: +1 905 816-4500, Fax: +1 905 858-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8526  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: info@kseam.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Edificio Atlantic, Oficina 3G  
Calle D e/ 1ra y 3ra  
10400 Vedado, Havana, Cuba  
Phone: +537 836 95 06, Fax: +537 836 97 76  
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.  
Lago Constanza No 326,  
Col. Granada Del. Miguel Hidalgo  
C.P. 11520 México D.F.  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Østensjøveien 15B  
0661 Oslo, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
12739 Skärholmen, Sweden  
14108 Kungens Kurva, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representation Office  
Kestucio g. 59  
08124 Vilnius, Lithuania  
Phone: +370 68567000  
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S  
Skovlyoftens 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)31 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +385 1 620 5880, Fax: + 386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE  
East Mediterranean and Gulf S.A.L.  
Block M, 3rd Floor  
Beirut Souks, Weygand Street  
2012 3301 Beirut, Lebanon  
Phone: +961 1 999390, Fax +961 1 999391  
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.  
P.O. 6061, Roggebaai 8012  
Cape Town, South Africa  
Phone: +27 21 417 2600, Fax: +27 21 421 5103  
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan  
Saryarka, 6, BC "Arman", off. 910  
010000 Astana, Republic of Kazakhstan  
Phone: +7 7172 552-549, 552-788, Fax: -444  
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE  
East Mediterranean & Gulf (branch)  
Building West Side 7A – Unit 7WA – 3008  
Dubai Airport Free Zone, P.O. Box 54983  
Dubai - United Arab Emirates  
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282  
Service Hotline: +971 (0)4 3415882  
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited  
11th Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG  
Resident Representative Office  
14th Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch,  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8 Jianguomenwai Street, Chaoyang District,  
Beijing, 100022, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shanghai Branch  
36F, Tower 1 Grand Gateway  
No. 1 Hongqiao Road, Xuhui District,  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
Chengdu, 6100414, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
Shenyang, 110014, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
Guangzhou, 510620, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.  
I-Park Bundang1, 102dong, 8F  
239 Jeongjail-Ro, Bundang-Gu  
Seongnam, 13556, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4th Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp  
www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)