





The **EASY CHECK** from Ritter at a glance:


- For the detection of SARS-CoV-2 in saliva
- Patient-friendly due to non-invasive sample collection
- Perfect for children, the elderly and disabled people
- Rapid test results after 4 minutes (depending on room temperature)
- BfArM-listed and producer has a

"Prequalification" at the WHO




Bundesinstitut für Arzneimittel und Medizinprodukte












One **EASY CHECK** pack contains:

- 20 test cassettes
- 20 pipettes
- 20 sampling tubes
- 20 disposable bags for collecting of sample material (saliva)
- Manual

(Excluded from the right of return like every test)

EASY CHECK TEST KIT PACKAGE







The EASY CHECK COVID-19 rapid antigen test can be performed with posterior oropharyngeal saliva samples.

For posterior oropharyngeal saliva samples: wash your hands with soap and water- or alcohol-based solutions.

1. Open the container. Clear your throat, loosen the saliva from the throat and spit it into the container (about 2ml). Avoid contamination of the outer surface of the container with saliva. The best time to collect the sample is after getting up before the patient has brushed, eaten or drunk their teeth.

HOW TO USE EASY CHECK Instructions for use:



2. Collect approximately 200µL of the saliva sample with the pipette.

3. Open the sample tube, add the sample and shake it to mix everything well.

4. Take the test device out of the pouch and place it on a table and cut off the protruding part of the sample tube lid. Then add 3 drops of the sample to the sample well and read Depending on the room temperature, you can see the result after 4-15 minutes. Fastest results at room temperature bet ween 18-30°C.

















MANUFACTURER INFORMATION




The company is a research and development focused Chinese biotechnology company that produces high quality medical In Vitro Diagnostic (IVD) rapid test kits as well as revolutionary customerspecific reagent kits developed, manufactured and delivered to all parts of the world.

The company was founded by a team of professionals with many years of combined engineering, marketing / sales, operations and manufacturing experience in the industry. Their in vitro diagnostic lateral flow kits screen for a wide range of targets including infectious diseases, tumors, cardiac abnormalities, substance abuse, and fertility.

Thanks to the comprehensive quality management system, which applies international standards (EN ISO 13485), a high quality of the test results and accuracy is guaranteed. Most of their products are CE and CFDA certified.















www.easycheck-test.de

Declaration of conformity

	EC Dec	laration of C	onformit	y				
Manufacturer: Name: JOYSBO (T Address: Tianjis Inte Medicine 9th floor, N China. Telt =88-022-053784 Essail: moltytikjoyd	mational Joint Acad io.220, Dongting Ra 115	Whene Authorized Representative: Name: Lona NJ, B, V Address: Koningin Julianaptein 10,1a Vord, 2595AA, The Hugar, Netherlands E-mail: petre@iotaard.com						
We, the manufactures	here with declars	that the product(s)						
Product Name	SARS-CoV-2 (Colloidal Gol	Antigen Rapid Test Kit di	Specification	20Tests/box (Test/hi +20 Bags), 40 Tests /box (Test/hig =40 Bags)				
Intended Use	denotity them and days of the unset	For in this qualitative detect of SARS-CoV-2 nucleoraged antipes in neuroball with specificant depecty from individuals who are neglected of COVED-19 by their heatmane prevaler within the first 3 days of the seast of the symptoms. This test is only provided the use by clinical laboratories or to heathbase workers for paint-of-care specing, will can for at home testing.						
Classification	Othen			-				
80 11493 80 14925 80 80 14925 80 80 1411 80 80 14113	1.2011	157 150 1813-1 201 157 13641 2002 150 11723-1 2018	- 15(13672-2082 3-23640-2013 62386-7-3013				
2		CE						
		sole responsibility that our pr the European Partiament and						
Medical Devices.								
We agete to develop.i	replement and main	ntain a documented post-proc	luction monitoring p	nocen.				
Name of Gene	ral Manager		10					
Signature		2.54						
Date		200	AL AND	and all				
Place		1.27	inim Chief					
Seal (Manufac	iarer)	120	出现市場公共 15	1				

EU registration "In Vitro Diagnostic Products



Ministry of Health, wellbeing and sports

> Return address PO Box 16114 2500 BC The Hague

Lotus NL B.V. Attn. Mr. X. Wei Koningin Julianaplein 10 2595 AA The Hague

Date : Aug 18, 2020 Subject : Notification In-vitro diagnostics

Dear Mr. Wei

I hereby acknowledge receipt on 29 April 2020 of the Article 4. 1st paragraph of Attachments the Dutch Decree in vitro diagnostics (BIVD) that company name JOYSBIO (Tianjin) Biotechnology Co., Ltd with European authorized Lotus NL B.V. market the product below as an in vitro diagnostic product on the European market.

The product is registered as an in vitro diagnostic under number:

SARS-CoV-2 IgG/Neutralizing antibody Rapid Test Kit(Colloidal Gold), SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold), Immunochromatography analyzer (no brand name) (NL-CA002-2020-53008)

Tuberculosis Antibody Test Kit (Colloidal Gold), Mycoplasma Pneumonia IgM Antibody Test Kit (Colloidal Gold), Treponema Pallidum Antibody Test Kit (Colloidal Gold), Morphine/Methamphetamine/Ketamine Test Kit (Colloidal Gold) (no brand name) (NL-CA002-2020-53009)

It means that you have fulfilled your obligation under Article 4 of the BIVD.

In all further correspondence regarding the above-mentioned product, I request that you state this number. No further rights can be derived from this number, it only serves to facilitate the administrative notification.

The registration of in vitro diagnostics as a medical device under the Classification Criteria (Annex II) to Directive 98/79/EC on medical devices for in vitro diagnostics is subject to possible revisions of European regulations on the classification of medical devices and to advancing scientific understanding (see Article 10 (1) of Directive 98/79/EC).

Farmate Visiting address: Hoftoren Riinstraat 50 2515 XP The Hague T 070 340 6161

http//hulpmiddelen.farmatec.nl

Information about: M.P. Meijer - Michiels

Medische hulpmiddelen@minvws.nl registration number: CIBG-20204011

Date of Application

Aug 13, 2020

letter

Correspondence should only be addressed to the return address stating the date and reference of this



CERTIFICATES / LISTINGS

FDA submission

FOM U.S. FOOD & DRUG

Acknowledgment Letter

9/11/2020

Hongyan Li JOYSBIO (Tianjin) Biotechnology Co., Ltd. Tianjin Tianjin TEDA 300457 CHENA

Dear Hongyan Li:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has received your submission. This submission has been assigned the unique document control number below. All future correspondence regarding this submission should be identified prominently with the number assigned and should be submitted to the Document Control Center at the above letterhead address. Failure to do so may result in processing delays. If you believe the information identified below is incorrect, please contact the Office of Product Evaluation and Quality (OPEQ) submission support at (301) 796-5640 or OPEOS/abmissionSupport/u fda hhs gov.

> Submission Number: EUA202733 Received: 9/11/2020 Applicant: JOYSBIO (Tianjin) Biotechnology Co., Ltd. Device: SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

We will notify you when the review of this document has been completed or if any additional information is required. For information about CDRH review regulations and policies, please refer to http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm.

Sincerely yours,

Center for Devices and Radiological Health

U.S. Fast & Drug Administration. 10902 New Hampshire Avenue Silver Spring, MD 20983 www.htia.gov

EN ISO 13485: 2016 Proof of quality

Certificate

The Certification Body of **TÜV Rheinland LGA Products GmbH**

hereby certifies that the organization JOYSBIO (Tianjin) Biotechnology Co., Ltd. Tianjin International Joint Academy of Biotechnology & Medicine 9th Floor No.220, Dongting Road, TEDA 300457 Tianjin P.R. China

has established and applies a quality management system for medical devices for the following scope:

(see attachment for scope)

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date 2020-06-07

Certificate Registration No SX 60143180 0001

An audit was performed Report No 16806278 004

This Certificate is valid until 2022-10-12


Certification Body





and the full the state of the end ware contract of the other war

TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg



CERTIFICATES / LISTINGS



Page 1 of 7

WHO Prequalification



20, AVENUE APPIA - CH-1211 GENEVA 27 - SWITZERLAND - TEL CENTRAL +41 22 791 2111 - FAX CENTRAL +41 22 791 3111 - WWW.WHO.IN

Prequalification Unit Inspection services WHO DESK ASSESSMENT REPORT **Emergency Use Listing (EUL)**

Review of Quality Management System Documentation

Part 1	General information					
Company information	1					
Name of	Joysbio (Tianjin) Biotechnology Co., Ltd					
manufacturer						
Corporate address	Tianjin International Joint Academy Biotechnology & Medicine					
of manufacturer	9th Floor, No. 220, Dongting Road, TEDA 300457 Tianjin, China					
Contact person	Ms Yang Man					
	Director Registration Department					
	Email: molly@joysbio.com					
	Tel: _+86-13821759311					
Manufacturing site(under assessment					
Address of						
manufacturing site	Same as above					
if different from						
that given above						
Desk assessment details						
Date of review	18-23 November 2020					
EUL number(s)	EUL 0582-223-00					
Inspector(s)	Conrad Mark					
Products covered	SARSV-2 Antigen Rapid Test Kit (Colloidal Gold)					
by this desk						
assessment						
List of documents	WHO-EUL Quality System Information 312 pages					
submitted						
Any documents						
missing?						
Abbreviations	Meaning					
NC	Non-conformity					
QC	Quality control					
QMS	Quality management system					

Part 2 Summary of the assessment of supporting documentation

1. Certification and audit reports:

ISO 13485:2016 certificate number SX 60143180 0001 was provided.

Organization:F

JOYSBIO (Tianjin) Biotechnology CO., Ltd.

Joysbio Biotech Co., Ltd, Tianjin, China-Dx-COVID 18-23 November 2020 This audit report is the property of the WHO Contact: prequalinspection@who.int

WHO Prequalification



20, AVENUE APPIA - CH-1211 GENEVA 27 - SWITZERLAND - TEL CENTRAL +41 22 791 2111 - FAX CENTRAL +41 22 791 3111 - WWW.WHO.INT

Tel. direct: +41 22 791 47 61 Fax direct: +41 22 791 47 30 Email: prequalinspection@who.int In reply please refer to: P5-447-3/KR/SL/1

Ms Yang Man Joysbio (Tianjin) Biotechnology Co. Ltd Tianjin International Joint Academy of Biotechnology & Medicine 9th floor, No. 220 Dongting Road TEDA 300457, Tianjin République Populaire de Chine

24 November 2020

Dear Ms Man,

Your reference

OUTCOME OF DESK ASSESSMENT EUL Emergency Use Listing WHO Prequalification Unit - Inspection Services Joysbio (Tianjin) Biotechnology Co. Ltd

Thank you for your email correspondence dated 12 October 2020 and the documents that were sent to the WHO PQT: Inspections Team for the Emergency Use Listing of SARS-CoV-2 (EUL) desk assessment of the Quality Management System of Joysbio (Tianjin) Biotechnology Co. Ltd. Kindly be advised that your application for a desk assessment was reviewed as described in the desk assessment report (enclosed). These related to the site, indicated as:

Name:	Joysbio (Tianjin) Biotechnology Co. Ltd
Address:	Tianjin International Joint Academy of Biotechnology & Medicine, 9th floor,
	No. 220 Dongting Road, TEDA 300457, Tianjin, China

The documents submitted for the desk assessment were found to be satisfactory and are considered to constitute adequate evidence of compliance with ISO 13485 and the requirements described in the "Instructions and requirements for Emergency Use Listing (EUL) submission: In vitro diagnostics detecting SARS-CoV-2 nucleic acid and rapid diagnostics tests detecting SARS-CoV-2 antigens (PQDx_ 347 version 4: 09 June 2020)"

Furthermore, this desk assessment allows Prequalification Inspection Team to recommend to the Prequalification Assessment Team that the site may be named as a manufacturing site in the dossier for the following product:

PQT Number	Product
EUL 0582-223-00	SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Please do not hesitate to send an email to prequalinspection@who.int should you require any further information regarding the closure of this inspection.

Yours sincerely,

Dr Joey Gouws Team Lead, Inspection Services Prequalification Unit Regulation and Prequalification Department Access to Medicines and Health Products Division





BfArM-Listing



Address: Tianjin International Joint Academy of Biotechnology& Medicine 9th floor, No.220, Dongting Road, TEDA 300457 Tianjin China TEL: 86-22-65378415 Web: en.joysbio.com E-mail: bd@joysbio.com

To: The Federal Institute for Drugs and Medical Devices 1/20/2020

Letter of Declaration

Confirmation of efficacy in novel mutations

We, JOYSBIO (Tianjin) Biotechnology Co., Ltd. (hereinafter "JOYSBIO"), with the address of Tianjin International Joint Academy of Biotechnology& Medicine 9th floor No.220, Dongting Road TEDA 300457 Tianjin, China, hereby declare that our product "JOYSBIO SARS-COV-2 Antigen Rapid Test Kit (Colloidal Gold)" is compatible with the new virus strain VUI - 202012/01.

This product is for qualitative detection of SARS-CoV-2 nucleocapsid antigen. We confirm this information is available on the IFU under section [Verwendungszweck], described as "Dieses Produkt wird zum extrakorporalen qualitativen Test der Infektion mit neuartiger Coronavirus-Pneumonie (COVID-19) oder des Proteins aus Nucleocapsid des neuartigen Coronavirus (SARS-CoV-2)".

If you have any questions or concerns, please feel free to contact us.

Sincerely Yours,

JOYSBIO (Tianjin) Biotechnology Co., Ltd.



) IL	iste der Antigen-	Tests zum d	rekteri E	Inegemachwek	e des Coronavirus	SARS-CoV-	2.							
					Hordhung bur Anden riardhung – TestiV/* s		trung av	m Ampruch auf beatme	the Territors	gen für den N	achwes	ava Vorie	qens sine	e:
					and a source of		eveninge	n Debrauchsinformation	en.:					
	ne Hirwetse zur vo entests auf SARS-		rigenterto	en Little sowie zu	den der Listung und	pph. auch Sim	khung	van del Latte zugnundele	egenden k	Unternen tinder	Sie auf	arsiener P	Vetrue-te z	R)
ş	v geens Los Attanet v												Elate	si s
	10 TR NACH 10	and selen												
	Matstader Deutscheit Verbeider			Europäischer Deväilmächtigter						Sevelovitat		Sportfall		
5	Name 1.	Stat	к.	Notes	Natur	514-8	4	Handkillsname des Teats	Testor	Arthum.	8	99%L. Ventr Mear	-	95%) Verter Johan
	BELING KEWEI GLINICAL DIAGNOETIC REAGENT INC.	BELING	01	GESING Gram	Lotus NE, R. V.	Chen Haulig	ni.	COVID-13 Antigen Abpid Test KatColoster Gotti		C18G- 20205003	w. 18	91,21- 85,24	100.00	100 100
	Hangphea Clangene Botech Co Litt	Heights	91	GENING GHIM	Sharighai Imternacional Holang Corporation Gintary (Europe)	Handurg	DE	DOVID-19 Antigen Rapid Sear Cassette		1009/5002- 1009/48	91,40	823- 960	100.00	94.J 100
	Hangston Really Tech Co. At	Hangshou	DE	GESING GMBH	Lione Lebenswet Gritter	WOLF	DE	Notel Spranavnuk (SARS-Colu-2) Antgen Mapid Test Casaudte (Swiet)		\$5714/00	95.17	9231- 98.07	98,90	98,88 500
	JOYSERO (Tianjin) Biolechnology 00., 670	Tianjin	CN	GILSING Gram	Loton HE, B, Y.	Пе марие	NL.	SARS-CW-2 Arragen Report Sect Kor (Collockel Gold) / Test Kor Kar meuartiges: Cononconse- Antigen devisable Cost- Methods:	POC Johne Gentti	009-46- 201 81073	95,10	819- 10.4	101.00	94.8 105.5

für Arzneimittel und Medizinprodukte



IMPORTER

GeSino GmbH Riedbergallee 38, 60438 Frankfurt am Main Germany



DISTRIBUTOR

Texolution Medical GmbH & Co. KG

Agentschap Benelux Heringsbosch 7-a 6451 PH Schinveld Nederland

Telefoon:+31 (0)45 - 5255117Mobiel:+31 (0)6 - 83565291E-Mail:bram@texolution-germany.de

