STORZ

KARL STORZ — ENDOSKOPE

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UH 400/400 U/401/401 U Hochfrequenz-Chirurgiegerät AUTOCON® III 400

INSTRUCTION MANUAL

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1 Important information for users of KARL STORZ devices

It is recommended to check the suitability of the products for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual forms part of the device. KARL STORZ SE & Co. KG, referred to in the following simply as KARL STORZ, assumes no liability nor provides any warranty whatsoever for damage and consequential damages that arise due to non-compliance with the instruction manual.

- Read the instruction manual carefully before using this device.
- Store the instruction manual in a safe place throughout the service life of the device.
- Keep the instruction manual accessible to surgical staff.
- Give the instruction manual to each successive owner and/or user of this device.

1.1 Revision index

Software version	Last revised
Valid from version 2.2.1	Version 2.0 – 2017/11

1.2 Validity

This instruction manual Art. no. 96206584D applies only to the devices UH 400, UH 400U, UH 401 and UH 401U.

1.3 Other applicable documents

- ▶ Comply with other applicable documents in the appendix.

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Important information for users of KARL STORZ devices

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Risk type, source and consequences thereof (Personal injury)!

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Risk type, source and consequences thereof (Material damage)!

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Symbol	Risk level	Probability of occurrence	Consequences of non-compliance
	DANGER	Immediate risk	Death, serious injuries
	WARNING	Possible risk	Death, serious injuries
	CAUTION	Possible risk	Minor injuries
	NOTE	Possible risk	Material damage

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Tips and additional information to facilitate tasks

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1.4.4 Other symbols and labeling

Symbol/ labeling	Meaning
☑	Prerequisite for an action
▷	Action with one step
1. 2. 3.	Action with several steps in a binding sequence
☛	Result of preceding action
•	List (first level)
•	List (second level)
Emphasis	Emphasis
..., see section xxx page xxx	Cross reference

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2 Safety

2.1 Intended use

The HF device is intended exclusively for the generation of electrical power for unipolar and bipolar cutting and coagulation on tissue structures in surgical operations.

It is used in the following areas:

- General surgery
- Arthroscopy
- Endoscopy
- Gastroenterology
- Gynecology
- Hand surgery
- ENT
- Cardiac surgery (including open-heart surgery)
- Minimally invasive surgery (MIS)
- Neurosurgery
- Orthopedics
- Plastic surgery and dermatology
- Pneumology
- Thoracic surgery
- Urology, including transurethral resection (TUR)
- Oral and maxillofacial surgery

2.1.1 Patient profile

The use of AUTOCON® III 400 is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

2.2 Contraindications

Not intended for the preparation of metallic implants such as stents.

- Not suitable for surgical procedures which require continuous activation of the HF generator.

The following body parts and tissue types must not be treated:

- Not suitable for working on bone tissue.

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The HF device must not be used if, in the opinion of an experienced physician or according to current specialist literature, such an application would endanger the patient, e.g. due to the patient's general condition, or if other contraindications are present.

-  KARL STORZ requires that the HF device is operated under the supervision of qualified and authorized personnel. The surgeon and medical staff must be trained in the fundamental principles, rules for use and risks of HF surgery and must be familiar with these in order to safely and reliably prevent putting patients, staff and equipment at risk.
-  Contact your KARL STORZ specialist dealer about training and training material.
-  Any other use is neither intended nor proper and must be effectively prevented.

2.3 General safety information

- ▷ Ensure that no electronic devices that are subject to interference from electromagnetic fields are set up in the vicinity of the HF device.
- ▷ Observe the instructions on electromagnetic compatibility provided in section EMC, page 217.

WARNING: Danger of electric shock!

- ▷ Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.
- ▷ Ensure constant access to the power plug at the site of installation. The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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
Additional equipment connected to medical electrical equipment must verifiably comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please contact your local specialist dealer or the Technical Support, see section Technical Support, page 100.

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i To protect personnel, KARL STORZ recommends the use of a smoke evacuator to extract electrosurgical smoke. To this end we recommend the KARL STORZ S-PILOT® UP 501 which, in combination with the AUTOCON® III 400 (plus connecting cable UP 004), guarantees automatic smoke extraction.

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2.4 Personal safety instructions

2.4.1 Ambient conditions

▷ Do not use the HF device in the immediate vicinity of the patient. Observe the minimum distances recommended by KARL STORZ (see Fig. 1).

2.4.2 Patients with pacemakers

Malfunction or destruction of the pacemaker can endanger the life of the patient or result in irreversible injuries to the patient.

▷ In the case of patients with pacemakers, consult the cardiologist before carrying out HF surgery.

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- ▷ Use bipolar HF methods.
- ▷ Attach the HF neutral electrodes close to the operating field.
- ▷ Set the demand pacemaker to a fixed frequency.
- ▷ Ensure that the pacemaker does not come into contact with the HF electrode.
- ▷ Keep a fully operational defibrillator within reach.
- ▷ Carry out a postoperative pacemaker check.

2.4.3 Hazard-free patient positioning

- ▷ Position patients so that they are not touching any metal parts that are grounded or have considerable capacitance relative to ground (e.g. operating table brackets). If necessary, place anti-static cloths between the patient and the underlay.
- ▷ Ensure that the patient does not touch any wet clothes or underlay.
- ▷ Place anti-static cloths between areas with heavy sweating and skin-to-skin contact areas on the patient's torso.
- ▷ Ensure that the patient is resting on a suitable surface in order to prevent pressure necrosis.
- ▷ Drain urine via a catheter.

2.4.4 Correct connection of the HF device

- ▷ Always ground the HF device via the potential equalization connector. Also observe the requirements in section 8.6.7 of IEC 60601-1 regarding medical electrical systems.

i The device features a potential equalization connector which must be connected via a potential equalization line with a connector in the room where it is set up. This additional potential equalization measure means that all conductive parts in the vicinity of the patient are connected electrically to each other and small contact voltages are compensated.

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As such, both transient disturbances from the switching of adjacent devices and the risk to patients from low-frequency electric current are reduced. The connection should be as low-resistance as possible and should not exceed a resistance value of 0.2 Ohm between the central potential equalization of the operating room and the potential equalization connector AUTOCON® III 400 and every other metallic touchable surface. If the design of the installation does not include an equipotential bonding bar, the connection to the PE conductor should be as low-resistance as possible.

- ▷ With simultaneous use of HF surgical devices and physiological monitoring devices on a patient, the monitoring electrodes should be positioned as far away as possible from the surgical electrodes.
- ▷ Needle electrodes for monitoring are not recommended.
- ▷ In all cases monitoring systems are recommended which include mechanisms for limiting the high-frequency current such as protective resistors and HF chokes.
- ▷ Place lines from monitoring devices so that they do not lie on the patient's skin.
- ▷ Keep the leads to the HF electrodes as short as possible and position them so that they do not touch the patient or other leads.
- ▷ Do not place any objects on the HF device.

2.4.5 Correct use of the HF device

Inadvertent activation of the HF device outside the user's field of vision can injure the patient.

- ▷ Activate the HF device only when the electrode is in your field of vision and you can quickly deactivate the HF device at all times.
- ▷ If the HF device is activated inadvertently, switch it off immediately using the Standby button.
- ▷ Take particular care when using a footswitch or manual switch.

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Improper preparation or errors in the HF device can cause damage to the HF device.

- ▷ Use the automatic monitoring functions to ensure that the HF device is working properly. See the section Test for proper functioning, page 26 for information on the automatic test functions.
- ▷ Ensure that no conductive fluids (e.g. blood or amniotic fluid) have penetrated the footswitch or the manual switch.
- ▷ Ensure that the cables for the footswitch and the manual switch are free from short circuits and broken leads.

2.4.6 Configuring HF device settings and using accessories

Setting the output power too high can injure the patient! Therefore, before you increase the output power, ensure that:

- the neutral electrode is attached properly;
- the working electrodes are clean;
- the plug connections are all correct.

A clearly lower output value or functional failure of the HF surgical unit during normal operation can be due to insufficient attachment of the neutral electrode or insufficient contact in the connections. In this case, check attachment of the neutral electrode and the connections before selecting a higher output power.

Setting the HF device correctly

- ▷ To prevent inadvertent (thermal) tissue damage during operations on body parts with small cross sections and in areas with high resistance (bones or joints), use the bipolar method in these areas.
- ▷ Set the level of the acoustic signal that sounds when the electrode is activated so that it is always clearly audible.

Nerve and muscle excitation due to low-frequency currents!

During HF surgical operations (especially when an arc is formed), part of the HF current is converted

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into a low-frequency current. This current can trigger muscle contractions in the patient:

- ▷ To minimize the risk of injury to the patient, set the power and the effect as low as possible.

Correct use of accessories

- ▷ Only use insulated accessories.
- ▷ Check all electrodes for sharp edges and projecting parts before use.
- ▷ Only use electrodes that are free of defects and in good working order.
- ▷ Never place active electrodes on or near the patient. Ensure that the electrodes cannot come into direct or indirect contact (through electrically conductive objects) with the patient.
- ▷ Electrodes which are not currently in use must not come into contact with the patient.
- ▷ Place the instruments in a safe place: sterile, dry, non-conductive, clearly visible. Instruments which have been placed down must not come into contact with the patient, the medical personnel, or flammable materials.
- ▷ Do not remove hot electrodes from the patient's body directly after cutting or coagulation.
- ▷ Ensure that there is sufficient distance between the patient cables and the cables of the HF device.
- ▷ Do not run the patient cables across the patient.

2.5 Product-related safety instructions

Devices manufactured by KARL STORZ are developed in accordance with the current state of technology and generally accepted safety rules. Despite this, using these products can lead to risks to the life and health of the user or third parties and/or damage to the device or other objects.

- ▷ Only use accessories approved by KARL STORZ, see section Accessories/Replacement parts, page 216.

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- ▷ Only use the device if it is free of technical defects and in good working order and only for the intended purpose, always remaining aware of safety requirements and risks and complying with this instruction manual.
- ▷ Have faults that can adversely affect safety (e.g. deviations from the permissible operating conditions) rectified without delay.
- ▷ Wipe down the HF device only with cleaning agents and disinfectants that are approved in the country of use for surface cleaning. Other reprocessing methods are not permitted! (See section Cleaning and disinfection, page 96.)
- ▷ If any fluids penetrate the device, drain them immediately.

Damage to the device can lead to an undesirable increase in output power due to a device malfunction.

Certain devices or accessories can cause danger in lower power settings. For example, the risk of gas embolism in argon assisted coagulation rises if the HF power is insufficient for the fast creation of an impenetrable eschar layer on the target tissue.

2.6 Safe handling (general instructions)

- ▷ Before each use of the device, check to ensure that it is functioning properly and is in good working order and connected properly.
- ▷ Observe the instructions on intended use in conformance with standards, see section Error indications for EASY monitoring, page 93.
- ▷ During use, always observe and comply with the acoustic signals and/or error indications of the HF device, see section Error indications for EASY monitoring, page 93.
- ▷ The device and accessories may be operated and used only by people who have the necessary training, knowledge and experience.

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- ▷ Regularly inspect the accessories, especially the electrode cables, endoscopic accessories and neutral electrodes, for damage to the insulation, proper operation and expiration date.
- ▷ Never place the instruments on the patient or on the devices.
- ▷ Ensure that no instruments are being cleaned when AUTOSTART is activated.
- ▷ Wear suitable gloves during operations.

2.6.1 Operation area: Avoiding ignition and explosions

Sparks are generated when the HF device is used as intended!


- ▷ Do not use the HF device in areas where there is a risk of explosion.
- ▷ Do not use any flammable or explosive liquids.
- ▷ If the display fails, do not use the HF device any longer!
- ▷ During operations (e.g. in the head or thoracic regions), avoid using ignitable anesthetics and gases that support combustion (e.g. nitrous oxide or oxygen) or extract them using a vacuum system.
- ▷ Use exclusively non-flammable cleaning agents, disinfectants and solvents (for adhesives). If you cannot avoid the use of flammable cleaning agents, disinfectants or solvents, ensure that they have fully evaporated before using HF surgical equipment.
- ▷ Ensure that no flammable liquids collect beneath the patient or in body cavities (e.g. the vagina). Suction or flush body cavities before activating the device.
- ▷ Wipe off all liquids before using the HF device.
- ▷ Ensure that no ignitable endogenous gases are present.
- ▷ Ensure that all flammable materials saturated with oxygen (e.g. cotton or gauze) are kept far enough away from the HF environment that they cannot ignite.

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Application point for neutral electrode

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2.6.2 Applying the neutral electrode

i Observe the instructions on the use of the neutral electrode in the instruction manual and the information on the package of the neutral electrode.

In the unipolar HF method, the neutral electrode feeds the current introduced into the patient's body at the surgical site back to the HF device.

▷ To prevent a rise in temperature at the current exit point, the following conditions must be ensured:

- Sufficiently large contact surface between the neutral electrode and the patient's body;
- High electrical conductivity between the neutral electrode and the patient's body.

▷ To prevent the patient being burned by the neutral electrode, you must comply with the following conditions:

- Select the application point for the neutral electrode so that the current paths between the active and neutral electrodes are as short as possible and run longitudinally or diagonally to the patient's body (because muscles are more conductive in the direction of the fibrils, see Fig. 2).
- For surgery in the thoracic region, do not run the current path transversely across the patient's body and ensure that the patient's heart is never in the current path.
- Depending on the surgical site, apply the neutral electrode to the nearest upper arm or thigh if possible, but never closer than 20 cm.
- In the case of self-adhesive disposable electrodes, comply with any further manufacturer instructions regarding the point of application.
- Ensure that the application area is free of scar tissue, bony protuberances, surface hair and ECG electrodes.

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- Ensure that there are no implants (e.g. bone nails, bone plates or endoprotheses) in the current path.
- Ensure that no short circuits can occur at the neutral electrode connection.
- Avoid areas where fluids may collect.
- As far as possible, use split neutral electrodes with a sufficiently large surface area (patient age and max. output power during operation have to be considered).

Before applying the neutral electrode

- ▷ Remove excess body hair.
- ▷ Clean the application site, but do not use any alcohol, since it dries out the skin and increases the transition resistance.
- ▷ If the patient has poor circulation, massage or brush the application site.
- ▷ Attach the neutral electrode over the entire contact surface evenly. Secure reusable neutral electrodes with rubber bands or elastic straps so that they do not loosen when the patient moves. Ensure that the patient's circulation is not impaired (risk of necrosis).
- ▷ Never use wet cloths or conductive pastes.
- ▷ Ensure that no liquids (e.g. cleaning fluids, disinfectants, blood or urine) penetrate between the patient and the neutral electrode.
- ▷ Do not place the neutral electrode under the patient's buttocks or back.
- ▷ Ensure that there are no ECG electrodes in the current path of the HF device.
- ▷ Check the neutral electrode for damage and to ensure that it is working properly.
- ▷ Replace defective accessories immediately.

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Example application using a disposable electrode

- ▷ Remove the protective film and attach the self-adhesive disposable electrode to the patient. Ensure that the long edge of the disposable electrode faces the operation site and the electrode is fully in contact with the skin. This avoids excessive current concentration on the short edge.
- ▷ Using both hands, press the self-adhesive disposable electrode firmly against the patient's skin.
- ▷ Clamp the electrode tab to the electrode cable.
- ▷ After the operation, remove the disposable electrode carefully to avoid skin damage.

With a one-piece neutral electrode

- ▷ Check the one-piece neutral electrode during surgery.
- ▷ Ensure that the one-piece electrode is not blocked at the device.

With a split neutral electrode

- ▷ Apply the split neutral electrode correctly and without any additional objects, as the HF device does not recognize the bridging of the split surfaces with other objects.
- ▷ Ensure that the HF current flows equally to both parts of the split neutral electrode.

i For monitoring the neutral electrode, see section EASY neutral electrode monitoring (EASY monitoring), page 28.

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3 Descripción

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3.1 Display and control elements

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3.1.1 Control elements on the front panel

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

i While activating an instrument, the activation bar (7 - 10) of the corresponding socket illuminates yellow or blue.

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- ① Standby button (ON – surrounded by a white light)
- ② Symbol 'Standby button'
- ③ Neutral electrode isolated from ground for HF
- ④ Symbol 'CF type applied part with defibrillation protection'
- ⑤ Symbol 'Observe instruction manual'
- ⑥ Touch screen with mode selection buttons
- ⑦ Activation bar upper unipolar socket
- ⑧ Activation bar lower unipolar socket
- ⑨ Activation bar upper bipolar socket
- ⑩ Activation bar lower bipolar socket

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i Disconnect patient-related sockets (11, 12, 13, 17, and 18) from the patient before commencing service and maintenance work, see section 7 Reprocessing and section 8 Maintenance and repair.

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3.1.2 Unipolar socket module (left)

- 11 Connection socket for unipolar instruments with hand or footswitch* (activated state – white surrounding light)
- 12 Connection socket for unipolar instruments with hand or footswitch* (activated state – white surrounding light)
- 13 Connection socket for neutral electrode* (activated state – white surrounding light)

*Applied part type CF according to IEC 60601-1, with defibrillation protection

Unipolar connection socket 11 (see Fig. 3)

- Connection socket 11 'top':
- 14 3-pin US type
- 15 5 mm KARL STORZ/Erbe VIO
- 16 4 mm socket (foot-switched)

Unipolar connection socket 12 (see Fig. 4)

- Connection socket 12 'bottom':
- 14 3-pin US type
- 15 BOVIE (foot-switched)
- 16 4 mm socket (foot-switched)

Connection socket for neutral electrode (see Fig. 5)

3.1.3 Bipolar socket module (right)

- 17 Connection socket for bipolar instruments with footswitch, fingerswitch or AUTOSTART* (activated state – white surrounding light)
- 18 Connection socket for bipolar instruments with footswitch, fingerswitch or AUTOSTART* (activated state – white surrounding light)

Bipolar connection sockets

- Upper bipolar socket 17 (see Fig. 6)
- 19 2-pin US type (28.58 mm)
- 20 1x KARL STORZ/Erbe VIO
- Lower bipolar socket 18 (see Fig. 7)
- 19 2-pin US type (28.58 mm)
- 20 2x KARL STORZ/Erbe VIO

* Applied part type CF according to IEC 60601-1, with defibrillation protection

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- ②① Connection socket 1 for footswitch
- ②② Connection socket 2 for footswitch
- ②③ Connection for potential equalization
- ②④ Power supply for IEC connector
- ②⑤ Fiber-optic signal input socket
- ②⑥ Fiber-optic signal output socket

Use the following connections only for service and training purposes:

- ②⑦ Ethernet connector (only for connection to KARL STORZ OR1™ control NEO)
- ②⑧ USB connector
- ②⑨ Connector S-PILOT® UP 501 with UP 004 S-PILOT® connecting cable, dia. 3.5 mm, 300 cm (max. 25 V/1 A)
- ③① UART communication interface
- ③② Power switch
- ③③ Line fuses (see page 102)

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-  The Ethernet interface is only suitable for connection of the OR1™ control NEO for a dedicated network.
-  Only use the USB interface as described in this manual. To this end, only use the KARL STORZ USB stick 20 0402 82.

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




Symbol	Designation
	Footswitch connection
	Neutral electrode isolated from ground for HF
	CF type applied part with defibrillation protection (cardiac floating)
	Alternating current
	Standby
	During activation (of the HF device) HF energy in the radio frequency range 9 kHz to 400 GHz is applied, which produces electromagnetic radiation. The device uses close-range radio frequency identification technology at the HF connection lines (RFID).
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE), see Disposal.
	Labeling of the (active) HF output; Caution: Hazardous voltage.

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


Symbol	Designation
	Manufacturer
	Date of manufacture
	Follow instructions for use
	Potential equalization connector
	Fiber-optic signal input
	Fiber-optic signal output
	Ethernet connector (connection to KARL STORZ OR1™ control NEO system).
	USB connector*
	S-PILOT® connector* (surgical smoke extraction)
	UART communication interface for service purposes*
	Reference number
	Serial number
	The device satisfies the Medical Device Directive (MDD) 93/42/EEC and 2007/47/EC (NB 0123) as well as the R & TTE Directive 1999/5/EC.
	Consult instructions for use
	General warning of physiological effect: The device generates physiological effects which may endanger the patient and/or user.
	'Rx only'

* In the event of queries contact Customer Service (see section 2.3).

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see Fig. 8: AUTOCON® III 400 UH 400.

3.2.2 Symbols on the packaging

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no.96216316DF. It can be downloaded from www.karlstorz.com.

3.3 Scope of supply*

- 1 AUTOCON® III 400 UH 400 or UH 401 (UH 401 incl. option BiVascularSafe)
 - 1 Power cord 400 A
 - 1 Instruction manual with inspection protocol
- * also see page 106 (specifications).

3.4 Components required for operation

- Power cord
- Footswitch
- Neutral electrode for unipolar applications
- Connecting cable for neutral electrode or instrument
- Instrument (unipolar or bipolar)

3.5 Operating conditions

Temperature: +10°C to +40°C
 Relative humidity: 30 to 75%, non-condensing
 Air pressure: 700 to 1 060 hPa
 Operating altitude (max.): 3000 m above sea level

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4 Preparation

4.1 Setting up the HF device

NOTE: HF equipment generates high-frequency electromagnetic fields that can interfere with the operation of sensitive electronic equipment and, for example, lead to image distortion. HF equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference. Parallel cable routing with cables from other devices should be avoided as a general rule, as HF cables emit HF emissions when activated, which may interfere with other equipment.

Additional measures include:

- Connect the power connection of the AUTOCON® III 400 or the impaired device to another circuit.
- Increase the distance between the AUTOCON® III 400 and the impaired device.
- Lay instrument cables so that, as far as possible, they are not in the vicinity of other devices and their connecting cables.

WARNING: Risk of electric shock! Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.

DANGER: Risk of burns to patients due to excessive leakage current! Locate the HF device outside the immediate vicinity of the patient, see section Ambient conditions, page 6.

HF devices may be used only in rooms used for medical purposes that meet the requirements of DIN VDE 0100-710 or IEC 60364-7-710.

If the HF device was previously stored or transported at temperatures below +10°C or a relative humidity above 75%, non-condensing, it will take approximately three hours to adjust to room temperature.

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1. Observe the specified operating conditions, see section Operating conditions, page 21.
2. Place the HF device on one of the following platforms:
 - a table,
 - an equipment cart,
 - a console suspended from a ceiling support or wall-mounted brackets.
3. Place the HF device a sufficient distance away from other electronic equipment, see section EMC, page 217.
4. Position the HF device with the front of the device facing the patient and surgeon.
5. Do not place any other device on the HF device.
6. Do not place any other objects on or above the HF device.
7. If necessary, place the HF device on the Argon beamer.
8. Connect the potential equalization cable to the potential equalization connector. The device's ground line should be installed by a qualified electrician.
9. Connect the power cord.

4.2 Switching on the HF device

i Do not use the HF device if the display components are not working! For troubleshooting see section Recognizing and remedying errors, page 79.

1. Switch the HF device on using the power switch on the rear side of the device and then touch the permanently lit Standby button on the front panel.
 - ☛ The HF device performs a self-test: All display and control elements light up.
2. Check all display and control elements for proper operation:
 - Power switch
 - Touch screen
 - Unipolar connection sockets
 - Bipolar connection sockets
 - Activation bar for unipolar and bipolar sockets

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- ☛ The main screen appears, and the HF device is ready for use.
- ☛ The parameters of the most recently selected program appear on the display.
- i** If the device is switched off at the front panel, the Standby mode is activated.

4.3 Connecting instruments

▷ Before connecting instruments, ensure that the following conditions are met:

- Combinations of accessories not mentioned in the instruction manual may be used only if they are explicitly designed for the intended use. Always observe performance characteristics and safety requirements.
- The insulation of the accessories (e.g. HF cables and instruments) must be sufficient for the maximum peak output voltage (see IEC 60601-2-2 and IEC 60601-2-18).
- The HF device is designed for connection of HF cables measuring 3 m or 4 m in length.
- Do not use accessories with defective insulation.

4.3.1 Instruments for unipolar application

1. Plug the neutral electrode cable into the socket for the neutral electrode and choose the corresponding neutral electrode type, see chapter Selecting the neutral electrode, page 36.
 - ☛ The socket illumination goes dark.
 - ☛ The neutral electrode button changes from gray to the measurement color (green, yellow or red).
2. Connect the electrode handpiece to one of the two unipolar connection sockets.
 - or –
 - In the case of an accessory without a finger button: Connect the footswitch and unipolar connecting cable.
 - or –
 - Connect the unipolar cable to one of the two unipolar connection sockets for unipolar instruments.

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4.3.2 Instruments for bipolar application

1. Connect the bipolar cable to the instrument (e.g. forceps).
 2. Connect the bipolar cable to one of the two bipolar connection sockets.
 3. For bipolar use without AUTOSTART, connect the bipolar switch.
– or –
Select a mode with AUTOSTART function at the corresponding socket.
- ☛ Once the instrument is connected, the application starts after the configured reaction time.

i When introducing bipolar instruments into a metal trocar, the AUTOSTART function should be disabled beforehand!

4.3.3 Connecting the footswitch

In addition to the manual switch, a footswitch can be used to activate various operating modes.

- ▷ Connect the desired footswitch only during operation to one of the two connection sockets for footswitches.
- ☛ The HF device automatically detects the connected footswitch and indicates this on the front panel display, including the selected connection socket.

i One dual-pedal footswitch and one single-pedal footswitch can be connected. Footswitches without a changeover switch cannot be used (see table).
During the operation, in addition to the footswitch, only the connection to the Argon beamer with fiber-optic cables and, if relevant, the connection to the KARL STORZ OR1™ control NEO system may be plugged into the rear of the device.

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The following footswitches can be connected to the HF device:

Art. no.	Designation
UF 901	Single-pedal footswitch with changeover function switch
UF 902	Dual-pedal footswitch with changeover function switch

4.4 Test for proper functioning

4.4.1 Automatic test function

The HF device automatically performs cyclic testing during operation. If any errors occur, see section Recognizing and remedying errors, page 79.

4.4.2 Performing the test for proper functioning

Perform the following functional test before putting the device into service:

i The accessories must be designed for the specified maximum voltage.

- Connect the neutral electrode and attach it securely to the patient's arm.
 - The EASY neutral electrode indicator changes to green.
- Remove the neutral electrode.
 - The indicator changes to red, acoustic signals sound.
- i** The neutral electrode used for this test may not later be used for an operation.
- Connect a unipolar HF handpiece to a unipolar socket if there is a green EASY indicator and use the manual switch and footswitch to individually activate 'Cut' and 'Coag'.
- Check the settings on the display.
- Now change to the bipolar output and connect bipolar forceps.
- Select a mode with AUTOSTART, grasp moist gauze with the forceps, and check the display.

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7. Now change to a mode without AUTOSTART and use the footswitch to activate the bipolar output. Check the settings and indicators in the bipolar section.

4.4.3 Actions in case of problems

Proceed as follows in case of functional problems:

1. Immediately disconnect the patient from the HF device.
2. Perform a technical inspection of the HF device.
3. Observe the relevant national regulations. For example, in Germany report incidents and near-accidents to the German Federal Institute for Drugs and Medical Devices in accordance with Section 3 of the German Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV). Observe the provisions of the in-house reporting system in this regard.
4. Consult the Technical Support, see section Technical Support, page. 100.

i The HF device can be switched off at any time using the power switch  on the rear of the device.

4.5 Neutral electrode monitoring

i Always use the largest possible electrode when attaching a neutral electrode.

4.5.1 General information

i KARL STORZ recommends using split neutral electrodes, since only this type of electrode allows the HF device to detect detachment of the neutral electrode if this occurs.

Monitoring of the neutral electrode minimizes the risk of burns at the site where the neutral electrode is attached.

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Two types of neutral electrodes can be monitored:

- Split electrodes for infants (for use with reduced power)
- Split neutral electrodes

The type of neutral electrode and its contact quality are selected and/or shown in the Neutral electrodes menu.

Errors related to the neutral electrode and possible remedies are shown on the display, see section Recognizing and remedying errors, page 79 (see Fig. 9, neutral electrode modes).

4.5.2 EASY neutral electrode monitoring (EASY monitoring)

- i** The maximum output power of unipolar current types is reduced to 50 W when a children's electrode is selected.

The EASY monitoring function measures changes in the resistance between the patient and the high frequency surgical unit before and during HF activation. If necessary, it generates visual and audible alarms to request staff intervention. This requires using a split neutral electrode with appropriate contact areas and suitable transition resistance, attached to the patient according to the manufacturer's instructions. The EASY system does not monitor the currents through the individual contact surfaces of the split neutral electrode.

A KARL STORZ electrode with a surface area of at least 90 cm² must be used for the unipolar program 'Resection' and the 'Moderate Coagulation' mode.

If an error message is generated, the display changes from green via yellow to red, depending on the type of error.

- !** **NOTE:** Incorrect application of the neutral electrode!
Ensure compliance with the specifications for correct attachment of the neutral electrode with regard to size, adhesive properties and full-surface contact of the complete electrode.

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5 Operation

5.1 Connecting the power supply

- ☑ The line voltage must match the voltage specified on the type plate.
 - ☑ The potential equalization cable must be connected to the potential equalization connector (see section 4.1). The device's ground line should be installed by a qualified electrician.
1. Connect the power cord to the generator and the plug to the grounded wall socket.
 2. Use the power switch 31 on the rear to switch the device on (see Fig. 10).
 - ☛ The Standby button on the front panel is permanently illuminated.
 3. Switch the device on at the Standby button.
 - ☛ The device performs a self-test.
 - ☛ The bars of the activation displays light up as do the white socket surrounds 11, 12, 13, 17 and 18 (see page 17).
 - ☛ Full functionality of the loudspeaker is indicated by the start melody.
 - ☛ The individually configurable start screen appears if it has been set up.

5.2 Program overview

5.2.1 Display

The status bar is located at the top of the display (see Fig. 11). The settings of the four sockets are shown below the status bar.

Settings can be configured for each of the sockets.

1. The 'Effect' button is used to set the effect of electrosurgical cutting or coagulation.
2. Assign the footswitch using the 'Pedal' button.
3. The 'Mode' button allows the desired type of current to be selected.
4. Click on 'max. Watt' to adjust the maximum output power.

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

5.2.2 Status bar

The program name and the buttons 'Favorite', 'Key lock', 'EASY', 'Information' and 'Menu' are arranged in the status bar (see Fig. 12).

In combination with an Argon beamer and selection of an Argon mode, the 'Argon' button is also additionally displayed (see Fig. 13).

5.3 Activating and deactivating sockets

▷ To activate a deactivated socket, plug a connecting cable into the connection socket.

– or –

Press the button 'activate socket' (see Fig. 14).

☛ An overview of the socket settings appears.

The overview is grayed out if no instrument is connected to the socket (see Fig. 15).

The socket illumination extinguishes and the socket selection window lights up when an instrument is plugged in.

▷ To hide sockets, click on the 'off' button next to the socket setting overview (see Fig. 16).

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5.4 Unlocking the screen

The device screen locks automatically. To unlock it, click on any control element and push the slider which then appears from left to right. An open lock then appears in the status bar (see Fig. 17, Unlocking the screen).

To deactivate the key lock or to change the time up to automatic screen lock, see section Menu 'System settings', page 67.

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
5.5 Configuring output currents

- i** All selection windows are closed after 30 seconds without assumption of the modification.
- i** If selection windows are open, all buttons outside this field are inactive. In this case, activations are not possible.
- i** A change to the currently loaded program, e.g. by adjusting the output, is indicated by the information 'changed' below the program name.

5.5.1 Selecting the mode


1. To select the type of current for unipolar cutting, first select one of the two sockets on the left side (see Fig. 18).
2. Press the yellow icon under the 'Mode' button.
 - A selection screen appears for the available modes, and the rim of the associated socket starts flashing.
3. Select the desired mode by pressing the corresponding button.
 - or -
 - Deactivate the mode by pressing the button 'Cutting off (X)'.
4. Press the '?' button for more information on this selection.
5. For additional selection options press ◀/▶ inside the selection window.
6. Confirm your selection by pressing the 'OK' button.
 - The main screen will be displayed.
 - or -
 - Press the 'Back' button to return to the main screen without changing the selection.
- i** If a mode is changed within a socket, the set parameters, e.g. effect and max. Watt, remain the same for the respective mode. However, if, for example, the factory default setting of a mode is adapted and subsequently changed to a different mode and then back again, the user changes are not undone.

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
20

5.5.2 Specifying power limits

1. To select the maximum output for unipolar cutting, first select one of the two unipolar sockets and press the yellow icon under the 'max. Watt' button (see Fig. 19).
2. Use the '+' and '-' buttons to adjust the power level in single steps.
- or -
Use the slider to set the power level in steps of 10.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
- or -
Press the 'Back' button to return to the main screen without changing the selection.

5.5.3 Selecting the effect

1. To select the effect for unipolar cutting, first select one of the two unipolar sockets and press the yellow icon under the 'Effect' button (see Fig. 20).
2. Use the '+' and '-' buttons to adjust the effect in individual steps.
- or -
Use the slider to set the effect.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
- or -
Press the 'Back' button to return to the main screen without changing the selection.

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22

5.5.4 Assigning the footswitch

i Handpieces and instruments with manual switches can be activated without a configuration setting.

A dual-pedal footswitch and/or single-pedal footswitch, each with a changeover switch, can be connected.

The changeover switch enables switching between pedal levels.

1. Call up the footswitch menu by pressing the button 'Pedal' (see Fig. 21).
2. Select the desired footswitch by pressing the corresponding button. For example, choose the active pedal level for cutting and coagulation for the upper left socket.
- or -
Deactivate the footswitch by pressing the button marked with an 'x'.

☛ The edge of the selected button lights up green (see arrow in Fig. 21).

3. Confirm the selection by pressing the 'OK' button.
- or -
Press the 'Back' button to return to the main screen without changing the selection.

☛ The socket is assigned to the active pedal level.

4. Pedal levels can be changed using the footswitches. Press the black changeover switch to change the socket.

☛ The orange background of the changeover switch in the 'Pedal' button indicates that the lower left-hand socket is activated (see Fig. 22).

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i If two footswitches are connected, either a single-pedal footswitch or a dual-pedal footswitch can be selected for coagulation (see Fig. 23).

The following pedal icons can be distinguished between:

Icon/ Button	Description	Icon/ Button	Description
	Dual-pedal footswitch CUT active		Dual-pedal footswitch COAG active
	Dual-pedal footswitch CUT inactive		Dual-pedal footswitch COAG inactive
	Dual-pedal footswitch not connected		Dual-pedal footswitch not connected
	Single-pedal footswitch COAG active		Single-pedal footswitch not connected
	Single-pedal footswitch COAG inactive		Deactivate footswitch CUT
	Deactivate footswitch coagulation		

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i The maximum power output of unipolar current types is reduced to 50 W when a children's electrode is selected.

2. Select the type of connected neutral electrode by pressing the corresponding symbol.
EASY: for monitoring split neutral electrodes (see Fig. 24)
BABY: for monitoring split neutral electrodes for infants (see Fig. 25)
MONO: to select a one-piece neutral electrode (see Fig. 26)

3. Press the '?' button for more information on this selection.

4. Confirm your selection by pressing the 'OK' button.
- or -
Press the 'Back' button to return to the main screen without changing the selection.

☛ The selected type of neutral electrode in connection with a color indicator for the contact quality is shown in the status bar.

i When using the 'EASY' and 'BABY' mode, no unsplit electrodes are accepted. Using the 'MONO' mode, no split electrodes are accepted.

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According to the contact quality, the icons are shown in color for neutral electrodes:

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Icon/Button	Description
	Split neutral electrode contact quality OK
	Split neutral electrode contact quality not optimum
	Split neutral electrode contact quality insufficient
	Split neutral electrode not connected
	Split baby neutral electrode contact quality OK
	Split baby neutral electrode contact quality not optimum
	Split baby neutral electrode contact quality insufficient
	Split baby neutral electrode not connected
	Non split neutral electrode contact quality OK
	Non split neutral electrode not detected or contact quality insufficient
	Non split neutral electrode not connected
	Display contact quality

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5.5.6 Coding system

The coding system enables automatic instrument identification. The coding system recognizes the connected coded KARL STORZ instrument and automatically selects the preferred parameters.

1. Insert the coded instrument in an AUTOCON® III 400 socket.

☛ The instrument data is read (see Fig. 27).

☛ A description of the instrument appears:

- Instrument name
- Recognized socket
- Article number
- Batch number
- Remaining usages

☛ The parameters are transferred automatically after 5 seconds and shown in the main screen.

• If the coded instrument is connected to a socket without preset parameters, with the coded system the ideal settings are loaded for the coded KARL STORZ instrument.

• If the coded instrument is connected to a socket with preset parameters, a plausibility check is performed. If the preset values for the coded instrument are within the permitted range, these are not overwritten. The coded instrument can now be used with the preset parameters.

- or -

The selection is saved by pressing the 'OK' button.

To return to the main screen without changing the selection, press the 'back' button.

☛ The coded instrument can now be used.

☛ The permitted parameters for the coded KARL STORZ instrument can still be selected, all other modes are gray.

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5.6 Mode overview

An overview of the current types that can be executed with the HF device is shown below.

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
Cutting mode symbol	Designation
	Gastro Knife 2
	Gastro Knife 3
	Argon*

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Coagulation mode symbol	Designation
	Laparoscopy
	Moderate
	Forced coag
	Resection
	Spray
	Forced mixed
	Forced cutting
	Gastro Coag
	Argon flexible*

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Coagulation mode symbol	Designation
	Argon flex. pulse*
	Argon open*
	Cardiac Thorax
	Cardiac Mammary
	SimCoag

i * These modes can only be activated in combination with an Argon beamer. Our service team will be happy to provide you with information on the available combination options.

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5.6.2 Bipolar modes


Cutting mode symbol	Designation
	Bip. cutting
	Bipolar scissors
	Bip. resection
	Bip. vaporization

Coagulation mode symbol	Designation
	RoBi®
	Laparoscopy
	Standard
	Bip. resection
	Bip. vaporization
	BiVascularSafe*
	BiVascularSafe Open Surgery*

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


Coagulation mode symbol	Designation
	Standard AUTO
	Bipolar scissors
	Micro
	Forced

i * These modes are available for the device versions UH 401 and UH 401U.

The information and data regarding settings, application points, application duration and instrument use are based on clinical practice. However, these are only basic guidelines which must be tested for suitability by the operator. Depending on individual conditions, it may be necessary to deviate from the provided data. Medical practice is continuously evolving as a result of R&D and clinical experience. This may also make deviations from the provided data necessary.

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5.7 Unipolar cutting modes

5.7.1 Laparoscopy

This mode is used in laparoscopy and arthroscopy for unipolar cutting.

Application areas

Laparoscopy, arthroscopy

Suitable instruments

- Arthroscopy electrodes
- Laparoscopy electrodes

5.7.2 Standard

In this mode a high-performance HF current with a low crest factor is used for cutting biological tissue.

ARC Control adjusts the power output to the minimum required level in response to variations in tissue type and changes in the cutting area or speed.

Application areas

Cutting tissue with low electrical resistance, such as muscle tissue or vascular tissue.

Cutting or preparing fine structures

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- Resektoskop (unipolar)
- Resektionsschlinge
- Vaporisations-Elektrode

5.7.3 Micro

This mode is used for electro-surgical cutting using micro-electrodes. It enables extremely fine control of the power level and precise work.

Application areas

Pediatric surgery, neurosurgery, plastic surgery

Suitable instruments

- Micro needle electrodes

5.7.4 Resection

This mode is used in gynecology and urology. ARC Control generates the cutting effect with simultaneously minimized output power. ARC Control facilitates direct cutting and prevents electrode adhesion.

- i** Use non-conductive irrigation fluids.

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- Gynäkologische Schlingen

5.7.5 Resection C-Cut®

This mode is used for electro-surgical cutting using unipolar cutting instruments and electrodes.

The cut features a particularly pronounced level of hemostasis.

Different hemostasis effects can be set.

i Use non-conductive irrigation fluids.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), vaporization of prostate tissue (TUR-VAP).

Suitable instruments

- Resectoscope (unipolar)
- Resection loop

5.7.6 SupraLOOP

This mode is used in gynecology for laparoscopic removal of the uterus. Removal can be achieved by emitting a unipolar cutting current and, at the same time, tensioning the loop.

Application areas

Gynecological, laparoscopic removal of the uterus.

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- Polypektomieschlingen

5.7.7 Dry

This mode is used for unipolar dry cutting. A large, controlled arc is generated, which allows significantly deeper coagulation to be obtained.


Application areas

Cardiac surgery and blood coagulation in retracting blood vessels in the sternum region.

Suitable instruments

- Knife electrodes

5.7.8 Gastro Loop 1

This mode is used in gastroenterology. Polypectomy snares are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a relatively slow pulse rate of 1 cutting pulse per second, this mode is suitable for especially cautious work.

Application areas

Polyp removal using polypectomy snares and flexible endoscopy.

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

- Polypektomieschlingen

5.7.9 Gastro Loop 2

This mode is used in gastroenterology.

Polypectomy snares are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic pulse rate of 1.5 cutting pulses per second this mode is suited to experienced users.

Application areas

Polyp removal using polypectomy snares and flexible endoscopy, with dynamic pulse rate for experienced users.

Suitable instruments

- Polypectomy snares

5.7.10 Gastro Loop 3

This mode is used in gastroenterology.

Polypectomy snares are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic and fast pulse rate of 2.2 cutting pulses per second, this mode is suitable for advanced users.

Application areas

Polyp removal using polypectomy snares and flexible endoscopy, with dynamic fast pulse rate for advanced users.

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
5.7.11 Gastro Knife 1

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a relatively slow pulse rate of 1.3 cutting pulses per second, this mode is suitable for especially cautious work.

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; slow pulse rate for cautious work.

Suitable instruments

- Papillotomes
- Needle knives

5.7.12 Gastro Knife 2

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic pulse rate of 1.8 cutting pulses per second, this mode is suitable for experienced users.

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- Papillotome
- Nadelmesser

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; dynamic pulse rate for experienced users.

Suitable instruments

- Papillotomes
- Needle knives

5.7.13 Gastro Knife 3

This mode is used in gastroenterology. Instruments for papillotomy and endoscopic resections are used for cutting and coagulation. ARC Control generates the cutting effect with simultaneously minimized output power. This mode consists of a series of cutting current pulses followed by a coagulation phase. With a dynamic and fast pulse rate of 2.2 cutting pulses per second, this mode is suitable for advanced users.

Application areas

Papilla incision using a papillotome and flexible endoscopy, resection with needle knives; dynamic fast pulse rate for advanced users.

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

5.7.14 Argon

This mode is used to perform open surgical interventions in combination* with an Argon beamer for argon-assisted cutting. With suitable instruments connected, argon-assisted cutting can be performed using rigid electrodes.

Application areas

Visceral surgery

Suitable instruments

- Rigid argon electrodes
- Argon handpiece

i * Our service team will be happy to provide you with information on the available combination options.

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5.8 Unipolar coagulation modes

5.8.1 Laparoscopy

This mode is used in laparoscopy and arthroscopy for unipolar coagulation.

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- Kugelelektroden
- Messerelektroden
- Spatelektroden

5.8.2 Moderate

This mode is used with contact coagulation to stop hemorrhagic oozing, for hemostasis of relatively large tissue areas, and for small-area coagulation. Tissue carbonization is prevented and electrode adhesion to the tissue is strongly reduced. Greater coagulation depth is achieved in this mode than in other coagulation modes. The degree of surface scabbing can be controlled by adjusting the 'Effect' setting in the range of 1 to 3.

Application areas

Coagulation with high penetration depth; low electrode adhesion to tissue.

Suitable instruments

- Electrodes with large contact areas, such as ball electrodes

5.8.3 Forced coag

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves a high degree of coagulation with low cutting tendency.


Application areas

Fast coagulation with small penetration depth

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- Messerelektroden
- Spatelektroden
- Nadelelektroden

5.8.6 Forced mixed

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves a high degree of coagulation with moderate cutting tendency.

Application areas

Fast coagulation with small penetration depth and moderate cutting tendency.

Suitable instruments

- Knife electrodes
- Spatula electrodes
- Insulated unipolar forceps

5.8.7 Forced cutting

This mode is used for contact coagulation with low tissue penetration, preferably using fine electrodes and electrodes with small contact areas. It achieves good hemostasis with very good cutting tendency.

Application areas

Fast coagulation with small penetration depth and very good cutting tendency.

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
5.8.10 Argon flex. pulse

This mode is used for argon-assisted electrosurgery in combination* with an Argon beamer.

This is the current type Spray. For argon-assisted coagulation, flexible probes are used in combination with endoscopes. The pulse frequency changes with the effect setting. The higher the effect level, the faster the pulse sequence.

Application areas

Gastroenterology, homogeneous surface coagulation

Suitable instruments

- Flexible argon probes

5.8.11 Argon open

This mode is used to perform open surgical interventions in combination* with an Argon beamer for argon-assisted electrocoagulation.

This is the current type Spray. With suitable instruments connected, argon-assisted coagulation can be performed using rigid electrodes.

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- Messerelektroden

5.8.12 Cardiac Thorax

This mode is used in thoracic surgery. It produces forced coagulation.


Application areas

Thoracic surgery

Suitable instruments

- Knife electrodes

5.8.13 Cardiac Mammary

This mode is used in mammary surgery and cardiac surgery. It produces forced coagulation.

Application areas

Mammary surgery and cardiac surgery

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
- Kugelelektroden
- Messerelektroden
- Spatelektroden

5.8.14 SimCoag

This mode is used for simultaneous preparation. Two unipolar sockets can be activated at the same time to enable the simultaneous use of two manually switched instruments. Both handpieces can be switched on and off independently.

This mode is used with non-contact surface coagulation using an arc.

i The output power setting applies to both output sockets, and the power is distributed according to the tissue structure.

! WARNUNG: The output power of all active electrodes can change during use.

Application areas

Simultaneous coagulation and preparation, e.g. for cardiac or mammary surgery

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5.9.3 Bip. resection

This bipolar mode is used in gynecology and in urology for resection with loop electrodes under conductive irrigant solution (saline). ARC Control generates the cutting effect with simultaneously minimized output power. ARC Control facilitates direct cutting and prevents electrode adhesion.

i Make sure that NaCl is used as an irrigation medium.

Secure a continuous irrigation during the application.

Always use conductive lubricants to avoid damages of the urethra.

Avoid continuous activations.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT).

Suitable instruments

- Resectoscope (bipolar)
- Resection loop

i Optimum results are possible only when using the coded KARL STORZ resection cable.

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5.9.4 Bip. vaporization

This bipolar mode is used in gynecology and urology with vaporization. Upon tissue contact the arc is ignited directly resulting in rapid tissue vaporization with minimal heat emission to the surrounding area.

- i** Make sure that NaCl is used as an irrigation medium. Secure a continuous irrigation during the application. Always use conductive lubricants to avoid damages of the urethra. Avoid continuous activations.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT), vaporization of prostate tissue (TUR-VAP).

Suitable instruments

- Resectoscope (bipolar)
- Vaporization electrode

Bipolar scissors should only be operated with the current type bipolar scissors for cutting or coagulation.


5.10 Bipolar coagulation modes

5.10.1 RoBi®

This mode designed specially for RoBi® is used for coagulation in combination with bipolar laparoscopic instruments.

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- i** Make sure that NaCl is used as an irrigation medium.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT)

Suitable instruments

- Resectoscope
- Resection loop
- Vaporization electrode

- i** Make sure that the instrument has contact with the tissue while activating bipolar coagulation to avoid an unintended heating of the irrigation fluid.

5.10.5 Bip. vaporization

This mode is used for bipolar hemostasis in gynecology and in urology with vaporization.

- i** Make sure that NaCl is used as an irrigation medium.

Application areas

Hysteroscopy, transurethral prostate resection (TUR-P), surgical treatment of bladder tumors (TUR-BT), vaporization of prostate tissue (TUR-VAP)

Suitable instruments

- Resectoscope
- Vaporization electrode

- i** Make sure that the instrument has contact with the tissue while activating bipolar coagulation to avoid an unintended heating of the irrigation fluid.

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**5.10.6 BiVascularSafe
(only UH 401/UH 401U)**

This mode is used for the permanent sealing of veins, arteries and tissue bundles. Conventional instruments cannot be used in this mode. All parameters are preconfigured and automatically regulated.

Application areas

Vessel sealing open and laparoscopic


Suitable instruments

- Sealing instruments*

Suitable electrodes*

i *Please contact your KARL STORZ authorized dealer for further advice and information.

**5.10.7 BiVascularSafe Open Surgery
(only UH 401/UH 401U)**

This mode is used for the permanent sealing of veins, arteries and tissue bundles with open surgical sealing instruments. Conventional instruments cannot be used in this mode. All parameters are preconfigured and automatically regulated.

Application areas

Vessel sealing in open surgery

Suitable instrument

- Open surgery sealing instruments*

i *Please contact your KARL STORZ authorized dealer for further advice and information.

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5.10.8 Standard AUTO

This mode is used for arcless contact coagulation with forceps. Activation starts automatically on contact with tissue. The adjustable delay time can be set under MENU – SYSTEM SETTINGS – AUTOSTART DELAY (page 67).

- i** Setting the AUTOSTART mode can result in unintentional coagulations, e.g. when bipolar forceps are used for gripping while the AUTOSTART mode is on.
- i** When introducing bipolar instruments into a metal trocar, the AUTOSTART function should be disabled beforehand!

Application areas

Bipolar coagulation with AUTOSTART

Suitable instruments

- Bipolar forceps

5.10.9 Bipolar scissors

This mode is used with bipolar scissors. It can be used for coagulation before or during mechanical cutting, point coagulation and surface coagulation.

Application areas

Preparation, coagulation and cutting of tissue

Suitable instruments

- Bipolar scissors

- i** Bipolar scissors should only be operated with the current type bipolar scissors for cutting or coagulation.

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- bipolare Pinzetten

5.10.10 Micro

This mode is used for arcless contact coagulation with micro forceps. It enables extremely fine control of power output down to 0.1 W and precise work for tightly restricted bipolar contact coagulation.

Application areas

Bipolar coagulation in pediatric surgery, neurosurgery, plastic surgery, etc.

Suitable instruments

- Bipolar forceps
- Micro forceps

5.10.11 Forced

This mode is used for forced coagulation with forceps.

i Given the higher voltage, we do not recommend the use of this current type for minimally invasive surgery. Note the permissible instrument voltage.

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

5.11 Menü

i Basic settings can be changed, e.g., operational language, audio, display and saving options.

5.11.1 Overview

The functions shown in the figure on the left (Fig. 28) are available.

Selecting a menu

▷ Press the corresponding selection button to open the menu.

Exiting a menu

▷ Press the 'back' button to return to the main screen.

5.11.2 Menu 'System settings'

In the 'System settings' menu the settings shown in the figure on the left (Fig. 29) can be changed.

Autostart delay

The time can be set under Autostart delay after which the HF generator is automatically activated impedance-dependent following tissue contact (see Fig. 30).

Incremental setting

▷ Press the buttons '+' and '-'.

Fast setting

▷ Move the slider to and fro.

Language selection

The following languages can be selected under 'Language selection':

German, English, Spanish, French, Portuguese, Italian, Russian, Polish, Turkish, Japanese, Korean, Thai, Indonesian, Chinese, Czech, Arabic, Hungarian, Danish, Finnish, Vietnamese, Swedish, Norwegian, Dutch, Bulgarian, Serbian, Lithuanian, Romanian, Slovakian

Key lock

In order to switch off the automatic screen lock or to set its duration press 'Key lock'. The duration can be set between 30 seconds and 5 minutes.

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OR system selection

The connection to the KARL STORZ OR system OR1™ can be created under the menu item 'OR system selection' (see Fig. 31). Upon delivery the KARL STORZ OR system OR1™ is in Standby mode. This is evident by the symbol

To create a connection to your OR1™ a LAN cable specified by KARL STORZ must be used. You can get more detailed information on this from our service team. If the OR system OR1™ is in Standby mode, a connection is established automatically when the connecting cable is inserted.

A detailed description can be found in the OR1™ manual SCB control NEO art. no. 96206572 D.

Deactivating the OR system (operating room without OR1 system):

1. Click on 'off'.
 2. Press OK to confirm.
- The OR1™ symbol is hidden. An automatic connection upon connecting the OR1™ system is no longer possible.

Returning the OR system (after deactivation) to Standby mode

1. Select 'OR1, Karl Storz'.
 2. Press OK to confirm.
- The OR1™ symbol is displayed

OR1 configuration (only for trained service personnel)

Special network configurations are possible by clicking on the symbol (see Fig. 31).

A detailed description can be found in the OR1™ manual SCB control NEO art. no.96206572 D.

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5.11.3 Menu 'Volume'

Use the 'Volume' menu to set the volume of the individual acoustic signals (see Fig. 32).

Incremental setting

▷ Press the buttons '+' and '-'.

Fast setting

▷ Move the slider to and fro.

i The volume of the activation signal should be increased as necessary for use in relatively noisy surroundings. The alarm tones have a minimum volume and limited changeability.

Mode	Category	Frequency (Hz)	Signal type
Unipolar Cut	Activation tones	635	Continuous sound
Unipolar Coag	Activation tones	475	Continuous sound
Bipolar Cut	Activation tones	565	Continuous sound
Bipolar Coag	Activation tones	505	Continuous sound
Sim Coag	Activation tones	755	Continuous sound
BiVascular-Safe end	Activation tones	-	Alternating sound
Footswitch changeover	Activation tones	-	Signal tone
ZAP mode	Activation tones	-	Signal tone
Fault	Alarm tones	-	Signal tone
Warning	Alarm tones	-	Signal tone
Note	Alarm tones	-	Signal tone

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5.11.4 Menu 'Service'

Service contact data can be called up in the 'Service' menu (see Fig. 33).

After entering a password, you can make further settings, e.g. resetting to the delivery state. To access the service level ('Service Tools') (see Fig. 34), enter the password 001224.

Saving device settings

With the function '1 Backup device' device settings can be saved on the KARL STORZ USB stick 20 0402 82. This includes all saved programs and system settings.

Installing device settings

Use '2 Restore device' to transfer saved device settings on the KARL STORZ USB stick to an AUTOCON® III 400.

Adding a start screen

Select the function '3 Add logo' to display a user-defined logo when starting up.

Deleting a start screen

Select the function '4 Remove logo' to remove a user-defined logo.

Program settings

Select the function '5 Program settings' to go to the setting level of the program/procedure configuration (see Fig. 35). In this menu restrictions ('Program restrictions') for the program handling and visibility of the procedures defined by the manufacturer can be set. Upon delivery, programs can be saved by the user (cannot be deleted or overwritten) and all KARL STORZ procedures are selected.

Saving logfiles

Select the function '6 Save logfiles' to save all logfiles onto the KARL STORZ USB stick.

Resetting to default setting

Select the function '7 Reset to default' to reset all settings and programs to the default setting.

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5.11.5 Menu 'System information'

In the 'System information' menu various system parameters such as version, serial number, TSI date for AUTOCON® III 400 and, if applicable, Argon-Beamer are displayed, as well as options (see Fig. 36).

For more details see the section Safety inspection, page 97.

5.11.6 Menu 'Programs'

1. In the 'Programs' menu you can select programs from a list and assign to or delete from favorites.
- or -
Click on the current program name in the main screen to access this menu quickly.

2. Press the relevant program name to select the program.
3. To navigate vertically in the program list, press the corresponding button ▲/▼. The programs are always sorted in alphabetical order.

4. Assignment to favorites is possible using the 'star buttons' at the bottom of the screen. Add programs to favorites using the 'green arrow' star button  (see Fig. 37 or 38 – left column) and remove these using the 'red arrow' star button  (see Fig. 39). The button  appears in place of the  button if a program is marked in the left-hand column.

5. Sort your favorites using the  /  star buttons.

6. Load the selected program by pressing 'OK'.
- or -
Return to the main screen by pressing 'back'.

 A maximum of 200 favorites can be selected.

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Saved programs can be deleted in the 'Programs' menu (if the function is enabled under 'Program settings' in section 5.11.4).

1. To this end, select the program to be deleted in the program list by clicking on the program name (see Fig. 36 or 37). To navigate vertically in the program list, press the corresponding button ▲/▼.
 2. To delete the selected program permanently, select the 'recycle bin' button.
- ☛ The selected program is deleted upon confirmation of a prompt.

5.11.7 Menu 'Favorites'

The defined favorites can be selected in the 'Favorites' menu. A fast selection of the favorites is possible using the star symbol in the main screen.

Access the next page of the favorites list using the buttons ◀/▶ in the lower part of the screen (see Fig. 40).

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

5.11.8 Menu 'Save program'

In the 'Save program' menu (see overview under 5.11.1) the current setting can be saved under the same or another program name (see Fig. 41). This function is enabled as standard (upon delivery); see 'Program settings' in section 5.11.4.

Program names can be created using the keyboard. Symbols, capital or small letters or numbers are selection options.

Program names spread over two lines can be assigned using the 'Enter' button.

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

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5.11.9 Socket settings

In the 'Socket configuration' menu the settings shown in the left-hand figure (Fig. 42) can be changed.

Expanding sockets

1. Call up 'Expand socket' in the 'Socket configuration' menu (see Fig. 43).
2. Click on  next to the two connectors to change the current setting (to divide the second bipolar socket on the bottom of the device into 2; see page 17, number .

☛ With active socket expansion the  on the left of the two connectors changes to .

3. Confirm with 'OK' to accept the selection.

4. Press the 'back' button to return to the main screen.

☛ You can now use the bipolar socket on the bottom of the device as two sockets.

ZAP mode

With the ZAP mode you can switch between two preset current forms for the same instrument.

1. Call up the ZAP mode in the 'Socket configuration' menu (see Fig. 44).
2. Activate or deactivate the changeover for the individual sockets by pressing the 'tick' next to the relevant socket (see Fig. 45).

If the ZAP mode is active, the 'tick' is displayed filled in.

Confirm with 'OK' to accept the selection.

Press the 'back' button to return to the main screen.

In the main screen a second level of the socket settings is shown for the active ZAP mode (see Fig. 46).

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You can switch in the main screen using either the handpiece or foot pedal.

For manual switching in the main screen tap on the ZAP mode symbol next to the socket settings (see Fig. 45).

On changing the ZAP mode level, the color of the symbol changes from *black* to *white* (see Fig. 46). Another way of changing the ZAP mode level is by actuating both buttons on the handpiece.

You can also change levels via the black changeover switch on the footswitch.

▷ To do this, select the ZAP mode under the 'Pedal' icon (see Fig. 47).

☛ The ZAP mode symbol now appears in the main screen below the 'Pedal' icon (see Fig. 48).

▷ To change the ZAP mode levels, actuate the black switch on the footswitch.

5.11.10 Menu 'System messages'

The 'System messages' menu can be accessed in the main screen via the 'Information' symbol.

Here the errors saved in the HF device which have occurred since switching on the HF device are called up (see Fig. 49).

These messages are not saved when switching off the HF device.

Calling up saved system messages:

1. Select a system message.
2. The selected system message is displayed again with '?'.
3. Press 'OK' to return to the overview.

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Argon flexibel:
COAG Flussrate: 0,4 l/min

5.11.11 Menu 'Argon'

In case of the successful recognition of an Argon beamer device, this menu is selectable in the status bar.

The 'Argon' menu enables the selection of argon flow rates for cutting and coagulation, as well as the selection of argon bottles and the display of filling levels.

The Argon menu can be selected quickly via the 'Argon' button in the main screen (see Fig. 50).

1. Use the '+' and '-' buttons to adjust the argon flow rate for cutting (Cut) or coagulation (Coag).
- or -
Use the slider to set the flow rate.
2. In case of two connected argon bottles use the buttons for output '1' or '2' to select the desired gas source.

Pressure reducers with electronic pressure sensor enable the display of exact filling levels of argon bottles.

3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.
- or -
Press the 'Back' button to return to the main screen without changing the selection.

i Default settings for argon flow rates according to the different modes are selected automatically:

Argon open:
CUT flow rate: 4.0 l/min
COAG flow rate: 3.0 l/min

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
i The procedure name is displayed in the status bar (see section 5.2.2).

5.11.12 Basic programs (procedures)

Procedures are programs in which parameters such as voltage and output are determined for specialist medical interventions. You can either create your own procedures/ programs (see section 5.11.8) or use procedures defined by KARL STORZ. These procedures are assigned to different disciplines, e.g. urology, laparoscopy and start with the number 0100. All specialist KARL STORZ procedures/programs and a standard program are enabled (visible) as standard (upon delivery). These programs and newly created ones cannot be changed or deleted (unless the corresponding setting in the menu Service – Program Restrictions, section 5.11.4, has been changed). In this menu you can display/ hide the specialist procedures to suit your needs. For selection see section 5.11.6 Menu ‘Programs’.

i The procedure name is displayed in the status bar (see section 5.2.2).

The predefined procedures/programs set in the factory for the AUTOCON® III 400 are divided into the following user groups:

- Specialization urology (0100-0105; see Fig. 51 to 56)
- Specialization gynecology (0200-0208; see Fig. 57 to 64)
- Specialization laparoscopy (0300-0303; see Fig. 65 to 68); Procedure 0303: BiVascularSafe only possible with models UH 401 and UH 401U
- Specialization gastroenterology (0400-0401; see Fig. 69 to 70)
- Application open surgery (0500-0501; see Fig. 71 to 72) Procedure 0501: BiVascularSafe only possible with models UH 401 and UH 401U
- Specialization pediatrics (0600; see Fig. 73)

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




- Specialization orthopedics and sports medicine (0700; see Fig. 74)
- Specialization ENT (0800; see Fig. 75)
- Specialization visceral surgery (0900; see Fig. 76); BiVascularSafe only possible with device models UH 401 and UH 401U

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


i An overview (screenshots) of the predefined procedures can be found on this page and the following page.

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5.11.13 Switching off the HF device (decommissioning)

1. Switch off the device at the Standby button on the front panel.
2. Switch off the device at the power switch ① on the rear (see Fig. 10).
3. Remove connected accessories from the device and, if necessary, reprocess.

i When the power switch on the rear is turned off, the Standby button (surrounding light) remains white for a few seconds. Once the surrounding light has gone off, the device is voltage-free.

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- System error
- EASY monitoring error

6.1 System information

System information is displayed as a message on the screen (see Fig. 77).

For further information on the cause and how to remedy this message, press the 'i' button in the status bar.

System information can be divided into three categories:

- Notices (gray display)
- Warnings (orange display)
- Faults (red display)

Notices are shown for five seconds on the screen. Warnings and faults are shown for ten seconds.

While a fault is present, activations are prohibited.

This information can be called up again under the now orange-framed 'i' in the status bar.

These messages are also available in the menu system information and will be deleted when switching off the device.

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The following table describes the cause of the error and the appropriate corrective action.

Header	Message text
Confirmation of AUTOSTART mode	You have selected an AUTOSTART mode. Setting the AUTOSTART mode can result in unintentional coagulations, e.g. when bipolar forceps are used for gripping while the AUTOSTART mode is on.
AUTOSTART Fault	The instrument is in contact with tissue. AUTOSTART cannot be selected when the instrument is in contact with tissue. Open up the instrument.
TSI Warning	The annual Technical Safety Inspection (TSI) is due.
Neutral Electrode Fault	No neutral electrode connected. No neutral electrode is connected. Connect a neutral electrode.
Neutral Electrode Fault	Wrong neutral electrode connected. The selection does not match the neutral electrode connected. Connect the neutral electrode which matches the selected mode, or change the mode to match the neutral electrode.
Neutral Electrode Fault	Poor contact with the patient. The resistance between the neutral electrode and the patient is too high. Establish better contact of the neutral electrode.
Neutral Electrode Warning	Poor contact with the patient. The contact resistance between the neutral electrode and the patient is getting worse. Establish better contact of the neutral electrode.
Neutral Electrode Fault	No cable for neutral electrode connected. No neutral electrode cable is connected. Connect a neutral electrode.
Neutral Electrode Notice	No cable for neutral electrode connected. The cable for neutral electrode has been removed. Unipolar activation is not possible.
Mode Fault	No mode selected. No mode was selected for this type of activation. Select the desired mode or change the footswitch assignment.
Mode Fault	This mode is not allowed for baby neutral electrodes. Use split neutral electrodes with a large surface and the suitable mode for this.
Mode Fault	This mode is not allowed for this socket. The current mode remains active. Choose another socket for this mode.
Footswitch Fault	No compatible footswitch. The connected footswitch is not compatible with this device. Connect a compatible footswitch with a changeover switch.
Footswitch Fault	Footswitch not assigned to a socket. The footswitch has not been assigned to a socket yet. Assign a socket to the footswitch using the 'Pedal' button.
Footswitch Fault	Error on footswitch connection. Check the footswitch. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.

Header	Message text
Fingerswitch Fault	Fault on fingerswitch connection. Check the handpiece and the connection cable. Please replace them if damaged. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.
Temperature Warning	The temperature of the device is higher than normal. The temperature of the device is elevated. This leads to a reduction of the maximum power.
Limitation of Continuous Activation	The maximum activation time has been exceeded. Please only activate the generator in short intervals, in order to avoid harming the patient and damaging the connected instruments or the generator.
Activation Fault	While switching on the device, there is an activation by footswitch, fingerswitch or AUTOSTART. Check the handpieces or footswitches for malfunctions. Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.
Activation Fault	There is an activation while connecting the footswitch or fingerswitch. Check the handpieces or footswitches for malfunctions. Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.
Activation Fault	There is no instrument connected on the activated socket. Connect an instrument on the designated socket.
Activation Warning	The mode for safety inspections is active. Activation is not possible. Quit this mode before activating again.
Bipolar Resection Warning	Use a coded KARL STORZ resection cable. Make sure that NaCl is used as an irrigation medium. Secure a continuous irrigation during the application. Always use conductive lubricants to avoid damages of the urethra. Avoid continuous activations.
GastroCut Warning	Polypectomy snare not in contact with tissue, or check connection cable at snare or generator. Please apply the snare and reactivate. First of all establish contact between tissue and polypectomy snare, or check the connection cable at the snare or the generator. Then activate with the yellow foot pedal.
BiVascularSafe Notice	Sealing has been completed.
BiVascularSafe Warning	There is a short-circuit in the area of the sealing instrument. Check the instrument for foreign bodies and whether there is contact with other objects. The sealing point must be free from foreign bodies.

Header	Message text
BiVascularSafe Warning	The sealing instrument is not in contact with tissue. Grasp the tissue and activate BiVascularSafe again. Check the connection to the generator. To seal the vascular and tissue bundles, it is necessary to grasp tissue and close the sealing instrument prior to activation.
Argon-Beamer Fault	Please connect the argon unit to the generator and switch it on. The argon device is connected to the generator with fiber optic cables. The active argon unit is automatically activated by the generator when switching to an argon mode.
Internal Error Argon-Beamer 5100	Argon-Beamer is inoperable. Connect an operational Argon unit to the generator. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Fault	The argon inlet pressure is too high. Max. inlet pressure: <4.5 bar Connect a source of argon gas in the appropriate pressure range. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
Argon-Beamer Fault	The argon inlet pressure has exceeded the permissible limits. Inlet pressure range: 2 – 4.5 bar Connect a source of argon gas in the appropriate pressure range. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
Argon-Beamer Warning	Mixed operation of argon bottles with and without an electric bottle pressure gauge is not recommended. Connect two identical pressure reducers.
Argon-Beamer Warning	Please check if the instrument is free of adhesions, and purge it with argon. If repeated purging does not solve the problem, the instrument and cable must be replaced.
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Warning	The filling level of the argon bottle is low. Please make sure that a replacement is available. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Fault	The argon bottle is empty. Connect a replacement bottle to enable activation. You can connect two argon bottles. The unit shifts automatically to the second bottle.
Argon-Beamer Notice	One argon bottle is empty. The unit has shifted automatically to the replacement bottle. Please make sure that a replacement is available.

Header	Message text
Argon-Beamer Fault	Please check if the argon bottles are connected and open. Empty bottles should be replaced. Subsequently restart the Argon-Beamer by activating the flashing 'Purge' button.
TSI Argon-Beamer Warning	The annual Technical Safety Inspection (TSI) for the Argon-Beamer is due.
Coded System Notice	The lifetime of the instrument is ending soon. Please order a replacement in good time. Any use of the instrument beyond its lifetime is not covered by warranty. Please contact your dealer in good time to purchase a new instrument.
Coded System Warning	The maximum lifetime of the instrument has been reached. Any further use is not covered by warranty. The maximum service lifetime of the instruments must not be exceeded, in order to guarantee safe usage. Any further use is at the user's risk.
Coded System Warning	A software update is necessary to use the Coded System with this instrument. Only carry out manual settings at this instrument. Please contact the Technical Support. Contact: MENU – SERVICE.
Coded System Warning	The preferred parameters for this coded instrument cannot be loaded. Only carry out manual settings at this instrument. Please contact the Technical Support. Contact: MENU – SERVICE.
Internal Error xxxx (e.g. with xxxx=4183)	If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.

With internal errors a number is shown in the header.

Please quote this number to the Technical Support.

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6.2 Error indications for EASY monitoring

Error indications are displayed in three stages (green, yellow and red) when problems occur.

When working with a split neutral electrode, the following errors may occur:

EASY monitoring	Cause	Indication	Remedy
Lights up yellow	Significant increase in resistance Depending on the indication, there may be heating under the neutral electrode	–	Stopping the application is not necessary. Check the positioning of the neutral electrode.
Switches from green to continuous red	A significant problem occurs when the unipolar current is activated.	An acoustic signal sounds. An error message appears on the display.	Check the neutral electrode and neutral electrode cable, see section EASY monitoring, page 28. ▷ Check the neutral electrode cable for proper connection and external damage.
	Loosened electrode	An acoustic signal sounds. A warning message appears on the display	▷ Correct the positioning of the neutral electrode. If the error persists, replace the neutral electrode.

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7 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: National laws and regulations must be observed.

The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

7.1 Reprocessing of accessories

▷ The reprocessing of accessories is described in the respective manual.

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7.2 Cleaning and disinfection

NOTE: Damage to the HF device due to incorrect reprocessing: No method other than the wipe-down disinfection described here is permitted.

WARNING: Fire hazard and risk from electric shock! Disconnect the plug before cleaning. Only use approved cleaning/disinfection agents for surface cleaning as per the manufacturer's instructions. It is essential to ensure that no liquids can penetrate the housing. Make sure that the AUTOSTART function is deactivated.

Wipe-down disinfection of devices and footswitches

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

KARL STORZ has qualified the resistance of the surfaces vis-à-vis wipe-down disinfection with the agent *Dr. Schumacher – Cleanisept Wipes*.

Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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8 Maintenance and repair

8.1 Maintenance

DANGER: Risk of infection! Carry out a surface disinfection and wrap the device in addition to the shipping packaging material before allowing the device to leave the hospital or practice to avoid spreading germs and infections.

- ▷ Check the device, the equipment cart and the accessories (e.g. footswitch, cable) after each use for damage or defects. In particular, make sure that the insulation is intact on all cables.
- ▷ Do not use any damaged device, damaged equipment cart or damaged accessories.
- ▷ Replace defective accessories immediately.
- ▷ Have the safety inspection for the device performed once a year. Please consult and comply with the respective service manual for additional technical information.

8.1.1 Safety inspection

Safety inspections must be performed once a year.

- ▷ The next safety inspection date of the AUTOCON® III 400 can be displayed in the menu, see section Menu 'System settings', page 71.

☛ A warning message appears during system start-up if a safety inspection is due. Press OK to confirm this message.

Any shorter safety inspection cycles specified in national regulations must be observed.

- ▷ The device and accessories may be inspected only by persons who have the required training, knowledge or experience and who can perform the inspection independently.
- ▷ With regard to the safety inspection, you must comply with the country-specific rules and regulations.

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The inspector documents the inspection results and measured values corresponding to the printed inspection protocol.

In the case of severe deviations from the values of the attached final acceptance report, or if the specified maximum values were exceeded:

- ▷ Send the HF device to the service center, see section 9.1 – Technical Support, page 100.

8.2 Repair

NOTE: You can damage the HF device by doing your own repairs and modifications of medical equipment! If a repair is necessary, only allow this to be done by the service center specified (see section 9.1). Never carry out any repairs yourself.

KARL STORZ is liable for safety, reliability and performance of the HF device under the following conditions:

- Full compliance with all instructions regarding the installation and proper use for the intended purpose contained in this instruction manual was maintained.
- Changes, repairs, new settings and similar procedures were carried out only by persons authorized to do this work by KARL STORZ.
- The electrical installations in the relevant room meet the local requirements and statutory provisions.

Fast and satisfactory repairs can only be guaranteed when all required data have been supplied in full.

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The following information is required for returning the device:

- Full address
 - Model number
 - Serial number
 - Software version
- ▷ Describe the problem, the corresponding application and the accessories used.
– or –
- ▷ Describe the repairs to be made.

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9 Storage

- ▷ If you store the HF device for longer than one year, pay specific attention to the indicators during automatic functional testing, see section Test for proper functioning, page. 26.
- ▷ Clean the HF device thoroughly before you put it into storage.
- ▷ Store the HF device in a clean, dry place in accordance with the storage conditions.

Storage conditions:

- Temperature: -20 °C to +50 °C
- Relative humidity: 0 to 75 %, non-condensing
- Air pressure: 500 to 1060 hPa

9.1 Technical Support

In Germany contact the following service center for maintenance and repair work:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

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10 Technical data
10.1 Technical data for AUTOCON® III 400

Insulation type / Classification	
EMC	IEC 60601-1-2
Protection class acc. to IEC 60601-1	I
Applied part type according to IEC 60601-1	CF
Standard compliance	IEC 60601-1: 2005+A1:2012, IEC 60601-1-2: 2007, IEC 60601-2-2: 2009, ISO 14971: 2007, ISO 13485: 2003 + Cor.1 2009
Classification according to Annex IX of the EC Directive 93/42/EEC	IIb
This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.	CE 0123

i The code number after the CE mark indicates the responsible notified body.

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Dimensions and weight	
Device dimensions (HxWxD)	177 x 447 x 457 mm
Net weight	12.5 kg
Packaging information/dimensions (HxWxD)	498 x 530 x 650 mm
Gross weight	18.2 kg

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Programs	
Number of programs in the device	300
Default programs, factory set	Yes
Individually programmable	Yes
Information shown on the display	Yes

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Neutral electrode monitoring	
EASY: Electrode Application System	Yes
Display of one-piece or split or Baby electrode	Main menu and neutral electrode selection
Contact resistance between individual sections of split neutral electrodes shown on display	Using color and contact indicator
Lead resistance shown on the display when a non-split neutral electrode is used	Yes
Maximum permissible resistance between the sections of a split electrode	300 Ω
Warning signal for hazardous conditions in connection with neutral electrodes	Visual, acoustic
Tones	Warning, activation and button tones, start melody
Warning message on the display	Text message with further information

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Safety features	
ISSys (Integrated Safety System)	Yes
ARC Control	ARC CONTROL
Continuous monitoring of HF leakage current and error message	Text message with further information
Dosage monitoring with error message on the display	Yes
Continuous self-test	Yes
Continuous status indication on the display	Yes
Operating errors shown on the display	Text message with further information
System errors shown on the display	Text message with further information
Safety inspection	Automatic reminder function (optional)
Instruction manual	Hard copy

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Documentation	
Data acquisition and storage in the device	System information
Documentation of error states	Yes
Documentation of operating errors	Yes
Retrieval of system messages via the display	Text message with further information

Communication	
Display	Capacitive touch screen 9"
External interface for communication between HF generator and an Argon beamer	Fiber-optic cables
USB interface for software updates	Yes
External PC interface for service support using KARL STORZ software	UART
External interface for communication between the HF generator and the KARL STORZ OR1™ system.	Ethernet connection

Service support	
Network connection for service support	Yes
Service support by service programs integrated in the device	Yes
Service support via ISSys	Yes

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Cooling	
Convection	Yes
Temperature-controlled fan	Yes

Operating mode	
Operating mode	Intermittent 10 s/30 s (on/off)

Specifications	
Max. monopolar power	400 W (at 200 Ω)
Max. bipolar power	400 W (at 75 Ω)
Output frequency	350 kHz / 1 MHz
Unipolar sockets	2x (footswitch and fingerswitch)
Bipolar sockets	3x (footswitch and fingerswitch)
Connection for footswitch	2x
AUTOSTART	Yes
Options	BiVascularSafe (UH 401 and UH 401U)
Scope of supply (see section 3.3)	Inspection protocol

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



Conditions of operation, transport and storage	Operation	Transport and storage
Temperature	+10 °C to +40 °C	-20 °C to +60 °C
Relative humidity	30 to 75 %, non-condensing	0 to 90 %, non-condensing
Air pressure	700 to 1 060 hPa	500 to 1 060 hPa
Operating altitude (max.)	3000 m above sea level	

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











i These max. values are not necessarily created at rated load.
The HF power limitation is subject to a tolerance of $\pm 20\%$.

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


10.2 Output, voltage and current diagrams
Unipolar cutting – laparoscopy

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Laparoscopy' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Laparoscopy'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Laparoscopy' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Laparoscopy' (idle mode)

Efecto	U (Vp)
1	400
2	450
3	560
4	650
5	650
6	700
7	700
8	700
9	750

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Standard' = 200 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Standard'. Rated load resistance = 200 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Standard' = 400 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Standard' (idle mode)

Efecto	U (Vp)
1	400
2	450
3	560
4	650
5	650
6	700
7	700
8	700
9	750

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

Measurement at ohmic resistances

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Micro' = 25 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Micro'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Micro' = 50 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Micro' (idle mode)

Efecto	U (Vp)
1	280
2	340
3	380
4	400
5	400
6	400
7	450
8	450
9	450

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
- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting Resection'. Rated load resistance = 500 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Resection' (idle mode)

Efecto	U (Vp)
1	650
2	700
3	700
4	700
5	750

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

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection C-Cut®' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Resection C-Cut®'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection C-Cut®' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Resection C-Cut®' (idle mode)

Efecto	U (Vp)
1	1400
2	1400
3	1400
4	1400
5	1500
6	1600
7	1600
8	1600
9	1600

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting SupraLoop'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting SupraLoop'. Rated load resistance = 100 Ω

Efecto	P (W)
1	300
2	350
3	400

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting SupraLoop' (idle mode)

Efecto	U (Vp)
1	650
2	650
3	650

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Dry' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Dry'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Dry' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Dry' (idle mode)

Efecto	U (Vp)
1	1400
2	1400
3	1400
4	1400
5	1500
6	1600
7	1600
8	1600
9	1600

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 1' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 1'. Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	400
2	21	400
3	26	400
4	30	400
5	35	400

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 1' cut phase

Effect	U (Vp)
1	750
2	750
3	750
4	750
5	750

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 1' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 2' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 2'. Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	400
2	21	400
3	26	400
4	30	400
5	35	400

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 2' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 2' (idle mode)

Efecto	U (Vp)
1	750
2	750
3	750
4	750
5	750

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 3' coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Gastro Loop 3'. Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	400
2	21	400
3	26	400
4	30	400
5	35	400

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Loop 3' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Loop 3' (idle mode)

Efecto	U (Vp)
1	750
2	750
3	750
4	750
5	750

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 1' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 1'.
Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	300
2	21	300
3	26	300
4	32	300
5	40	300

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 1' cut phase

Effect	U (Vp)
1	650
2	650
3	650
4	700
5	750

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 1' (idle mode)

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 2' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 2'. Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	300
2	21	300
3	26	300
4	32	300
5	40	300

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Measurement at ohmic resistances

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 2' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 2' (idle mode)

Efecto	U (Vp)
1	650
2	650
3	650
4	700
5	750

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 3' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Gastro Knife 3'.
Rated load resistance = 500 Ω

Efecto	P (W) Fase Coag	P (W) Fase Cut
1	17	300
2	21	300
3	26	300
4	32	300
5	40	300

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Gastro Knife 3' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Gastro Knife 3' (idle mode)

Efecto	U (Vp)
1	650
2	650
3	650
4	700
5	750

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Argon' = 150 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Argon'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Argon' = 300 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Argon' (idle mode)

Efecto	U (Vp)
1	400
2	450
3	560
4	650
5	650
6	700
7	700
8	700
9	750

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Unipolare Koagulation – Laparoskopie

Unipolar coagulation – Laparoscopy

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Laparoscopy' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Laparoscopy', Rated load resistance = 500 Ω

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Laparoscopy' = 120 W
- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Laparoscopy' (idle mode) = 1800 Vp

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

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Moderate' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Moderate'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Moderate' = 120 W

Effect	U (Vp)
1	250
2	250
3	250

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Moderate' (idle mode)

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced coag' = 40 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced coag'. Rated load resistance = 1000 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced coag' = 80 W
- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Forced coag' (idle mode) = 3500 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Resection'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 120 W
- HF output voltage U [Vp] for the setting 'Unipolar Coagulation Resection' (idle mode) = 2200 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Spray'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 120 W

Effect	U (Vp)
1	3000
2	3800
3	4600
4	5000

- Diagram of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Spray' (idle mode)

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Unipolare Koagulation – Forciert mixed

Unipolar coagulation – Forced mixed

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced mixed' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced mixed'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced mixed' = 120 W

Effect	U (Vp)
1	1500
2	2000
3	2500

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced mixed' (idle mode)

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
Unipolare Koagulation – Forciert cutting

Unipolar coagulation – Forced cutting

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
- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced cutting' = 125 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced cutting'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced cutting' = 250 W

Effect	U (Vp)
1	1500
2	1500
3	1300
4	1300

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced cutting' (idle mode)

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Unipolare Koagulation – Gastro Coag

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Gastro Coag' = 25 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Gastro Coag'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Gastro Coag' = 50 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Gastro Coag' (idle mode)

Efecto	U (Vp)
1	1800
2	2200
3	2800

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Unipolar coagulation – Argon flexible
(only in combination with an Argon beamer)

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flexible' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon flexible'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flexible' = 120 W

- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Argon flexible' (idle mode) = 4400 Vp

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Unipolar coagulation – Argon flex. pulse
(only in combination with an Argon beamer)

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flex. pulse' = 40 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon flex. pulse'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon flex. pulse' = 80 W

Effect	U (Vp)
1	4400
2	4400
3	4400

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Argon flex. pulse' (idle mode)

i The pulse frequency changes with the effect setting. The higher the effect level, the faster the pulse sequence.
Effect 1: 1 Hz, effect 2: 5 Hz, effect 3: 10 Hz
The mode 'Argon flexible' is paused due to these pulse sequences.

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Unipolar coagulation – Argon open
(only in combination with an Argon beamer)

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon open' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Argon open'. Rated load resistance = 500 Ω

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Argon open' = 120 W
- HF output voltage U [Vp] for the setting 'Unipolar Coagulation Argon open' (idle mode) = 4600 Vp

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Unipolare Koagulation – Cardiac Thorax

Unipolar coagulation – Cardiac Thorax

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Thorax' = 50 W

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- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Cardiac Thorax'. Rated load resistance = 500 Ω

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Thorax' = 100 W
- HF output voltage U [Vp] for the setting 'Unipolar Coagulation Cardiac Thorax' (idle mode) = 1800 Vp

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Unipolare Koagulation – Cardiac Mammaria

Unipolar coagulation – Cardiac Mammary

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Mammary' = 30 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Cardiac Mammary'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Cardiac Mammary' = 60 W
- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Cardiac Mammary' (idle mode) = 1800 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation SimCoag' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation SimCoag'.
Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation SimCoag' = 120 W

Effect	U (Vp)
1	2000
2	2500
3	4600

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation SimCoag' (idle mode)

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bip. Cutting'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 200 W
- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bip. Cutting' (idle mode) = 400 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bipolar Scissors'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bipolar Scissors' (idle mode) = 200 Vp

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Resection' cut phase

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
- Table of output power P [W] as a function of the setting 'Bipolar Cutting Bip. Resection'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Resection' phase after cutting

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Cutting Bip. Resection' (idle mode)

Efecto	U (Vp)
1	500
2	500
3	500
4	500

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Vaporization'

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- Table of power output P [W] as a function of the setting 'Bipolar Cutting Bip.Vaporization'.
Rated load resistance = 75 Ω

Efecto	P (W)
1	300
2	300
3	400

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Cutting Bip. Vaporization' (idle mode)

Efecto	U (Vp)
1	350
2	400
3	450

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation RoBi®'. Rated load resistance = 25 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 100 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation RoBi®' (idle mode) = 110 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Laparoscopy' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Laparoscopy'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Laparoscopy' = 120 W
- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Laparoscopy' (idle mode) = 150 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard' (idle mode) = 150 Vp

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
Bipolare Koagulation – Bip. Resektion

Bipolar coagulation – Bip. resection

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bip. Resection' = 175 W

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- Table of power output P [W] as a function of the setting 'Bipolar Coagulation Bip. Resection'. Rated load resistance = 75 Ω

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Coagulation Bip. Resection' (idle mode)

Efecto	U (Vp)
1	190
2	190
3	190
4	190

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
Bipolare Koagulation – Bip. Vaporisation

Bipolar coagulation – Bip. vaporization

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bip. Vaporization'

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- Table of power output P [W] as a function of the setting 'Bipolar Coagulation Bip. Vaporization'. Rated load resistance = 75 Ω

Efecto	P (W)
1	250
2	250
3	250

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
- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Coagulation Bip. Vaporization' (idle mode)

Efecto	U (Vp)
1	190
2	400
3	500

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
- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery'
- Power output P [W] with the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery' (Rated load resistance = 25 Ω) = 200 W
- HF output voltage U [Vp] with the setting 'BiVascularSafe'/'BiVascularSafe Open Surgery' (idle mode) = 190 Vp

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Bipolare Koagulation – Standard AUTO

Bipolar coagulation – Standard AUTO

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard AUTO'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 120 W
- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard AUTO' (idle mode) = 150 Vp

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Bipolare Koagulation – Bipolare Schere

Bipolar coagulation – Bipolar scissors

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Bipolar Scissors'. Rated load resistance = 75 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 120 W
- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Bipolar Scissors' (idle mode) = 200 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 10 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Micro'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 20 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Micro' (idle mode) = 90 Vp

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
Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Forced'. Rated load resistance = 50 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 100 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Forced' (idle mode) = 550 Vp

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10.3 Technical documentation

On request, the manufacturer will provide those service manuals, itemized parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

 We reserve the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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**11 Accessories/
Replacement parts**

Original KARL STORZ accessories are suitable for use with the AUTOCON® III 400 devices. When using accessories made by other manufacturers, the user must ensure that they are designed for and compatible with the maximum HF peak voltage of the HF device.

For the use and correct reprocessing of the autoclavable devices, compliance with the relevant instruction manuals accompanying these devices is required.

Detailed information on accessories and replacement parts is available in the corresponding KARL STORZ catalogs.

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12 EMC

12.1 Guidance and manufacturer's declaration as per IEC 60601-1-2, section 6.8.3.201

WARNING: Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this section during installation and commissioning.

WARNING: The AUTOCON® III 400 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile HF communication equipment may have an impact on this or other medical electrical devices.

WARNING: If KARL STORZ original accessories are not used, this may result in increased emissions or decreased immunity of the AUTOCON® III 400. Some KARL STORZ original accessories are also suitable for high frequency surgical units from other manufacturers. In use, the instruction manual for the high frequency surgical unit must be observed and the compatibility checked. The operation of KARL STORZ original accessories with high frequency surgical units from other manufacturers may result in increased interference emission or decreased immunity of the high frequency surgical unit.

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i The tables and guidelines that are included in this section provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If electromagnetic interference arises while using the device, the user can eliminate it by taking the following actions:

- ▷ reorient or relocate the receiving device
- ▷ increase the separation between the equipment
- ▷ connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Guidance and manufacturer's declaration – electromagnetic emissions (IEC 60601-1-2, Table 1)

AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.

Interference emissions measurements	Compliance	Electromagnetic environment – guidelines
HF emissions acc. to CISPR 11	Group 2	The AUTOCON® III 400 must emit electromagnetic energy in order to perform its intended function. Nearby electronic devices may be affected.
HF emissions acc. to CISPR 11	Class B	The AUTOCON® III 400 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity (IEC 60601-1-2, Table 2)

The AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	±6 kV contact discharge	±6 kV contact discharge	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30%.
	±8 kV air discharge	±8 kV air discharge	
Rapid transient electric interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV for input/output lines	
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
	± 2 kV line(s) to earth	± 2 kV line(s) to earth	
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	< 5% U_T (> 95% dip in U_T) for ½ cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	< 5% U_T (> 95% dip in U_T) for ½ cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AUTOCON® III 400 requests continued functioning even when interruptions of the power supply occur, it is recommended that the AUTOCON® III 400 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the a. c. mains voltage prior to application of the test level.

**Guidance and manufacturer's declaration – electromagnetic immunity
– for equipment and systems that are not life-supporting (IEC 60601-1-2, Table 4)**

The AUTOCON® III 400 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 400 must make sure that it is used in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
<p>Conducted HF disturbances acc. to IEC 61000-4-6</p> <p>Radiated HF disturbances acc. to IEC 61000-4-3</p>	<p>3 V_{effective value} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V_{effective value}</p> <p>3 V/m</p>	<p>Portable and mobile HF communications equipment should be used no closer to any part of the AUTOCON® III 400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1.2 \sqrt{P}$</p> <p>$d = 1.2 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz</p> <p>Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m].</p> <p>The field strength of stationary transmitters at all frequencies on site ^a should be, according to a study, less than the conformity level ^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed HF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AUTOCON® III 400 is used exceeds the applicable HF compliance level above, the AUTOCON® III 400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AUTOCON® III 400.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile HF communications equipment and the AUTOCON® III 400 (IEC 60601-1-2, Table 6)

The AUTOCON® III 400 is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the AUTOCON® III 400 can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications equipment (transmitters) and the AUTOCON® III 400 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance d according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
2	1.7	1.7	3.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

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13 Disposal

 The national regulations must be observed when disposing of or recycling the device or its components.

Symbol	Designation
	Devices bearing this symbol must be collected separately as electrical and electronic devices. Within the EU, disposal is handled free of charge by the manufacturer.

▷ For questions on the disposal of the device, please contact the Technical Support, see section Technical Support on page 100.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skårholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjial-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

CE 0123

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com