# STORZ

## KARL STORZ — ENDOSKOPE

**GEBRAUCHSANWEISUNG**  
TH 200 – VITOM® 3D-Exoskop  
TC 014 – IMAGE1 PILOT

**INSTRUCTION MANUAL**  
TH 200 – VITOM® 3D Exoscope  
TC 014 – IMAGE1 PILOT

**MANUAL DE INSTRUCCIONES**  
TH 200 – Exoscopio VITOM® 3D  
TC 014 – IMAGE1 PILOT

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**1 Important information  
for users of  
KARL STORZ devices**

**It is recommended to check the suitability of the products for the intended procedure prior to use.**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to aid you with proper connection and operation of the VITOM® 3D in conjunction with the instruction manual for the IMAGE1 S. All the necessary details and actions are clearly explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment. Insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the IMAGE1 S instruction manual for the individual functions of the CCU that can be used with the VITOM® 3D exoscope.

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

**2.1.1 Controls, displays, connectors and their uses**

- ① Rotary knob for horizontal image adjustment
- ② Sheath support for connection to holding system
- ③ Rocker switch – Scroll up/down in menu/programmed functions
- ④ Control button: Open menu, selection/programmed functions
- ⑤ Handle
- ⑥ Connection for fiberoptic light cable
- ⑦ Stereo lens
- ⑧ Light outlet – integrated illumination unit
- ⑨ Video connection cable

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2.2 IMAGE1 PILOT

**2.2.1 Controls, displays, connectors and their uses**

- ① Control buttons 1 – 4 for programmed functions
- ② LED for operating status display
- ③ 3D wheel for control of zoom, focus and ROI as well as navigation through the menu
- ④ Right control button: Cancel, jump back from a menu and Esc
- ⑤ Left control button: Open the camera menu, Setup menu and selection
- ⑥ KSLOCK quick-release coupling (male) for connection to a KARL STORZ holding system

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	Follow instructions for use
	Type CF device/ Defib. resistant
	CE mark

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
	Consult instructions for use
	Catalogue number
	Serial number
	Number of products in the product packaging
	Fragile, handle with care
	Keep dry
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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## **4 General information**

### **4.1 Description of the device**

The VITOM® 3D exoscope is a high-resolution stereoscopic visualization system with integrated zoom, focusing and image alignment mechanism for visualization of microsurgical interventions. The applications are comparable to those of a surgical microscope. With the VITOM® 3D exoscope, the displayed site is enlarged on a 3D screen.

The VITOM® 3D exoscope is compatible with the IMAGE1 S via the D3-LINK™.

The VITOM® 3D exoscope is positioned for use with the aid of a holding arm in such a way that the user can visualize the site.

The VITOM® 3D exoscope is controlled using the control unit supplied (IMAGE1 PILOT). The control unit is secured to the operating table using a second holding arm in such a way that the user can reach and operate it with ease.

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## 5 Safety instructions

### 5.1 Warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in the instruction manual, please read the text that follows carefully, in order to ensure safe and efficient operation of the device.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** Notes contain additional information on the safe operation of the device.

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**WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

**WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

**WARNING:** Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

**WARNING:** Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments which transmit energy pointed at the target in the field of vision during use.

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**WARNING:** Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Never allow the exoscope and the fiberoptic light cable to rest on the patient or come in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

**WARNING:** Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected fiberoptic light cable or an exoscope. Always wear standardized protective equipment when using lasers.

**WARNING:** Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

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**WARNING:** Risk of injury: Patient leakage currents may add if powered endotherapy devices are used simultaneously. This is particularly important if an endoscope with a Type CF applied part is used. In that case, an endotherapy device with a Type CF applied part must be used to minimize the total patient leakage current.

**WARNING:** The video system may stop working during use. An adequate replacement system must be kept available.

**WARNING:** Risk of electric shock: Do not open the device! Have servicing carried out by authorized personnel only. The removal of covers by unauthorized staff voids the warranty.

**WARNING:** Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the telescope is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the setting ring which can cause cardiac pacemakers or other electrical devices to malfunction.

**WARNING:** Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the product must not be used.

**CAUTION:** Protect the device from water splashes and other liquids.

**CAUTION:** Do not store liquids on or above the device.

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## 5.2 Applications

### 5.2.1 Intended use

The VITOM® 3D is used for the short-term, non-invasive visualization of the operating field during microsurgical, open and endoscopic interventions.

### 5.2.2 Indication

Extracorporeal visualization of the operating field during microsurgical, open and endoscopic interventions.

### 5.2.3 Contraindications

Use of the VITOM® 3D exoscope is contraindicated, if, in the opinion of the attending physician, the VITOM® 3D system is not suitable for visualization of the operating field, or the general condition of the patient does not permit the use of the VITOM® 3D system.

As KARL STORZ light sources are not intended for use in ophthalmology, the use of the VITOM® 3D system is also contraindicated in this field.

## 5.3 Patient profile

The use of the VITOM® 3D is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

## 5.4 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been trained thoroughly at least once in the use of the device

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- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means.

## 5.5 Intended conditions of use

### 5.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operating properly before using the device. During treatment with the VITOM® 3D exoscope the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

**i NOTE:** Any treatment may only be performed if visual observation of the action of the device is ensured.

### 5.5.2 Other intended conditions

Frequency of use:	one or more times a day
Length of use:	from a few minutes to several hours a day
Place of installation:	Positioning on a suitable holding system
Mobility:	can be positioned and moved in combination with a holding arm
Combination:	can be used on the patient at the same time as other equipment required for the operation
Control:	can be controlled via the IMAGE1 PILOT

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**5.6 Safety precautions when operating the instrument**

It is the user's responsibility to make sure the instrument is safe and operates properly before using it.

During treatment in conjunction with the VITOM® 3D exoscope, the patient must be treated and kept under observation using the usual medical care. This includes monitoring sterile application conditions where required by the type of intervention.

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#### 6.2.2 Optionales Zubehör

- 1 VITOM® 3D-Illuminator (TH 002)

## 6 Installation and operating instructions

### 6.1 Unpacking the equipment

Carefully remove the VITOM® 3D exoscope from its packaging. Check for missing items and evidence of damage during shipment.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

### 6.2 Components of the VITOM® 3D system

#### 6.2.1 Necessary basic equipment Camera system and monitor:

- 1 IMAGE1 S CONNECT module (TC 200 EN)
- 1 IMAGE1 S D3-LINK™ module (TC 302)
- 1 32" 3D monitor (TM 323)
- 1 3D polarization glasses (9800 GF)

#### VITOM® 3D exoscope:

- 1 VITOM® 3D exoscope (TH 200)
- 1 IMAGE1 PILOT (TC 014)
- 1 Sterile cover for VITOM® 3D exoscope (TH 001)
- 1 Sterile cover for IMAGE1 PILOT (041150-20)

#### Recommended light source:

- 1 Power LED 300 (TL 300)
- 1 Light cable with a diameter of 4.8 mm (e.g., 495 TIP or 495 VIT)

#### Holding arms:

- 1 Suitable KARL STORZ holding system for VITOM® 3D exoscope (e.g., 28272 HS with 28272 HP)
- 1 VITOM® 3D clamping jaw with ball joint for KARL STORZ holding systems (28272 VTK)
- 1 Suitable KARL STORZ holding system for the IMAGE1 PILOT (e.g., 28272 HA with 28172 HR)

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- 1 AIDA™ documentation system (WD 200 or WD 250)

When using the POINT SETTER holding system:

- 1 VITOM® 3D clamping jaw for POINT SETTER holding system (28272 VTP)

### 6.3 Compatible camera control units (CCU)

#### CCU

IMAGE1 S CONNECT (TC 200)  
from SW version **2.6** in combination  
with IMAGE1 S D3-LINK™ (TC 302)  
from SW version **2.6**

Please check the software version as follows:

1. Open the Setup menu (Setup key on the dashboard or Shift+Enter key on the keyboard).
2. Select "Module Information".
3. Check the software version.

If the above-mentioned or a higher software version (TC 200 – 2.6, TC 302 – 2.6) is not installed on your IMAGE1 S camera control unit, you can install the software version using the USB stick optionally provided or contact your KARL STORZ service partner.

If it is necessary to upgrade the hardware of the compatible Link module (TC 302), the system will show you this. For more information on this, please contact your KARL STORZ service partner.

### 6.4 Before each use

Before each use, the instrument must be inspected for damage, especially around the distal end (cover glass, light outlet, ...). Before each use, the video connecting cable and light cable must be inspected for damage, especially breaks and cuts.

The use of instruments with damaged cable must be discontinued. Cables must be replaced (light cable to be replaced by the user or the Service Dept., and the video connection cable by the Service Dept.).

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
### 6.5 General information on commissioning of the VITOM® 3D exoscope

The VITOM® 3D exoscope is always held over the site by a holding system during a practical application. The holding arm should be positioned in such a way that the operating theater procedures are not interrupted and the user has an optimal and ergonomic view of the 3D monitor. If possible, the monitor should therefore be positioned directly opposite the user at eye level.

**i NOTE:** The VITOM® 3D exoscope including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

**i NOTE:** Always observe the instructions for the holding system used and its suitability for use with the VITOM® 3D exoscope.

### 6.6 Attaching the VITOM® 3D exoscope to the holding arm

For operation of the VITOM® 3D exoscope, a KARL STORZ holding arm (e.g., endoscopic holding system) which can carry a load of at least 2 kg is required.

Depending on the holding system, the VITOM® 3D exoscope is attached to the KARL STORZ holding system using a clamping jaw (28272 VTK or 28272 VTP).

Attach the clamping jaw to the holding system correctly and ensure that the clamping jaw is firmly connected to the holding system (see Fig. A: Clamping jaw correctly attached to the holding system and Fig. B: Clamping jaw incorrectly attached to the holding system).

**i NOTE:** Observe the relevant instruction manuals for use of the clamping jaw 28272 VTK or 28272 VTP on the holding system (Art. no.: 97000154 or 97000155).

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Then open the clamping jaw 28272 VTK or 28272 VTP by simultaneously pressing the unlock buttons. Place the VITOM® in the clamping jaw and pay particular attention to correct orientation (see Fig. C).

**NOTE:** Due to the circumferential groove, the VITOM® 3D can only be placed in the clamping jaw in one direction.

Hold the VITOM® 3D in position and close the clamping jaw with your other hand until it clicks into place audibly (see Fig. D). Check whether the VITOM® 3D is safely held in the clamping jaw before you let it go. The clamping jaw ensures that the VITOM® 3D can rotate along its axis. This permits a comfortable position during the intervention (see Fig. E).

Depending on the holding system and clamping jaw, the VITOM® 3D can be swiveled in all directions either via the ball joint on the clamping jaw (Art. no.: 28272 VTK) or via the integrated ball joint in the POINT SETTER.

**NOTE:** Before every intervention, check whether the VITOM® 3D exoscope can be positioned using the holding system selected in such a way that optimal visualization of the site is ensured.

### 6.7 Connecting the light cable to the VITOM® 3D exoscope

The light cable is mounted by fitting to the light port provided and turning clockwise. Check that it is firmly connected (see Fig. F). Light cables with a fiber diameter of 4.8 – 5.0 mm are recommended for the VITOM® 3D.

**CAUTION:** Please note that the circular symbol on the VITOM® 3D and the circular symbol on the light cable must correspond (see following table).

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**Overview: Recommended combinations\***

Light cable diameter	
4.8 – 5.0 mm	⊙
3.0 – 3.5 mm	⊗
2.0 – 2.5 mm	⊚
Endoscope/exoscope diameter	
6.6 – 12.0 mm	⊙
3.0 – 6.5 mm	⊗
0.8 – 2.9 mm	⊚

\*Special endoscopes/exoscopes may deviate slightly

**6.8 Fitting the VITOM® 3D and holding arm with a sterile cover**

The VITOM® 3D exoscope cannot be sterilized and must therefore be fitted with a sterile cover for the duration of the intervention. The cover TH 001 has been specially designed for the VITOM® 3D exoscope and features a plastic frame with a transparent screen in the distal area.

Remove the sterile cover from its packaging (see Fig. G). Insert the plastic frame over the cover glass and light outlet on the distal end of the VITOM® 3D exoscope (see Fig. H). Then pull the cover over the VITOM® 3D exoscope and holding arm (see Fig. I).

**i NOTE:** The cover TH 001 has been specially designed for the VITOM® 3D. The use of other sterile covers is prohibited.

**i NOTE:** Observe the instruction manual with article number 97000169 for handling instructions for the sterile cover.

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### 6.9 Attaching the IMAGE1 PILOT

The VITOM® 3D exoscope system is always operated with the IMAGE1 PILOT (TC 014). It is attached to the operating table using a KARL STORZ articulated stand (e.g., base 28172 HR and articulated stand 28272 HA) and positioned in such a way that the user can reach the control unit with ease.

The holding arm is attached to the operating table and fixed in position as described in the corresponding instruction manual with Art. no. 96116335DF.

The IMAGE1 PILOT is equipped with a KSLOCK quick-release coupling (male) on its underside and can thus be easily attached to KARL STORZ articulated stands (see Fig. J).

To this end, open the KSLOCK quick-release coupling (female) on the articulated stand by turning in the 'OPEN' direction and insert the IMAGE1 PILOT so that it clicks into place audibly. Make sure that the IMAGE1 PILOT is correctly locked in (see Fig. L: fully locked in/ Fig. K: not fully locked in).

Then lock the KSLOCK quick-release coupling (female) on the articulated stand by turning in the 'CLOSE' direction (see Fig. M).

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Then connect the USB cable of the control unit to a USB socket on the IMAGE1 S CONNECT module (see Fig. N). The message 'IMAGE1 PILOT connected' appears on the screen. The LED for display of the function status must be lit up constantly.

**i NOTE:** The IMAGE1 PILOT calibrates the 3D Wheel position on power up. Therefore, do not press any keys or touch the 3D Wheel while attaching the USB connection, or powering up the IMAGE1 S system.

### 6.10 Fitting the IMAGE1 PILOT with a sterile cover

The IMAGE1 PILOT cannot be sterilized and is therefore fitted with a sterile cover during the intervention (041150-01).

Remove the sterile cover from its packaging and use it to cover the control unit and articulated stand (see Fig. O).

Make sure that the elastic insert on the tip of the sterile cover is placed directly over the IMAGE1 PILOT. The elastic insert increases the grip when operating the IMAGE1 PILOT (see Fig. P).

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**6.11 Connecting the VITOM® 3D to a light source**

Plug the light cable into the connecting socket on the light source as far as it will go. To disconnect the light cable pull on the connector (see Fig. Q), never on the cable itself.

**6.12 Connecting the VITOM® 3D to the D3-LINK™**

Insert the video connecting cable connector into the connection socket on the D3-LINK™.

**NOTE:** Do not use moist cables.

**NOTE:** The VITOM® 3D exoscope is first initialized after being connected. It takes approx. 20 seconds until an image is displayed.

**6.13 Switching on the CCU, light source and monitor**

**NOTE:** Heed the instruction manuals for the IMAGE1 S, the light source, the light cable and the monitor.

**6.14 Optimal light efficiency**

**NOTE:** The entrance and exit faces of the light cable and VITOM® 3D must be kept clean if optimal light transmission is to be maintained. The light cable must first be unscrewed (see Fig. R). Then clean the entrance and exit faces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70% isopropyl alcohol. Please note that both the glass surface inside the light port of the VITOM® 3D (see arrow, Fig. S) and on the light cable itself (see Fig. T) must be cleaned. Please observe the information on cleaning the light cable connection on page 32.

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### 6.15 Disassembly of the VITOM® 3D system after use

Switch off the equipment after using the VITOM® 3D exoscope. Remove and dispose of the sterile covers of the VITOM® 3D exoscope and IMAGE1 PILOT.

Remove the VITOM® 3D exoscope and IMAGE1 PILOT from the holding system, perform a wipe-down disinfection (also see chapter 8 'Reprocessing') and keep both safe.

**NOTE:** The handle on the VITOM® 3D exoscope can become very hot during operation, especially on the sides and on the light cable connection (over 50°C under certain circumstances). This is not a malfunction. It is recommended to allow the system to briefly cool down after use prior to disassembly of the holding system. To this end, please remove the camera head cable from the D3-LINK™ and switch the light source to standby mode.

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## 7 Operation of the VITOM® 3D exoscope system

### 7.1 Setting the horizon

There is a rotary knob on the distal end of the VITOM® 3D with which the stereo lens can be rotated through +/- 90° (see Fig. U). This enables flexible positioning of the 3D exoscope and correction of the horizontal image alignment during the application. To set the horizon, turn the black rotary knob in the desired direction (see Fig. V). If the image has to be rotated further, the VITOM® 3D image can be rotated digitally through 180° using the 'Image Flip' function. This makes it possible to rotate an image through 360°. If necessary, the 'Image Flip' function can be assigned to one of the 4 programmable buttons on the IMAGE1 PILOT or to the camera head buttons. Also see chapter 7.4 'Programmable function buttons on the IMAGE1 PILOT'.

**i NOTE:** To avoid damage to the VITOM® 3D exoscope, do not turn the knob beyond the limit stops.

### 7.2 Controlling the optical functions of the VITOM® 3D via the IMAGE1 PILOT

Most of the functions of the VITOM® 3D exoscope can be remote controlled easily via the IMAGE1 PILOT control unit. To do so, you do not have to hold the VITOM® 3D exoscope.

The IMAGE1 PILOT can be most easily operated using your finger tips (see Fig. W).

The most important functions are explained below.

**i NOTE:** For reasons of simplicity, a sterile cover is not displayed in this image. Always use a sterile cover for practical application in a sterile environment.

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### 7.2.1 Adjusting the magnification

The VITOM® 3D has an infinitely variable zoom which can be conveniently adjusted via the IMAGE1 PILOT.

To magnify the image, press the 3D Wheel of the IMAGE1 PILOT straight down (see Fig. X).

To reduce the image, pull on the 3D Wheel of the IMAGE1 PILOT.

The more you press the 3D Wheel, the faster you will be able to zoom in or out. This ensures that the desired magnification can be set precisely.

### 7.2.2 Adjusting the focus

The working distance of the VITOM® 3D can be variably adjusted from 20–50 cm. The 3D Wheel enables focusing of the image at any time.

To bring a lower level (greater working distance) into focus, turn the 3D Wheel clockwise.

To bring a higher level (shorter working distance) into focus, turn the 3D Wheel counterclockwise (see Fig. Y).

The more you turn the 3D Wheel, the faster you will be able to focus. This ensures that the required focus point can be set with great precision.

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### 7.2.3 Selectable region of interest during magnification (ROI)

Whenever a magnification greater than 1x is selected, the selected region of interest can be moved upward, downward, to the left or to the right. This facilitates centering of the magnified image without having to move the holding arm.

To move the region of interest, push the 3D Wheel in the required direction (see Fig. Z and a).

The more you press the 3D Wheel, the faster you will be able to move the region of interest.

**i NOTE:** The region of interest can only be moved in the image level. This does not change the angle of view of the VITOM® 3D relative to the object. The region of interest can only be moved with a magnification setting of >1x. The movement range is limited. The image cannot be moved beyond the margins of the fully zoomed out image.

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
By pressing the 3D Wheel on the IMAGE1 PILOT, an overview screen (minimap) is displayed in the top right-hand corner which provides the user with visual feedback on the zoom area and the position of the selectable region of interest (ROI) (see Fig. b).

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The large rectangle shows the available overview area. As soon as you use the zoom, a small filled rectangle which represents the selectable region of interest appears. The size of the rectangle varies with the selected zoom level. If the region of interest is moved, the rectangle also moves and shows the position within the available overview area. Under the overview area, a bar shows the relative zoom level. The bar is filled as the magnification factor increases.

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### 7.3 Control of the IMAGE1 S using the IMAGE1 PILOT

The IMAGE1 PILOT permits direct access to the menu of the IMAGE1 S so that the surgeon can operate all the functions from the sterile field as necessary (see Fig. c).

As soon as the menu of the IMAGE1 S is called up using the 'ENTER' button, you can navigate through the menu, adjust settings or call up functions by moving the 3D Wheel. You can exit the menu by pressing the ESC button.

The following table offers an overview of the functions.

**i NOTE:** As long as the IMAGE1 S menu is active, the VITOM® 3D functions (e.g., zoom and focus) cannot be controlled.

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



#### 7.4 Programmable function buttons on the IMAGE1 PILOT

The 4 buttons on the IMAGE1 PILOT can be programmed in a similar way as the camera head buttons and can be assigned functions (see Fig. d and e). In doing so, certain functions can be performed quickly or rapid access to menus can be provided.

The programmable buttons can be assigned the following functions:

- White Balance
- Light source intensity
- Capture
- Video
- S technologies
- Brightness
- Nothing
- Source swap
- AEL Activate
- Image Flip

#### 7.5 Advanced Exposure Levelling (AEL) function

The function **Advanced Exposure Levelling (AEL)** stabilizes the Exposure control in situations with extreme illumination conditions, and prevents the VITOM® 3D from under- or overexposing the image, for example when a shiny instrument is inserted.

The effect is shown schematically in the following images:

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Well exposed image through automatic exposure control.

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By inserting the shiny instrument, the automatic exposure control dims and the operating field is displayed too dark.

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The exposure is stabilized when the AEL function is activated. The reflection on the instrument may become larger, however the operating field is displayed well exposed.

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The AEL function can be assigned to one of the programmable buttons on the IMAGE1 PILOT. As soon as the button is pressed, the AEL function is active. The AEL logo appears on the screen. The exposure does not change even when inserting instruments.

The AEL function remains active until one of the following events takes place:

- AEL button is pressed again
- The region of interest (ROI) is considerably changed using the IMAGE1 PILOT
- The zoom is considerably changed using the IMAGE1 PILOT

After deactivating the AEL function, the VITOM® 3D returns to the normal exposure mode.

### 7.6 Extended PIP-Swap function

The PIP-Swap function can be assigned to one of the programmable buttons on the IMAGE1 PILOT. As soon as more than one Camera LINK module is connected to the IMAGE1 Platform, the PIP-Swap function permits switching between the cameras defined as Primary and Secondary Source. For example, this may be beneficial during endoscopically assisted microsurgery (EAM) to quickly switch from the VITOM® 3D to an endoscopic camera.

Pressing the function button to which PIP-Swap is assigned switches the image between the primary source and secondary source.

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**7.7 Operation of the VITOM® 3D  
exoscope without IMAGE1  
PILOT**

The VITOM® 3D exoscope can also be operated without the IMAGE1 PILOT. This is necessary if no IMAGE1 PILOT is connected or if the IMAGE1 PILOT fails.

If one of these events occurs, the VITOM® 3D exoscope can alternatively be operated via the touchpad.

**i NOTE:** For safety reasons, always connect the keyboard to the IMAGE1 S.

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	<b>Zoom</b>
	By moving your finger upward on the touchpad, the image view is magnified.
	By performing a downward movement on the touchpad, the image view is reduced.
	<b>Focus</b>
	Move to the left on the touchpad to focus in the distance.

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**8 Reprocessing not validated for US customers**

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**NOTE:** Instructions for US customers can be found in Chapter 9.

**8.1 General care instructions**

**NOTE:** Always disconnect the power cord from the Camera Control Unit before cleaning.

We recommend storing the VITOM® 3D in its original packaging when not in use. Do not coil the video connecting cable of the VITOM® 3D into tight loops or cause it to kink while storing. Do not expose the VITOM® 3D unprotected to direct sunlight or excessive heat.

When returning the VITOM® 3D for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The entrance and exit faces of the light cable and exoscope must always be kept clean (see corresponding instruction manual) to ensure optimal reproduction of the image.

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**8.1.1 Care and handling**

To maximize the life of the VITOM® 3D, please follow the guidelines given below for proper handling and care.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. The connector must always be completely dry and free of residue. Where necessary, remove any moisture with a low lint disposable cloth before plugging in the device. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. If necessary, remove remaining moisture with a low lint disposable cloth. The connectors of the video connecting cable and of the light cable must also be completely dry before connecting!

NEVER PICK UP THE VITOM® 3D BY ITS CABLE! Always pick up the set by the instrument handle first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE VITOM® 3D FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to the device. When transporting, use the original packaging to prevent damage to the device.

DO NOT STORE THE VITOM® 3D WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Store the cable in loops at least 15 cm in size.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful 'kinks' can result. Always allow the cable to uncoil gently and loosely.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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**8.2 General warnings**

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

**NOTE:** Risk of infection: Always use the medical device with the sterile cover (TH 001), which is available to order separately.

**WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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**CAUTION:** Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at [www.karlstorz.com](http://www.karlstorz.com).

**CAUTION:** National laws and regulations must be observed.

**NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

### 8.3 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility and necessary exposure time for the disinfectant must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth. Clean and wipe down the cover slips of the VITOM® 3D exoscope for the telescope and lighting, making sure that there are no stripes. Ensure that the lenses are not scratched.

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#### 8.4 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

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**9 Reprocessing for US customers only**

**WARNING:** The light cable must be removed for reprocessing of the VITOM® 3D. To do so, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the VITOM® 3D and must not be removed for reprocessing.

**9.1 General care instructions**

**NOTE:** Always disconnect the power cord from the Camera Control Unit before cleaning

We recommend storing the VITOM® 3D in its original packaging when not in use. Do not coil the video connecting cable of the VITOM® 3D exoscope into tight loops or cause it to kink.

Do not expose the VITOM® 3D unprotected to direct sunlight or excessive heat.

Prior to returning for service, the VITOM® 3D must be thoroughly cleaned per validated procedures. The entrance and exit faces of the light cable and exoscope must always be kept clean (see corresponding instruction manual) to ensure optimal reproduction of the image.

**CAUTION:** Chemical residue on the card-edge connector can cause failure of the camera electronics.

**CAUTION:** Heat from the light source can cause chemical residue on the light cable and on the light entrance of the exoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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### 9.3 Manual wipe-down disinfection

After each use, the device must be wiped down with a low level disinfection wipe with 0.1 – 1% Sodium Hypochlorite Wipe. In addition, clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol.

### 9.4 Inspection after cleaning

The cleaned medical device must be visually inspected to verify that it is clean, complete, free of damage, and dry.

- If residues or contamination are still present, the medical device must be wipe-down disinfected again.
- The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.

### 9.5 Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

**CAUTION:** Inspect the video connecting cable and the fiberoptic light cable for breaks or cracks before and after disinfection, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cables and cause damage. Where a VITOM® 3D has a damaged fiberoptic light cable or damaged video connecting cable, the relevant cable must be replaced.

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**10 Servicing****10.1 Maintenance and safety check****10.1.1 Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

**10.1.2 Safety check**

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such metrological checks and who are not subject to any instructions as far as such checking activities are concerned.

**Visual inspection:**

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the labels relevant to safety are legible.

**Electric measurements:**

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.

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**Test for proper functioning:**

- A functional check must be carried out in accordance with the instruction manual.
- The metrological check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

**10.2 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

**CAUTION:** Do not open the instrument or the lens. Opening these sealed assemblies affects the camera's tightness and voids all warranties. The covers on the camera processor must only be removed by authorized personnel.

**10.3 Disposal**

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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#### 10.4 Repair program

Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In the case of repairs, please contact your local KARL STORZ branch or authorized dealer.

In Germany, you can refer repairs direct to:

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Take-off-Gewerbepark 83  
78579 Neuhausen

#### Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### 10.5 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

### 10.6 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Dimensions (LxWxH):	approx. 351 x 45 x 80 mm
Weight:	approx. 950 g
Working distance:	20 – 50 cm
<b>Buttons/Input options</b>	
• 3 camera head buttons	
• Rotary knob for mechanical image alignment:	+/- 90°
<b>Interfaces</b>	
• Video interface for the IMAGE1 S D3-LINK™ module	
• Interface for fiberoptic light cable	
• Interface for attaching a holding system using the clamping jaw 28272 VTK or 28272 VTP	
Operating temperature:	5°C ... 30°C
<b>Storage</b>	
Ambient temperature:	-10°C ... +60°C
Humidity:	5 ... 95% (non-condensing)
Air pressure:	500 ... 1060 hPa

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**11.2 IMAGE1 PILOT technical data**

Dimensions (LxWxH):	approx. 139 x 93 x 77 mm
Weight:	approx. 354.6 g
<b>Buttons/Input options</b>	
<ul style="list-style-type: none"> <li>• 1 3D Wheel with 6 DoF (Degrees of Freedom)</li> <li>• 4 programmable buttons</li> <li>• 2 buttons for menu navigation</li> </ul>	
<b>Interfaces</b>	
<ul style="list-style-type: none"> <li>• USB interface for operation with the IMAGE1 S CONNECT</li> <li>• KSLOCK for fast and secure fixation on KARL STORZ holding systems</li> </ul>	
Operating temperature:	5°C ... 30°C
<b>Storage</b>	
Ambient temperature:	-10°C ... +60°C
Humidity:	5 ... 85% (non-condensing)
Air pressure:	500 ... 1060 hPa

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**11.2.1 Standard compliance**

for VITOM® 3D exoscope (TH200) in connection with the KARL STORZ IMAGE1 S CCU

**According to IEC 60601-1, UL 60601-1, CAN/GSA C22.2 No. 601.1-M90:**

- Type of protection against electric shock: Protection Class III
- Degree of protection against electric shock: Defibrillation-protected Type CF 

**11.2.2 Directive compliance (for VITOM® 3D)**

In accordance with the Medical Device Directive (MDD):

This medical device belongs to Class I  
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

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**12 Information on  
electromagnetic  
compatibility (EMC)**

**WARNING:** Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The VITOM® 3D exoscope model TH 200 corresponds to IEC 60601-1-2 [CISPR 11 Class A] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

The employed limit values provide a basic measure of safety from typical electromagnetic interference that is to be expected in a medical environment. The VITOM® 3D exoscope model TH 200 is a Group 1 device (according to CISPR 11). Group 1 contains all of the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If electromagnetic interference arises while using the device, the user can eliminate it by taking the following actions:

- reorient or relocate the receiving device;
- increase the separation between the equipment;
- connect the devices to different electrical circuits.

In case of additional questions, please contact your area representative or our Service Department.

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**Information on  
electromagnetic  
compatibility (EMC)**

**WARNING:** The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

**WARNING:** The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the device. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The VITOM® 3D exoscope model TH 200 is intended for use in the environment specified below. The customer or user of the VITOM® 3D exoscope model TH 200 should make sure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	The VITOM® 3D exoscope model TH 200 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class A	The VITOM® 3D exoscope model TH 200 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions acc. to IEC 61000-3-3	Complies	

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<b>Table 2</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
<p>The VITOM® 3D exoscope model TH 200 is intended for use in the electromagnetic environment specified below. The customer or user of the VITOM® 3D exoscope model TH 200 should make sure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5% $U_T^*$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T^*$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

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<b>Table 4</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting</b>			
The VITOM® 3D exoscope model TH200 is intended for use in the electromagnetic environment specified below. The customer or user of the device should make sure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbance acc. to IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz	3 V <sub>eff</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the VITOM® 3D exoscope model TH200, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5 / 3] \sqrt{P}$
Radiated RF disturbance acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5 / 3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7 / 3] \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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<b>Table 6</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the VITOM® 3D exoscope model TH200</b>			
<p>The VITOM® 3D exoscope model TH200 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter  [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz  $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz  $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz  $d = \left[ \frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance <b>d</b> in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where <b>P</b> is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.</p> <p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga  
Ontario, L5N 3R3, Canada  
Phone: +1 905 816-4500, Fax: +1 905 858-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8526  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
175 Cremona Drive  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia  
Miramar Trade Center,  
Edificio Jerusalem, Oficina 108,  
La Habana, Cuba  
Phone: +53 72041097, Fax: +53 72041098  
E-Mail: karlstorz@enet.cu

KARL STORZ Endoscopia México S.A de C.V.  
Lago Constanza No 326,  
Col. Granada Del. Miguel Hidalgo  
C.P. 11520 México D.F.  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Østensjøveien 15B  
0661 Oslo, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
12739 Skärholmen, Sweden  
14108 Kungens Kurva, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Valimotie 1 B, 4. kerros  
00380 Helsinki, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG  
Representation Office  
Kestucio g. 59  
08124 Vilnius, Lithuania  
Phone: +370 68567000  
E-Mail: ausra.kukeniene@karlstorz.com

KARL STORZ Endoskopi Denmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +32 473 810 451  
E-Mail: info@stopleer.be

KARL STORZ Endoscopia France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Verovškova c. 60A  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gilwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd., Hong Kong  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch,  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
Beijing, 100022, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shanghai Branch  
36F, Tower 1 Grand Gateway  
No. 1 Hongqiao Road, Xuhui District,  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
Chengdu, 610044, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
Shenyang, 110014, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
Guangzhou, 510620, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.  
I-Park Bundang1, 102dong, 8F  
239 Jeongjail-Ro, Bundang-Gu  
Seongnam, 463-859, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: mingwang@karlstorz.tw#

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

 **KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen  
Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)