

TruPort Support Head

Think laterally



The TruPort Support Head – now horizontal

The reliable TruPort modularity combined with space saving horizontal orientation makes TruPort Horizontal the ideal solution for the anesthesia workstation. Positioned directly above the anesthesia machine, TruPort Horizontal provides more freedom of movement in the OR and enables efficient workflow around the anesthesia workstation. Through modular architecture, supply modules for gases, power and data can be simply repositioned or even extended at any time according to your requirements – even post installation.


More flexibility for anesthesia

- 2 lengths available (500mm / 750mm)
- To be combined with all TruPort arm systems
- Up to 4 side rails laterally mounted
- With the M Rail and MPC Rail available
- Maximum pay load of the head is 80kg
- Suitable TruPort accessories available



Compact

More space in the OR through the compact architecture of this horizontal TruPort support head. Positioned directly below the ceiling allows maximum headroom and freedom of movement.



Adaptable

React quickly and efficiently to new requirements in the OR. The simple repositioning and possibility to retrofit supply modules for gases, power and data, TruPort Horizontal provides the flexible solution for the future.



Ergonomic

User-friendly control handles and an optional remote control allow ergonomic control for the user. Compatible with all arm systems, TruPort Horizontal can be easily controlled in all situation.