

## ...for Examination, Treatment - height adjustable short tables

Please choose - without surcharge – between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandora Plus (20 colours)



Picture with extras:  
Castors 100 mm diam.,  
each with total lock  
(4x Mod. 1120)

All tables are produced in our company in Bielefeld, Germany. On request the tables may be delivered in special measures (Length, Width, Height), additional functions or other frame colours. **Please contact us.**

### Model 2100...XL and Model 2105...XL and Model 2110...XL

**Alternative 1:** Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table.

**Alternative 2:** Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch)  
Time of electric level positioning from min. to max.: 22 sec.

#### Basic equipment:

- Lying surface 1-part: Length: 1.400 mm (other length on request)
- Not movable. Possibilities for mobility see list of extras, page 34.
- Standard-Thickness of upholstery: 6 cm.
- Upholstery base: Laminated wood.
- Load capacity: 250 kg (max. weight of patient)

Model	Width	Vertical height adjustment by	Frame colour:
2100-00XL/H	650 mm	Hydraulic pump: 500 - 920 mm	White powder painted
2100-04XL/H	650 mm	Hydraulic pump: 500 - 920 mm	Silver powder painted
2105-00XL/H	800 mm	Hydraulic pump: 500 - 920 mm	White powder painted
2105-04XL/H	800 mm	Hydraulic pump: 500 - 920 mm	Silver powder painted
2110-00XL/H	1.000 mm	Hydraulic pump: 500 - 920 mm	White powder painted
2110-04XL/H	1.000 mm	Hydraulic pump: 500 - 920 mm	Silver powder painted
2100-00XL	650 mm	Electric motor: 480 - 920 mm	White powder painted
2100-04XL	650 mm	Electric motor: 480 - 920 mm	Silver powder painted
2105-00XL	800 mm	Electric motor: 480 - 920 mm	White powder painted
2105-04XL	800 mm	Electric motor: 480 - 920 mm	Silver powder painted
2110-00XL	1.000 mm	Electric motor: 480 - 920 mm	White powder painted
2110-04XL	1.000 mm	Electric motor: 480 - 920 mm	Silver powder painted



Picture with extras:  
Castors 100 mm diam.,  
each with total lock  
(4x Mod. 1120)

All tables are produced in our company in Bielefeld, Germany. On request the tables may be delivered in special measures (Length, Width, Height), additional functions or other frame colours. **Please contact us.**

### Model 2114...XL and Model 2115...XL

**Alternative 1:** Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table.

**Alternative 2:** Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch)  
Time of electric level positioning from min. to max.: 22 sec.

#### Basic equipment:

- Lying surface 1-part: Length: 1.400 mm (other length on request)
- With 3-side upholstered border in covering of lying surface colour  
Border height: 150 mm
- Not movable. Possibilities for mobility see list of extras, page 34.
- Standard-Thickness of upholstery: 6 cm.
- Upholstery base: Laminated wood.
- Load capacity: 250 kg (max. weight of patient)

Model	Width	Vertical height adjustment by	Frame colour
2114-00XL/H	650 mm	Hydraulic pump: 500 - 920 mm	White powder painted
2114-04XL/H	650 mm	Hydraulic pump: 500 - 920 mm	Silver powder painted
2115-00XL/H	800 mm	Hydraulic pump: 500 - 920 mm	White powder painted
2115-04XL/H	800 mm	Hydraulic pump: 500 - 920 mm	Silver powder painted
2114-00XL	650 mm	Electric motor: 480 - 920 mm	White powder painted
2114-04XL	650 mm	Electric motor: 480 - 920 mm	Silver powder painted
2115-00XL	800 mm	Electric motor: 480 - 920 mm	White powder painted
2115-04XL	800 mm	Electric motor: 480 - 920 mm	Silver powder painted

**Tables with electric motor powered adjustment** must have a safety-lock-box to protect against unauthorized usage (regulations of BfArM based on EU-ordinance acc. EN 60601). Please choose between our 2 possibilities:

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)