## ...for Examination, Treatment, Massage

## - height adjustable manual by hydraulic lift





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)

Our Models 2000XL/H / 2001XL/H and 4000/H / 4001/H are high grade versatile tables.

Please consider for decision to buy a table with height adjustability by electric motor or by manual hydraulic lift:

- Tables with height adjustment by hydraulic pump we recommend by predominant mobility usage.
- The height adjustment is possible without electric current
- The height is manual adjustable by foot operated pump (repeated treading of the foot pedal).
- For operator convenience the foot-pedals always fitted (as standard) bothsides the table.



Picture with extras: - Castors 100 mm diam. with central locking system (Mod. 1180N), - Paper roll holder (Mod. 1270K)

Model 2000-00XL/H: Frame white powder coated Model 2000-04XL/H: Frame silver powder coated

Manual foot operated height adjustment by hydraulic pump from 500 mm to 920 mm.

Pedals for foot-controlled height positioning are bothsides the table. **Basic equipment:** 

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring
- Length head part: 550 mm
- Foot (back) part fixed. Length: 1.400 mm
- Not movable. Possibilities for mobility see list of extras, page 34.

Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 1.950 mm, Width: 700 mm Load capacity: 250 kg (max. weight of patient)

Model 2001-00XL/H: Frame white powder coated Model 2001-04XL/H: Frame silver powder coated

like Mod. 2000...XL/H, but difference: Width of lying surface: 800 mm



Picture with extras: Comfort-Castors 125 mm diam. with central locking system (Mod. 1185)

Model 4000-00/H: Frame white powder coated Model 4000-04/H: Frame silver powder coated

Manual foot operated height adjustment by hydraulic pump from 440 mm to 900 mm.

Pedals for foot-controlled height positioning are bothsides the table. **Basic equipment:** 

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring
- Length head part: 500 mm
- Foot (back) part fixed. Length: 1.500 mm
- Not movable. Possibilities for mobility see list of extras, page 34.

Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 2.000 mm, Width: 700 mm Load capacity: 250 kg (max. weight of patient)

Model 4001-00/H: Frame white powder coated Model 4001-04/H: Frame silver powder coated

like Mod. 4000.../H, but difference: Width of lying surface: 800 mm

### Protection for upholstery of Models 2000XL, 2001XL and 4000, 4001



Picture with extras: - Castors 100 mm diam. with central locking system (Mod. 1180N)

- Wall protection bumpers (Mod. 1175-00)
- Upholstery shock-protection (Mod. 1178-00)

Mod. 1175-00- Model-No.: Wall protection bumpers, holder white Mod. 1175-04- Model-No.: Wall protection bumpers, holder silver For Models 2000XL and 4000 (lying surface 700 mm)

4 wall protection bumpers to retrofit at the edges of head and foot part

Mod. 1176-00- Model-No.: Wall protection bumpers, holder white Mod. 1176-04- Model-No.: Wall protection bumpers, holder silver For Models 2001XL and 4001 (lying surface 800 mm)

4 wall protection bumpers to retrofit at the edges of head and foot part

Mod. 1178-00- Model-No.: Upholstery shock-protection, holder white Mod. 1178-04- Model-No.: Upholstery shock-protection, holder silver For Models 2000XL and 4000 (lying surface 700 mm)

2 longitudinal protection rods for foot upholstery to assemble at left and right side of the table

Mod. 1179-00- Model-No.: Upholstery shock-protection, holder white Mod. 1179-04- Model-No.: Upholstery shock-protection, holder silver For Models 2001XL and 4001 (lying surface 800 mm)

2 longitudinal protection rods for foot upholstery to assemble at left and right side of the table

# ...for Examination, Treatment, Massage - height adjustable by electric motor drive





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)

Our Models 2000XL / 2001XL and 4000 / 4001 are high grade versatile tables.

Please consider for decision to buy a table with height adjustability by electric motor or by manual hydraulic lift:

- Tables with electric motor must have a 230 V power connection for motor function.

  With lifting force of 7.000 to 10.000 N the motor is very powerful with constant speed burdened or
- unburdened. A safety-lock-box to protect against unauthorized usage is necessary.

   To handle the motor you can choose between hand or foot swich, as standard the table has a foot switch.
- Tables with electric motor lift must be checked-up under specifications of IEC 60601 / IEC 62353.



Picture with extras: - Castors 100 mm diam., each with total lock (4x Mod. 1120)

- Paper roll holder (Mod. 1270F)

## Model 2000-00XL: Frame white powder coated Model 2000-04XL: Frame silver powder coated

Height adjustment by electromotive drive from 480 mm to 920 mm.

Time of electric level positioning from min. to max.: 22 sec.

The electric motor is controlled by 1 foot switch

(possible option: Hand switch).

### Basic equipment:

Lying surface 2-sectioned:


- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring
- Length head part: 550 mm
- Foot (back) part fixed. Length: 1.400 mm
- Not movable. Possibilities for mobility see list of extras, page 34.

Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

**Measures:** Length: 1.950 mm, Width: 700 mm Load capacity: 280 kg (max. weight of patient)

### Model 2001-00XL: Frame white powder coated Model 2001-04XL: Frame silver powder coated

like Mod. 2000XL, but difference: Width of lying surface: 800 mm



Height adjustment by electromotive drive from 420 mm to 970 mm.

Time of electric level positioning from min. to max.: 18 sec.

The electric motor is controlled by 1 foot switch (possible option: Hand switch).

### Basic equipment:

Lying surface 2-sectioned:


- Head part stepless up- (+45°) and downwards (-35°)
- adjustable by gas spring
- Length head part: 500 mm
- Foot (back) part fixed. Length: 1.500 mm
- Not movable. Possibilities for mobility see list of extras, page 34.


Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 2.000 mm, Width: 700 mm Load capacity: 250 kg (max. weight of patient)

## Model 4001-00: Frame white powder coated Model 4001-04: Frame silver powder coated

like Mod. 4000, but difference: Width of lying surface: 800 mm





## Special equipment Model 1665 for tables Model 2000XL / 2001XL and 4000 / 4001:

Head part stepless up- and downwards adjustable by electric motor.

The motor is controlled by hand switch (standard).

## Some possibilities to provide the different alternatives of height adjustable tables series 2000 and 4000 with accessories:

 Side guards, - Storage baskets, - Upholstery in soft and elastic configuration, - Support for arm, - Stainless steel side rails, - Different possibilities of castors, - Paper roll holders under head part or at foot part, - Opening for nose and cushion to close....

Please have a look at our special equipment on page 34 in this catalogue.

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

## ...for Examination, Treatment

## - height adjustable





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)

The ease and convenience designed table Model 5000:

- Upholstery ergonomic configurated
- Upholstered in comfortable softness to support patients during treatment relaxed and ease
- Nice designed table, usage of large scale dimensioned oval and round pipe





Model 5000-04

### Model 5000-00, frame white powder coated Model 5000-04, frame silver powder coated

Height adjustment only by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch)

### **Basic equipment:**

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring. Length head part: 500 mm
- Foot (back) part fixed. Length: 1.500 mm
- Upholstery in high-quality comfortable softness Thickness of upholstery: 8 cm
- Upholstery ergonomic configurated, special: Widening in the middle of back part for the arms
- Not movable. Possibilities for mobility see list of extras, page 34. Measures:

Length: 2.000 mm, Height: Adjustable from 480 to 1.020 mm Width: Trapezoidal from 600 mm (begin head part) to 790 mm middle of back (foot) part to 650 mm end of back (foot) part Load capacity: 250 kg (max. weight of patient) Time of electric level positioning from min. to max.: 28 sec.



Model 2006-...XL and Model 2007-...XL

Measures: Length min.: 1.550 mm, Length max.: 1.950 mm Alternative 1: Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. Alternative 2: Height adjustment by electromotive drive.

The electric motor is controlled by 1 foot switch

(possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec.

### **Basic equipment:**

Lying surface 2-sectioned:

- Head part stepless adjustable from downwards vertical position (-90°) to upwards position (+30°) by gas spring. Length head part: 450 mm
- Foot (back) part fixed. Length: 1.500 mm
- Not movable. Possibilities for mobility see list of extras, page 34. Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood. Load capacity: 250 kg (max. weight of patient)

Both pictures with extras: Castors 100 mm diam., each with total lock (4x Mod. 1120)



Model 2006-00XL

Model	Width	Vertical height a	Frame colour	
2006-00XL/H	700 mm	Hydraulic pump:	540 - 960 mm	White powder painted
2006-04XL/H	700 mm	Hydraulic pump:	540 - 960 mm	Silver powder painted
2007-00XL/H	800 mm	Hydraulic pump:	540 - 960 mm	White powder painted
2007-04XL/H	800 mm	Hydraulic pump:	540 - 960 mm	Silver powder painted
2006-00XL	700 mm	Electric motor:	520 - 960 mm	White powder painted
2006-04XL	700 mm	Electric motor:	520 - 960 mm	Silver powder painted
2007-00XL	800 mm	Electric motor:	520 - 960 mm	White powder painted
2007-04XL	800 mm	Electric motor:	520 - 960 mm	Silver powder painted

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

### ...for Examination - height adjustable with variable length





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)



Both picture with extras: Castors 100 mm diam., each with total lock (4x Mod. 1120)

Model 2010...XL and Model 2011...XL: Length min.: 1.550 mm, Length max. 1.950 mm

Alternative 1: Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. Alternative 2: Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch). Time of electric level positioning from min. to max.: 22 sec.

#### **Basic equipment:**

Lying surface 3-sectioned:

- Head part stepless upwards adjustable by gas spring till +60°. Length head part: 450 mm
- Middle part fixed. Length: 1.050 mm
- Foot part stepless adjustable from downwards vertical position (-90°) to upwards position (+30°) by gas spring. Length foot part: 450 mm
- Not movable. Possibilities for mobility see list of extras, page 34. Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood. Load capacity (foot part vertical downwards): 250 kg (max. weight of

Model	Width	Vertical height a	Frame colour	
2010-00XL/H	700 mm	Hydraulic pump:	550 - 970 mm	White powder painted
2010-04XL/H	700 mm	Hydraulic pump:	550 - 970 mm	Silver powder painted
2011-00XL/H	800 mm	Hydraulic pump:	550 - 970 mm	White powder painted
2011-04XL/H	800 mm	Hydraulic pump:	550 - 970 mm	Silver powder painted
2010-00XL	700 mm	Electric motor:	530 - 970 mm	White powder painted
2010-04XL	700 mm	Electric motor:	530 - 970 mm	Silver powder painted
2011-00XL	800 mm	Electric motor:	530 - 970 mm	White powder painted
2011-04XL	800 mm	Electric motor:	530 - 970 mm	Silver powder painted



Model 2010-00XL/H with extras



Model 2008-...XL and Model 2009-...XL

Length min.: 1.790 mm, Length max. 2.190 mm

Alternative 1: Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. Alternative 2: Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec. **Basic equipment:** 

Lying surface 3-sectioned:

- Head part stepless up- (+50°) and downwards (-45°) adjustable by gas spring. Length head part: 450 mm
- Middle part fixed. Length: 1.290 mm
- Foot part stepless adjustable from downwards vertical position (-90°) to upwards position (+30°) by gas spring. Length foot part: 450 mm
- Not movable. Possibilities for mobility see list of extras, page 34. Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood. Load capacity (foot part vertical downwards): 250 kg (max. weight of

Model	Width	Vertical height a	Frame colour	
2008-00XL/H	700 mm	Hydraulic pump:	540 - 960 mm	White powder painted
2008-04XL/H	700 mm	Hydraulic pump:	540 - 960 mm	Silver powder painted
2009-00XL/H	800 mm	Hydraulic pump:	540 - 960 mm	White powder painted
2009-04XL/H	800 mm	Hydraulic pump:	540 - 960 mm	Silver powder painted
2008-00XL	700 mm	Electric motor:	520 - 960 mm	White powder painted
2008-04XL	700 mm	Electric motor:	520 - 960 mm	Silver powder painted
2009-00XL	800 mm	Electric motor:	520 - 960 mm	White powder painted
2009-04XL	800 mm	Electric motor:	520 - 960 mm	Silver powder painted



Model 2008-00XL

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

### ...for Examination, Treatment, - height adjustable





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)



Alternative 1: Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. Alternative 2: Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec.

### **Basic equipment:**

Lying surface 2-sectioned:

- Head part upwards adjustable +30° by 2 metal rachets, safety and for longevity. Length of head part: 550 mm
- Foot part fixed. Length: 1.400 mm
- Standard-Thickness of upholstery: 8 cm

Model 2600-...XL and Model 2605-...XL

- Movable on 4 castors, each with total lock
- (Other possibilities for mobility see list of extras, page 34)

Measures: L: 1.950 mm, Load capacity: 250 kg (max. weight of patient)

Picture with extras:

- 1 pce stainless steel side rail 1.300 mm long (Mod. 1420)

Model	Width	Vertical height a	adjustment by	Frame colour	
2600-01XL/H	650 mm	Hydraulic pump:	540 - 960 mm	Top frame chromed, underpart white	
2600-00XL/H	650 mm	Hydraulic pump:	520 - 940 mm	White powder painted	
2600-04XL/H	650 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted	
2605-01XL/H	800 mm	Hydraulic pump:	540 - 960 mm	Top frame chromed, underpart white	
2605-00XL/H	800 mm	Hydraulic pump:	520 - 940 mm	White powder painted	
2605-04XL/H	800 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted	
2600-01XL	650 mm	Electric motor:	520 - 960 mm	Top frame chromed, underpart white	
2600-00XL	650 mm	Electric motor:	500 - 940 mm	White powder painted	
2600-04XL	650 mm	Electric motor:	500 - 940 mm	Silver powder painted	
2605-01XL	800 mm	Electric motor:	520 - 960 mm	Top frame chromed, underpart white	
2605-00XL	800 mm	Electric motor:	500 - 940 mm	White powder painted	
2605-04XL	800 mm	Electric motor:	500 - 940 mm	Silver powder painted	



### Model 2150-...XL and Model 2155-...XL

Alternative 1: Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. Alternative 2: Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec.

### **Basic equipment:**

Lying surface 2-sectioned:

- Head part stepless upwards adjustable by gas spring till +75° as a back for seat positioning. Length: 750 mm
- Foot part fixed. Length: 1.200 mm
- Standard-Thickness of upholstery: 8 cm.
- Movable on 4 castors, each with total lock
- (Other possibilities for mobility see list of extras, page 34)

Measures: L: 1.950 mm, Load capacity: 250 kg (max. weight of patient)

Both pictures with extras: Castors 100 mm diam. with central locking system (Mod. 1180N)

Model	Width	Vertical height a	adjustment by	Frame colour
2150-00XL/H	650 mm	Hydraulic pump:	520 - 940 mm	White powder painted
2150-04XL/H	650 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted
2155-00XL/H	800 mm	Hydraulic pump:	520 - 940 mm	White powder painted
2155-04XL/H	800 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted
2150-00XL	650 mm	Electric motor:	500 - 940 mm	White powder painted
2150-04XL	650 mm	Electric motor:	500 - 940 mm	Silver powder painted
2155-00XL	800 mm	Electric motor:	500 - 940 mm	White powder painted
2155-04XL	800 mm	Electric motor:	500 - 940 mm	Silver powder painted

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

# ...for Examination, Treatment - height adjustable





Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)



Picture with extras: Castors 100 mm diam. with central locking system (Mod. 1180N)

#### Model 2610-...XL and Model 2615-...XL

**Alternative 1:** Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. **Alternative 2:** Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec.

### **Basic equipment:**

Lying surface 2-sectioned:

- Head part stepless up- (+40°) and downwards (-20°) adjustable by gas spring. Length head part: 550 mm
- Foot (back) part fixed. Length: 1.400 mm
- Standard-Thickness of upholstery: 8 cm
- Movable on 4 castors, each with total lock

(Other possibilities for mobility see list of extras, page 34)

Measures: L: 1.950 mm, Load capacity: 250 kg (max. weight of patient)

Model	Width	Vertical height adjustment by		Frame colour	
2610-01XL/H	650 mm	Hydraulic pump:	540 - 960 mm	Top frame chromed, underpart white	
2610-00XL/H	650 mm	Hydraulic pump:	520 - 940 mm	White powder painted	
2610-04XL/H	650 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted	
2615-01XL/H	800 mm	Hydraulic pump:	540 - 960 mm	Top frame chromed, underpart white	
2615-00XL/H	800 mm	Hydraulic pump:	520 - 940 mm	White powder painted	
2615-04XL/H	800 mm	Hydraulic pump:	520 - 940 mm	Silver powder painted	
2610-01XL	650 mm	Electric motor:	520 - 960 mm	Top frame chromed, underpart white	
2610-00XL	650 mm	Electric motor:	500 - 940 mm	White powder painted	
2610-04XL	650 mm	Electric motor:	500 - 940 mm	Silver powder painted	
2615-01XL	800 mm	Electric motor:	520 - 960 mm	Top frame chromed, underpart white	
2615-00XL	800 mm	Electric motor:	500 - 940 mm	White powder painted	
2615-04XL	800 mm	Electric motor:	500 - 940 mm	Silver powder painted	



Picture with extras: Castors 100 mm diam. with central locking system (Mod. 1180N) and Paper roll holder (Mod. 1265K)

### Model 2650-...XL and Model 2651-...XL

**Alternative 1:** Manual foot operated height adjustment by hydraulic pump. Pedals for foot-controlled height positioning are bothsides the table. **Alternative 2:** Height adjustment by electromotive drive. The electric motor is controlled by 1 foot switch (possible option: Hand switch) Time of electric level positioning from min. to max.: 22 sec.

### **Basic equipment:**

Lying surface 3-sectioned:

- Head part (L: 550 mm) upwards adjustable +30° by 2 metal rachets
- Middle part fixed. Length: 450 mm
- Foot part (L: 950 mm) stepless upwards adjustable by gas spring till 25°.
- Standard-Thickness of upholstery: 8 cm
- Movable on 4 castors, each with total lock
- (Other possibilities for mobility see list of extras, page 34)

Measures: L: 1.950 mm, Load capacity: 250 kg (max. weight of patient)

### Model 2655-...XL and Model 2656-...XL (without picture)

like Mod. 2650XL – 2651XL, but difference:

 Head part stepless up- (+40°) and downwards (-20°) adjustable by gas spring. Length head part: 550 mm

Model with head part adjust 0° bis + 30°	Model with head part adjust -20° bis + 40 °	Width	Vertical height a	adjustment by	Frame colour
2650-00XL/H	2655-00XL/H	650 mm	Hydraulic pump:	560 - 980 mm	White powder painted
2650-04XL/H	2655-04XL/H	650 mm	Hydraulic pump:	560 - 980 mm	Silver powder painted
2651-00XL/H	2656-00XL/H	800 mm	Hydraulic pump:	560 - 980 mm	White powder painted
2651-04XL/H	2656-04XL/H	800 mm	Hydraulic pump:	560 - 980 mm	Silver powder painted
2650-00XL	2655-00XL	650 mm	Electric motor:	540 - 980 mm	White powder painted
2650-04XL	2655-04XL	650 mm	Electric motor:	540 - 980 mm	Silver powder painted
2651-00XL	2656-00XL	800 mm	Electric motor:	540 - 980 mm	White powder painted
2651-04XL	2656-04XL	800 mm	Electric motor:	540 - 980 mm	Silver powder painted

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

### ...for Examination, Treatment

### - height adjustable by central positioned electric lifting column

Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or Skai® Pandoria Plus (20 colours)
- Height adjustable by central positioned electric lifting column (Therefore free space under top frame)
- Large cross section of column and precise guidance ensure tip-resistant operation
- Soft and smooth height adjustment (Almost constant speed burdened or unburdened)
- Covering of the column, made of 3 telescope sliding aluminium sections, silver anodized





Model 6000-04

### Model 6000-04

Underpart silver powder coated.

Lifting column aluminium, silver anodized.

### **Basic equipment:**

Height adjustable from 560 mm to 960 mm.

Time of electric level positioning from min. to max.: 24 sec.

**DEWERT** 

The electric lifting column is controlled by 1 foot switch (possible option: Hand switch)

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring.
- Length head part: 550 mm
- Foot (back) part fixed. Length: 1.400 mm
- Not movable

Possibilities for mobility see list of extras, page 34.

Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 1.950 mm, Width: 700 mm, Height: 560 – 960 mm

Load capacity: 250 kg (max. weight of patient)

#### Model 6001-04

like Mod. 6000-04, but difference: Width of lying surface: 800 mm

Both pictures with extras: Castors 75 mm diam., each with total lock (4x Mod. 1100) Note: Then height adjustment from 650 to 1.050 mm





Model 6050-04

### Model 6050-04

Underpart silver powder coated.

Lifting column aluminium, silver anodized.

### **Basic equipment:**

### Height adjustable from **560 mm to 960 mm**.

Time of electric level positioning from min. to max.: 24 sec.

The electric lifting column is controlled by 1 foot switch

(possible option: Hand switch)

Lying surface 2-sectioned:

- Head part (L: 750 mm) stepless up- +75° (as a back for seat positioning) and downwards (-25°) adjustable by gas spring.
- Head part in aesthetic half-round design
- Foot (back) part fixed. Length: 1.200 mm
- Complete lying surface is adjustable by electric motor to inclined position (Trendelenburg) till 20°
- Electric motor to adjust the lying surface is controlled by 1 hand switch Note: Adjustment of head part is independent of table height or positioning of lying surface

Possibilities for mobility as special equipment.

Standard-Thickness of upholstery: 6 cm.

Upholstery base: Laminated wood.

Measures: Length: 1.950 mm, Width: 700 mm, Height: 560 – 960 mm

Load capacity: 250 kg (max. weight of patient)

### Model 6051-04 (without picture)

like Mod. 6050-04, but difference: Width of lying surface: 800 mm

All tables are produced in our company in Bielefeld, Germany. On request the tables may be delivered in different measures of height adjustment, other measures of lying surface, other sectioned of upholstery or with additional functions.

Other frame colours, not column, are also possible. Please contact us.


- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

### ...for Examination, Treatment,

### - height adjustable, lying surface inclinable



- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)



Tables with electric motor must have a 230 V power connection for motor function.
 The motor is very powerful with constant speed burdened or unburdened.

A safety-lock-box to protect against unauthorized usage is necessary.

- Tables with height adjustment by hydraulic pump we recommend by predominant mobility usage

The height adjustment is possible without electric current

The height is manual adjustable by foot operated pump (repeated treading of the foot pedal)





Both pictures with extras:

- Castors 100 mm diam. with central locking system (Mod. 1180N) (Then height: 570 – 990 mm)
- Paper roll holder (Mod. 1270K)
- Special push handle (Mod. 1812)



Manual foot operated height adjustment by **hydraulic pump** from **550 mm to 970 mm**.

Pedals for foot-controlled height positioning are bothsides the table. **Basic equipment:** 

### sasic equipment:

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-35°) adjustable by gas spring. Length head part: 500 mm (Adjustment of head part is independent of table height or position of lying surface)
- Foot (back) part fixed. Length: 1.500 mm
- Complete lying surface is stepless adjustable by gas spring to inclined position (Trendelenburg) till 14°
- A push handle (= Option, Mod.-No. 1812) at foot part facilitate the operation of manual inclination (please see picture)
- Not movable. Possibilities for mobility see list of extras, page 34. Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 2.000 mm, Width: 700 mm Load capacity: 250 kg (max. weight of patient) Length overall (inclusive push handle): 2.150 mm

Model 2053-00XL/H: Frame white powder coated Model 2053-04XL/H: Frame silver powder coated

like Mod. 2052...XL/H, but difference: Width of lying surface: 800 mm



Both pictures with extras:
- Wheel-lifting/fixing-me-

- chanism (Mod. 1045)
   Holder for side guards stickable (Mod. 1030)
- Opening for nose and upholstery part to close (Mod. 1610 + 1620)



Additional: 2 side guards stickable (Mod. 1010)

Model 2052-00XL: Frame white powder coated Model 2052-04XL: Frame silver powder coated

Height adjustment by **electromotive drive** from **540 mm to 980 mm**.

Time of electric level positioning from min. to max.: 22 sec.

The electric motor is controlled by 1 foot switch

(possible option: Hand switch)

### Basic equipment:

Lying surface 2-sectioned:

- Head part stepless up- (+45°) and downwards (-20°) adjustable by gas spring. Length head part: 500 mm (Adjustment of head part is independent of table height or position of lying surface)
- Foot (back) part fixed. Length: 1.500 mm
- Complete lying surface is adjustable by electric motor to inclined position (Trendelenburg) till 30°
- The electric motor to adjust lying surface is controlled by 1 hand switch.
- Not movable. Possibilities for mobility see list of extras, page 34. Standard-Thickness of upholstery: 6 cm. Upholstery base: Laminated wood.

Measures: Length: 2.000 mm, Width: 700 mm Load capacity: 250 kg (max. weight of patient)

Model 2053-00XL: Frame white powder coated Model 2053-04XL: Frame silver powder coated like Mod. 2052, but difference:

Width of lying surface: 800 mm

- Mod. 1993: Safety box, fixed on the underpart, with 1 magnetic key (see picture on page 34)
- Mod. 1995: Emergency stop button, fixed on top frame, with 2 keys (see picture on page 34)

