...for Examination, Treatment

Please choose - without surcharge - between 2 different variants of leatherettes:

- Stamskin®Top (34 colours) or - Skai® Pandoria Plus (20 colours)





Model 16-01: Frame completely chromed Model 16-00: Frame white powder coated Model 16-04: Frame silver powder coated

Underpart of square tube 30x30x2 mm, not movable. Basic: With transverse bracing and levelling feet. Possibilities for mobility as special equipment. Note: Therefore we recommend longitudinal bracing. Lying surface 1-part:

Measures: L: 1.450 mm, W: 650 mm, H: 650 mm Standard-Thickness of upholstery: 6 cm Load capacity: 250 kg (max. weight of patient)

Model 17-01: Frame completely chromed Model 17-00: Frame white powder coated Model 17-04: Frame silver powder coated Like Model 16, but width of lying surface 800 mm


(Other measures on request)



Picture with extras: Paper roll holder (Mod. 1210)



Model 113-01: Frame completely chromed Model 113-00: Frame white powder coated Model 113-04: Frame silver powder coated Underpart of square tube 30x30x2 mm,



Possibilities for mobility as special equipment. Note: Therefore we recommend longitudinal bracing.

- Lying surface 3-sectioned:
 Head part (L: 550 mm) upwards adjustable (+30°) by 2 metal rachets.
- Middle part fixed. Length: 920 mm
- Foot part stepless from vertical downwards -90° till +10° upwards adjustable by gas spring. Length foot part: 480 mm

Measures: Length min.:1.600 mm / max.: 1.950 mm, W: 650 mm, H: 800 mm. Standard-Thickness of upholstery: 8 cm. Load capacity (foot part vertical downwards): 250 kg (max. weight of patient)

Model 114-01: Frame completely chromed Model 114-00: Frame white powder coated Model 114-04: Frame silver powder coated Like Model 113, but width of lying surface 800 mm

Model 121-01: Frame completely chromed Model 121-00: Frame white powder coated Model 121-04: Frame silver powder coated

Underpart of round pipe 32x2 mm, not movable. Basic: With transverse and longitudinal bracing and levelling feet.

Possibilities for mobility as special equipment. Lying surface 3-sectioned:

- Head part upwards adjustable (+30°) by 2 metal rachets, safety and for longevity. Length of head part: 550 mm
- Middle part fixed. Length: 900 mm
- Foot part stepless from vertical downwards -90° till +25° upwards adjustable by gas spring. Length foot part: 500 mm

Measures: Length min.:1.550 mm / max.: 1.970 mm, W: 650 mm, H: 800 mm. Standard-Thickness of upholstery: 8 cm.

Load capacity (foot part vertical downwards): 250 kg (max. weight of patient)



Model 122-01: Frame completely chromed Model 122-00: Frame white powder coated Model 122-04: Frame silver powder coated

Underpart of round pipe 32x2 mm, not

Basic: With transverse and longitudinal bracing and levelling feet. Possibilities for mobility as special equipment.

- Lying surface 3-sectioned: Head and foot part stepless from vertical downwards -90° till +25° upwards adjustable by gas spring.
- Length head / foot part: 450 mm (each)
- Middle part fixed. Length: 1.120 mm

Measures: Length min.: 1.300 mm / max.: 2.020 mm, W: 650 mm, H: 800 mm. Standard-Thickness of upholstery: 8 cm.

Load capacity (head and foot vertical part downwards): 250 kg (max. weight of patient)

