

KaWe COLORSCOP Duo

DOPPELKOPF-STETHOSKOP | DUAL-HEAD STETHOSCOPES


KaWe COLORSCOP Duo Stethoskope bringen Farbe in Ihren Alltag. Wir verwenden für die Herstellung dieser Stethoskope nur hochwertige Materialien. Durch die erstklassige Verarbeitung eignet sich dieses Stethoskop hervorragend für den täglichen Gebrauch in der Klinik oder im Pflegedienst für alle Untersuchungen. Das Doppelkopf-Bruststück kann von der Membran- auf die Trichterseite durch Drehen am Schlauch umgeschaltet werden. Die austauschbare Membrane ist glasfaserverstärkt und daher besonders sensibel. Die Schalltrichter sind aus Metall gefertigt und gewährleisten so eine gute Resonanz durch eine optimale Schallleitung. Die Auswahl an acht verschiedenen Farben schließt ein Verwechseln unter den Anwendern nahezu aus.

KaWe COLORSCOP Duo stethoscopes bring colour into your everyday routine. We use high-quality materials for the production of these stethoscopes. Due to its first-rate workmanship, this stethoscope is suited extremely well for daily use in clinics or by nurses for all types of examinations. The double-sided chest piece can be switched from diaphragm to bell mode by twisting the tube. The exchangeable diaphragm is reinforced with fibre-glass and is therefore especially sensitive. The sound bell is made of metal and therefore guarantees a good resonance because of optimal sound conduction. The available selection of eight different colours almost completely eliminates the possibility of a user mistaking their instrument for that of another.

gut | good ausgezeichnet | excellent



krankes Herz Sick heart (16-2000 Hz)	
Gefäßgeräusche Vascular sounds (100-2000 Hz)	
Blutdruckmessung Blood pressure (100-2000 Hz)	
Lungengeräusche Pulmonary sounds (100-2000 Hz)	
Darmgeräusche Bowel sounds (16-1000 Hz)	



134

KaWe COLORSCOP Duo

- das bunte Doppelkopf-Stethoskop, umschaltbar, leichtes Eigengewicht
- Membran Ø 43,8 mm
- Gesamtlänge ca. 78 cm
- Gewicht: ca. 110 gr.

- the colourful dual-head stethoscope, reversible, lightweight
- diaphragm Ø 43.8 mm
- total length approx. ca. 78 cm
- weight: approx. 110 gr.



Farbbeispiele | Example of available colours

rot red		REF 06.22100.014
schwarz black		REF 06.22100.024
blau blue		REF 06.22100.034
grün green		REF 06.22100.044



gelb yellow		REF 06.22100.054
lila violet		REF 06.22100.094
türkis turquoise		REF 06.22100.124
pink pink		REF 06.22100.134

ERSATZTEILE SPARE PARTS

Ohrlivern | Ear tips

- schwarz
- Ø 5,0 mm
- VE = 10 St.

- black
- Ø 5.0 mm
- PU = 10 items



REF 06.91120.021